EN 3 Cook® Chorionic Villus Sampling Set

Instructions For Use

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Bruksanvisning

COOK® CHORIONIC VILLUS SAMPLING SET

PROFESSIONAL LABELING AND DIRECTIONS FOR USE

COOK® CHORIONIC VILLUS SAMPLING SET STERILE SINGLE USE/DISPOSABLE

DESCRIPTION

The Cook Chorionic Villus Sampling (CVS) Set is designed to provide the physician with a means to obtain a sample of tissue from the chorion frondosum transcervically.

The Cook CVS Set is sterile and disposable and consists of a flexible plastic tube. The catheter has a formed distal tip and a proximally mounted Luer connector. Integral to the catheter is a removable stainless steel stylet which is visible under ultrasound.

INDICATIONS FOR USE

The Cook CVS Set is indicated for use in obtaining chorionic tissue samples transcervically for the purpose of prenatal diagnosis of genetic abnormalities **during weeks 10-12 of pregnancy.**

CAUTION: Sterile if the package is unopened or undamaged. Do not use if package is broken.

CONTRAINDICATIONS

The use of the Cook CVS Set is contraindicated when one or more of the following conditions exist:

- Signs of intrauterine infection or gonorrheal pelvic infection;
- Atypical anatomy of the cervix or uterus which makes a successful sampling unlikely, or situations requiring extreme flexion of the uterus;
- Known or suspected uterine or cervical malignancy including an unresolved abnormal pap smear.

ACUTE ADVERSE EFFECTS OF THE DEVICE

Adverse effects reported within the first week post procedure period with use of the Cook CVS Set during a study of 4537 patients included the following:

1. **Maternal Bleeding:** Maternal spotting and/or bleeding was reported in 439 of 4537 sampled patients (9.7%). This can be potentially minimized by passage of a catheter through the cervical os **no more than two times per pregnancy**, and by utilizing a new catheter for each attempt, and by careful passage of the catheter utilizing ultrasound guidance.
2. **Cramping:** 26 of 4537 patients (0.6%)
3. **Fluid Leakage:** 6 of 4537 patients (0.1%)
4. **Pain:** 43 of 4537 patients (0.9%)
5. **Flu-like Symptoms:** 3 of 4537 patients (0.07%)
6. **Placental Separation:** 3 of 4537 patients (0.07%)
7. **Hypoglycemic Seizure:** 1 of 4537 patients (0.02%) experienced an acute hypoglycemic seizure following the procedure. The patient was a known diabetic and the seizure was not felt to have been directly procedure related.

LONG TERM ADVERSE EFFECTS FOLLOWING CVS

Adverse effects occurring in the Cook CVS clinical study in the 4363 patients with long term follow-up data included the following:

1. **Spontaneous Abortion:** Spontaneous abortion was reported in 163 of the 4363 patients available for follow up (3.7%). This adverse effect can potentially be minimized by ultrasound examination immediately prior to the procedure to ensure viability of the pregnancy and by utilizing safe sampling techniques.
Additionally, this risk can be minimized by using a new sterile catheter for each cervical os pass and by making a maximum of two passages through the os per pregnancy in order to attempt to obtain the tissue sample.
2. **Oromandibular/Limb Dystrophies (OM/LD):** Oromandibular/Limb Dystrophy (OM/LD) is a condition in which a child is born with shortened and/or missing digits or limbs, and/or malformed lower jaw and mouth. Some authors have suggested a link between CVS and OM/LD (See warnings for further discussion).

OM/LD was observed in 3 of 4363 (0.068%) pregnancies during the clinical investigation.

POTENTIAL ADVERSE EFFECTS

The CVS procedure in general has the potential for producing the conditions listed below. However, none of these potential adverse reactions were observed during the clinical investigation of the Cook CVS Set in 4537 patients.

1. **Infection:** The potential risk for infection can be minimized by adherence to appropriate aseptic techniques, proper pre-procedural vaginal cleaning and sterilization of the equipment used for the procedure. Also, the risk of infection can be potentially decreased by using a new catheter for each cervical os pass and by making no more than two passes through the cervical os per pregnancy.
2. **Septic Shock:** There is a potential risk of morbidity, including death, from septic shock secondary to an infection introduced from the CVS procedure. This possibility can be minimized by adherence to appropriate aseptic techniques, proper pre-procedural vaginal cleaning and sterilization of the equipment used in the procedure. Also, using a new sterile catheter for each pass through the cervical os and making no more than two passes through the cervix per pregnancy will potentially minimize the risk.
3. **Rh Sensitization:** There is a risk of Rh sensitization in Rh negative women. This possibility can be minimized by administering a RhoGAM injection to all Rh negative women after the procedure.

4. **Fetal Bleeding:** This potential risk may be minimized by careful ultrasound examination prior to the procedure and by cautious passage of the catheter under ultrasound guidance during the procedure.
5. **Uterine Perforation:** This potential risk can be minimized by careful passage of the catheter under ultrasound guidance.
6. **Intrauterine Growth Retardation:** There is a theoretical risk of intrauterine growth retardation as a result of the CVS procedure.
7. **Intrauterine Death:** An increased risk of intrauterine death is known to be associated with significant interruption of the fetal blood supply by chorionic (placental) damage or intrauterine infection. Adherence to CVS safety procedures will minimize the potential risk in terms of the procedure.

WARNINGS

- Risk of Oromandibular/Limb Dystrophies (OM/LD)

Oromandibular/Limb Dystrophy (OM/LD) is a condition in which a child has shortened and/or missing digits or limbs, and/or malformed lower jaw and mouth. Some authors have suggested a link between CVS and OM/LD.

Data exist from clinical studies which suggest an increased risk of OM/LD associated with CVS procedures performed **prior to 10 weeks gestation**.¹ Incidents of OM/LD have also been reported when the CVS procedure was performed from 10-12 weeks gestation.² These reports describe CVS procedures done trans-abdominally and transcervically. One study of births in the general population showed OM/LD defects to occur at a rate of approximately 5.4 per 10,000 live births (0.054%).³

OM/LD was observed in 3 of 4363 (0.068%) pregnancies with known outcomes during the clinical investigation of the Cook CVS Set.

The three clinical study cases were first published by Burton et. al.⁴ A second article by the same authors⁵ stated that a fourth case of OM/LD had been observed following the use of the Cook CVS Set. However, after investigation, it was found that this particular patient had *not* had CVS performed using the Cook CVS Set and the published report was in error.⁶

Due to the risk of OM/LD, it is essential to perform a detailed ultrasound study of the gestational sac and uterus before the procedure in order to assure a correct assessment of the gestational age of the fetus, as well as fetal and uterine anatomical locations, and to restrict the number of sampling attempts per pregnancy.

- Number of Attempts to Obtain a Sample of Tissue

The clinical studies have shown that the total fetal loss rate increased directly with the number of passes made with a catheter through the cervical os in order to obtain a tissue sample. **It is strongly recommended that no more than two catheter passes through the cervical os be made per pregnancy in order to attempt to obtain tissue.** In the clinical study of the Cook CVS Set, the loss rate after the CVS procedure was 4.2% with one attempt, 7.0% after two attempts and 10.4% after three attempts.

- Single Use Catheter

Each sampling pass attempt through the cervical os must be done with a new sterile catheter.

- Post-Procedure Symptoms

In the acute post procedure period, the patient must be advised to report **immediately** to her health care provider all abnormal conditions such as abdominal cramping, flu-like symptoms, fever, vaginal discharge, or severe bleeding.

PRECAUTIONS

- Multiple Pregnancies

Insufficient clinical data are available concerning the safety and effectiveness of the procedure in multiple pregnancy situations (e.g., twins). The procedure is therefore not recommended in this circumstance.

- Intrauterine Devices

Insufficient clinical data are available concerning the safety and effectiveness of chorionic villus sampling in the presence of an intrauterine device. Therefore, the procedure is not recommended for patients with an intrauterine device.

PATIENT INFORMATION

- Patient Counseling

Prior to the procedure, the physician, nurse, genetic counselor, or other trained health care professional must provide the patient with the document entitled "Patient Information: Chorionic Villus Sampling During Pregnancy". *In addition*, the patient should be given thorough genetic counseling, be encouraged to ask questions, and review the procedure options available for fetal genetic diagnosis. The risk/benefit situation for individual patients must be thoroughly explained prior to the procedure.

- Patient Evaluation and Clinical Considerations

- a. The chorionic villus sampling technique requires a concerted team effort: a sonologist for ultrasound guidance, an obstetrician to perform the villus sampling, and a geneticist to perform chromosomal and biochemical analysis. Each have critical responsibilities in order to obtain satisfactory results. These individuals, along with the nursing and genetic counseling staff, form a team which must be capable of discussing the risks and benefits of the procedure with the patient.
- b. A physical examination for acute pelvic inflammatory disease is necessary prior to the procedure. A cervical culture for *Neisseria gonorrhoeae* is advisable. Any patient demonstrating signs of active pelvic infection must not receive the procedure. In addition, an initial ultrasound evaluation is advisable to determine viability of pregnancy, age of gestation, number of fetuses present, relationship of placenta to cervical canal, the presence of fibroids, adnexal masses, or any other abnormalities.

- c. Due to the clinical study data which suggest an increased risk of Oromandibular/Limb Dystrophies associated with CVS procedures performed **prior to 10 weeks gestation**, it is critical to perform a detailed ultrasound study of the gestational sac and uterus **on the day of the procedure** to recheck the assessment of the gestational age of the fetus as well as the fetal and uterine positions.
- d. The patient should be told to anticipate some minor bleeding and/or cramping during the two weeks following the procedure. But, **if these symptoms continue or increase in severity**, they should be reported **immediately** to her primary health care provider.
- e. Careful ultrasound screening for developmental anomalies (e.g. Oromandibular/Limb Dystrophies) should be performed at **16-18 weeks gestation for all patients**.

• Patient Information Literature

The "Patient Information: Chorionic Villus Sampling During Pregnancy" document must be provided to the patient. The literature is provided to help the patient and health care provider discuss the risks and benefits of the CVS procedure and review the clinical study data. **However, do not rely on the patient information literature alone to be able to explain all the risks and benefits of this complex issue to the patient.**

SUMMARY OF COOK® CVS SET CLINICAL FINDINGS

- Data from 4537 patients demonstrated greater overall success with increased experience. It is recommended that physicians with **suitable CVS training and experience perform the sampling procedure**.
- Clinical studies have shown that to minimize patient risks, **no more than two cervical os passes per pregnancy**, utilizing a **new sterile catheter** with each pass, should be performed due to the associated potential increased risks for infection and spontaneous abortion.
- The CVS procedure is **best performed between the 10th and 12th post-menstrual week**.
- The most common reasons for performing the CVS procedure in the Cook CVS Set study population (4537 patients) were:
 - a. Maternal age of 33 years, or older (in this study protocol, patients aged 33 or older were allowed to enroll);
 - b. Previous delivery of an aneuploid fetus or infant;
 - c. Parental carrier of chromosomal translocation; or
 - d. Maternal risk of X-linked or autosomal recessive disorder.

CAUTION: Federal law restricts this device to sale, distribution, use by, or on the order of a physician with appropriate training and experience.

REFERENCES

1. Firth, H.V., et.al., "Severe Limb Abnormalities After Chorion Villus Sampling at 55-66 Days' Gestation," *Lancet*, 1991;337:762-763.
2. Miny, P., et. al., "Limb-reduction Defects and Chorion Villus Sampling (Letter to Editor)," *Lancet*, 1991;337:1423-1424.
3. Froster-Iskenius, U.G., et. al., "Limb Reduction Defects in Over One Million Consecutive Live Births," *Teratology*, 1989;39:127-135.
4. Burton, B.K., et. al., "Limb Anomalies Associated With Chorionic Villus Sampling," *Obstetrics and Gynecology*, 1992 May;79(5):726-730.
5. Burton B.K., Schulz C.J., Burd L.I., "Spectrum of limb disruption defects associated with chorionic villus sampling," *Pediatrics*, 1993 May;91(5): 989-93.
6. Burton, B.K. Correspondence. June 11, 1993.
7. Blakemore K., Filkins K., Luthy D.A., Platt L.D., et. al. "Cook obstetrics and gynecology catheter multicenter chorionic villus sampling trial: comparison of birth defects with expected rates," *American Journal of Obstetrics and Gynecology*, 1993 Oct;169(4):1022-6.

DIRECTIONS FOR USE

STERILE TECHNIQUE SHOULD BE OBSERVED

1. **It is essential that patients be adequately informed about the potential risks and benefits of undergoing the CVS procedure.** *Prior to the procedure* the patient should have received extensive counseling, including genetic counseling, as well as a thorough explanation of the risks and benefits of the CVS procedure, and a full description of the events that are to be anticipated during and after the procedure.
2. A thorough ultrasound study of the placenta, fetus, and uterus should be performed immediately prior to the procedure in order to reassess the gestational age and the position of the uterus, cervix and placenta to facilitate catheter insertion. **The ultrasound measurement should provide the distance from the internal os to the sampling site and help determine the optimal curve, if any, required for the assembled catheter and stylet that will be used for sampling.**
3. After the patient has been positioned in the lithotomy position, cleanse the external genitalia with an appropriate antiseptic solution.
4. Partially fill a 20 ml syringe with 5 ml nutrient media.
5. Internally swab the vagina using ring forceps, cotton balls and an appropriate antiseptic solution.
6. Insert a sterile speculum and visualize the cervix.
7. Cleanse the cervix with an appropriate antiseptic solution and dry excess solution from cervix.
8. The cervix may be grasped with a sterile tenaculum in order to apply gentle traction as required.
9. Using sterile gloves and adhering to sterile technique, remove the Cook CVS Set from the package. Remove the protective sheath from the catheter. With stylet *in-situ*, manually curve the catheter as necessary to direct its insertion based on the ultrasound findings.
10. Slowly insert the catheter through the endocervix to the sampling site under ultrasound guidance. **Care must be taken to avoid entering the amniotic sac.**

11. Slowly withdraw the stylet without displacing the catheter.
12. Attach the prefilled syringe to the catheter Luer mount.
13. With the syringe, apply intermittent gentle suction to the biopsy site while slowly withdrawing the catheter. The degree of necessary negative pressure and motion should be determined by viewing the procedure on the ultrasound monitor. Very gentle suction should be used when the catheter is actually placed on the chorionic membrane. Slightly more suction (up to 5 ml of negative pressure) may be applied when the catheter is **clearly in the center of the chorion** according to the ultrasound view.
14. Withdraw the catheter while maintaining slight negative pressure on the syringe.
15. While adhering to sterile technique, flush the contents of the syringe into a sterile petri dish.
16. The specimen should be quickly evaluated under a dissecting microscope to determine the adequacy of the sample. A sample size of 10 mg to 30 mg (*approximately 0.1 cc of sample will contain 1 mg of tissue*) is generally adequate for subsequent genetic analysis.
17. A second sterile catheter should be introduced into the cervical os *only* if the prior sampling pass was judged to contain an insufficient quantity of villi to permit the desired testing.

NOTE: it is important that no catheter be passed through the cervical os more than one time and that no more than two catheter passes be made through the cervical os per pregnancy. Making three or more passes through the cervical os has been associated with a significantly increased risk of pregnancy loss.

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REFERENCE

- Firth, H.V., et al., "Severe Limb Abnormalities After Chorion Villus Sampling at 55-66 Days' Gestation," *Lancet*, 1991;337:762-763.
- Miny, P., et. al., "Limb-reduction Defects and Chorion Villus Sampling (Letter to Editor)," *Lancet*, 1991;337:1423-1424.
- Froster-Iskenius, U.G., et. al., "Limb Reduction Defects in Over One Million Consecutive Live Births," *Teratology*, 1989;39:127-135.
- Burton, B.K., et. al., "Limb Anomalies Associated With Chorionic Villus Sampling," *Obstetrics and Gynecology*, 1992 May;79(5):726-730.
- Burton B.K., Schulz C.J., Burd L.I., "Spectrum of limb disruption defects associated with chorionic villus sampling," *Pediatrics*, 1993 May;91(5): 989-93.
- Burton, B.K. Correspondence. June 11, 1993.
- Blakemore K., Filkins K., Luthy D.A., Platt L.D., et. al. "Cook obstetrics and gynecology catheter multicenter chorionic villus sampling trial: comparison of birth defects with expected rates," *American Journal of Obstetrics and Gynecology*, 1993 Oct;169(4):1022-6.

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LITTERATUR

1. Firth, H.V., et al., "Severe Limb Abnormalities After Chorion Villus Sampling at 55-66 Days' Gestation," *Lancet*, 1991;337:762-763.
2. Miny, P., et. al., "Limb-reduction Defects and Chorion Villus Sampling (Letter to Editor)," *Lancet*, 1991;337:1423-1424.
3. Froster-Iskenius, U.G., et. al., "Limb Reduction Defects in Over One Million Consecutive Live Births," *Teratology*, 1989;39:127-135.
4. Burton, B.K., et. al., "Limb Anomalies Associated With Chorionic Villus Sampling," *Obstetrics and Gynecology*, 1992 May;79(5):726-730.
5. Burton B.K., Schulz C.J., Burd L.I., "Spectrum of limb disruption defects associated with chorionic villus sampling," *Pediatrics*, 1993 May;91(5): 989-93.
6. Burton, B.K. Correspondence. June 11, 1993.
7. Blakemore K., Filkins K., Luthy D.A., Platt L.D., et. al. "Cook obstetrics and gynecology catheter multicenter chorionic villus sampling trial: comparison of birth defects with expected rates," *American Journal of Obstetrics and Gynecology*, 1993 Oct;169(4):1022-6.

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1. Firth, H.V., et.al., "Severe Limb Abnormalities After Chorion Villus Sampling at 55-66 Days' Gestation," *Lancet*, 1991;337:762-763.
2. Miny, P, et. al., "Limb-reduction Defects and Chorion Villus Sampling (Letter to Editor)," *Lancet*, 1991;337:1423-1424.
3. Froster-Iskenius, U.G., et. al., "Limb Reduction Defects in Over One Million Consecutive Live Births," *Teratology*, 1989;39:127-135.
4. Burton, B.K., et. al., "Limb Anomalies Associated With Chorionic Villus Sampling," *Obstetrics and Gynecology*, 1992 May;79(5):726-730.
5. Burton B.K., Schulz C.J., Burd L.I., "Spectrum of limb disruption defects associated with chorionic villus sampling," *Pediatrics*, 1993 May;91(5): 989-93.
6. Burton, B.K. Correspondence. June 11, 1993.
7. Blakemore K., Filkins K., Luthy D.A., Platt L.D., et. al. "Cook obstetrics and gynecology catheter multicenter chorionic villus sampling trial: comparison of birth defects with expected rates," *American Journal of Obstetrics and Gynecology*, 1993 Oct;169(4):1022-6.

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REFERENCIAS

1. Firth, H.V., et al., "Severe Limb Abnormalities After Chorion Villus Sampling at 55-66 Days' Gestation," *Lancet*, 1991;337:762-763.
2. Miny, P., et al., "Limb-reduction Defects and Chorion Villus Sampling (Letter to Editor)," *Lancet*, 1991;337:1423-1424.
3. Froster-Iskenius, U.G., et al., "Limb Reduction Defects in Over One Million Consecutive Live Births," *Teratology*, 1989;39:127-135.

4. Burton, B.K., et. al., "Limb Anomalies Associated With Chorionic Villus Sampling," *Obstetrics and Gynecology*, 1992 May;79(5):726-730.
5. Burton B.K., Schulz C.J., Burd L.I., "Spectrum of limb disruption defects associated with chorionic villus sampling," *Pediatrics*, 1993 May;91(5): 989-93.
6. Burton, B.K. Correspondence. June 11, 1993.
7. Blakemore K., Filkins K., Luthy D.A., Platt L.D., et. al. "Cook obstetrics and gynecology catheter multicenter chorionic villus sampling trial: comparison of birth defects with expected rates," *American Journal of Obstetrics and Gynecology*, 1993 Oct;169(4):1022-6.

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BIBLIOGRAPHIE

1. Firth, H.V., et al., "Severe Limb Abnormalities After Chorion Villus Sampling at 55-66 Days' Gestation," *Lancet*, 1991;337:762-763.
2. Miny, P., et al., "Limb-reduction Defects and Chorion Villus Sampling (Letter to Editor)," *Lancet*, 1991;337:1423-1424.
3. Froster-Iskenius, U.G., et al., "Limb Reduction Defects in Over One Million Consecutive Live Births," *Teratology*, 1989;39:127-135.
4. Burton, B.K., et al., "Limb Anomalies Associated With Chorionic Villus Sampling," *Obstetrics and Gynecology*, 1992 May;79(5):726-730.
5. Burton B.K., Schulz C.J., Burd L.I., "Spectrum of limb disruption defects associated with chorionic villus sampling," *Pediatrics*, 1993 May;91(5): 989-93.
6. Burton, B.K. Correspondence. June 11, 1993.
7. Blakemore K., Filkins K., Luthy D.A., Platt L.D., et al. "Cook obstetrics and gynecology catheter multicenter chorionic villus sampling trial: comparison of birth defects with expected rates," *American Journal of Obstetrics and Gynecology*, 1993 Oct;169(4):1022-6.

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1. Firth, H.V., et al., “Severe Limb Abnormalities After Chorion Villus Sampling at 55-66 Days’ Gestation,” *Lancet*, 1991;337:762-763.
2. Miny, P., et. al., “Limb-reduction Defects and Chorion Villus Sampling (Letter to Editor),” *Lancet*, 1991;337:1423-1424.

3. Froster-Iskenius, U.G., et. al., "Limb Reduction Defects in Over One Million Consecutive Live Births," *Teratology*, 1989;39:127-135.
4. Burton, B.K., et. al., "Limb Anomalies Associated With Chorionic Villus Sampling," *Obstetrics and Gynecology*, 1992 May;79(5):726-730.
5. Burton B.K., Schulz C.J., Burd L.I., "Spectrum of limb disruption defects associated with chorionic villus sampling," *Pediatrics*, 1993 May;91(5): 989-93.
6. Burton, B.K. Correspondence. June 11, 1993.
7. Blakemore K., Filkins K., Luthy D.A., Platt L.D., et. al. "Cook obstetrics and gynecology catheter multicenter chorionic villus sampling trial: comparison of birth defects with expected rates," *American Journal of Obstetrics and Gynecology*, 1993 Oct;169(4):1022-6.

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- Firth, H.V., et.al., "Severe Limb Abnormalities After Chorion Villus Sampling at 55-66 Days' Gestation," *Lancet*, 1991;337:762-763.
- Miny, P., et. al., "Limb-reduction Defects and Chorion Villus Sampling (Letter to Editor)," *Lancet*, 1991;337:1423-1424.
- Froster-Iskenius, U.G., et. al., "Limb Reduction Defects in Over One Million Consecutive Live Births," *Teratology*, 1989;39:127-135.
- Burton, B.K., et. al., "Limb Anomalies Associated With Chorionic Villus Sampling," *Obstetrics and Gynecology*, 1992 May;79(5):726-730.
- Burton B.K., Schulz C.J., Burd L.I., "Spectrum of limb disruption defects associated with chorionic villus sampling," *Pediatrics*, 1993 May;91(5): 989-93.
- Burton, B.K. Correspondence. June 11, 1993.
- Blakemore K., Filkins K., Luthy D.A., Platt L.D., et. al. "Cook obstetrics and gynecology catheter multicenter chorionic villus sampling trial: comparison of birth defects with expected rates," *American Journal of Obstetrics and Gynecology*, 1993 Oct;169(4):1022-6.

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MANUFACTURER
COOK INCORPORATED
750 Daniels Way
Bloomington, IN 47404 U.S.A.

www.cookmedical.com
© COOK 2019

EC REPRESENTATIVE
COOK IRELAND LTD.
O'Halloran Road
National Technology Park
Limerick, Ireland

2019-09
T_J-CVS_REV3

Chorionic Villus Sampling During Pregnancy	EN 3
Odběr choriových klků v průběhu těhotenství	CS 4
Chorion-villus-prøvetagning under graviditeten	DA 6
Chorionzottenbiopsie während der Schwangerschaft	DE 8
Δειγματοληψία χοριακής λάχνης κατά την κύηση	EL 10
Muestreo de vellosidades coriónicas durante el embarazo	ES 12
Prélèvement de villosités chorioniques pendant la grossesse	FR 14
Chorionboholy-mintavétel terhesség idején	HU 16
La villocentesi durante la gravidanza	IT 18
Chorionvillusbiospie tijdens de zwangerschap	NL 20
Morkakebiopsi (CVS) under graviditet	NO 22
Biopsja kosmówki w czasie ciąży	PL 24
Colheita de amostra de vilosidade coriônica durante a gravidez	PT 26
Provtagning av chorionvilli under graviditet	SV 28
孕期绒毛取样	ZH 30

CE
0123

CHORIONIC VILLUS SAMPLING DURING PREGNANCY

PATIENT INFORMATION

WHAT IS CHORIONIC VILLUS SAMPLING (CVS)?

Chorionic Villus Sampling (CVS) is a procedure used to get a sample of tissue from the placenta for genetic testing of the fetus. The sample is examined to help identify birth disorders such as Down's Syndrome.

CVS is performed by sampling the chorionic villi tissues in the womb of an expectant woman. Chorionic tissue surrounds part of the sac where the fetus develops. The diagram in this brochure shows the location of the chorionic villi tissue.

Chorionic tissue originally starts from the same cell as the fetus. So taking a sample of the tissue is similar to taking a sample cell from the fetus.

ARE YOU A CANDIDATE FOR CVS?

You should ask your doctor whether or not the procedure is needed in your particular case.

Candidates for CVS are women who are concerned about an inherited disorder in their family or have a child with a birth disorder, or are concerned for other medical reasons that the fetus may have a birth disorder.

Having a child with a birth disorder may be due to several reasons, but the most common ones include:

- mother's age of 35 or older (risk of having an affected child increases with age)
- a family history of mental retardation
- multiple miscarriages
- inherited (genetic) conditions
- multiple birth problems, or
- disorders in previous children.

WHEN IS CVS DONE?

The procedure is most safely done between the 10th and 12th weeks of pregnancy. Counting from the first day of the last menstrual period will give you the number of weeks you are in your pregnancy.

CVS must not be done on women with certain medical conditions. Discuss this issue with your doctor.

HOW IS THE CVS PROCEDURE PERFORMED?

Before and during the CVS procedure, you will have an ultrasound examination. Ultrasound is the use of sound waves to form a picture of a developing fetus and is painless.

Your doctor will confirm the age of the fetus before the CVS procedure is done. Then during the CVS procedure ultrasound is used to allow the doctor to locate and identify the safest site for the sampling.


You will be asked to lie on your back with your legs placed in the stirrups, just like a pelvic exam. The vagina is cleaned and held in an open position with an instrument known as a speculum.

A thin plastic tube or catheter is gently inserted through the cervix (opening of the womb) to reach the chorion tissue. Ultrasound helps the doctor guide the catheter to the proper position.

You may feel some cramping during the procedure. This is minor in nature and usually no anesthesia is needed.

When the catheter is properly located, a small sample of tissue is removed. The tissue is then sent to the laboratory for analysis. Results are usually available within a week.

You should be able to return home following the procedure.

TABLE

WHAT ARE THE RISKS ASSOCIATED WITH CVS?

Miscarriage

There is always a risk of miscarriage even if CVS is not done. On average about 1 in 20 (4.7%) pregnancies end in miscarriage by the 12th week. When CVS was done using the Cook catheter in a clinical study, the results showed 3.7% of pregnant women had a miscarriage. Thus, it is difficult to say if any of these miscarriages were due to the CVS procedure or not.

You should also know that miscarriage rates vary with the age of the mother, her past pregnancies and medical histories. For example, among mothers 35 years of age and older, the rate of miscarriage increases to 8.7% by the 12th week.

Oromandibular/Limb Dystrophies

Oromandibular/Limb Dystrophies (OM/LD) is a condition in which the child is born with shortened fingers, toes, arms or legs. The child may also have a malformation (abnormality) of the jaw or mouth.

OM/LD occurs in some patients when CVS is done during 10 to 12 weeks gestation. An increased risk of OM/LD may be linked to CVS when performed **before** 10 weeks of gestation.

OM/LD occurs in the general population about 5 to 6 times for every 10,000 live births. In a study of the Cook CVS Set, about 7 per 10,000 children were born with some shortened fingers, toes or both. None of these had malformations of the jaw or mouth.

To check for OM/LD, a follow-up ultrasound examination should be done. You should make an appointment with your doctor during the 16th to 18th week of pregnancy.

Other Complications

A clinical study using the Cook CVS Set found that complications included: maternal bleeding, cramping, fluid leakage, pain, flu-like symptoms, and placental separation.

IS THERE AN ALTERNATIVE TO CVS?

Amniocentesis is an alternative prenatal procedure. It checks for the same disorders as CVS plus disorders of the brain and spinal cord.

Amniocentesis is done between the 16th and 18th week when fewer miscarriages are seen. The risk of miscarriage following amniocentesis is 1% lower than CVS; this figure takes into account the ages of the mother and the fetus.

You may want to ask your doctor about amniocentesis as an alternative procedure.

DECIDING ABOUT HAVING THE CVS PROCEDURE

To undergo CVS is an important decision. The decision to have this or any test is up to you and your partner. You are not required to have this procedure. You must understand the importance of an early diagnosis, the risks of the procedure, and the potential risk of having an abnormal child.

When considering the CVS procedure, you should contact your health care provider as soon as possible. The procedure is done between the 10th and 12th week of pregnancy, so timing is critical.

If you are considering having CVS, you should discuss it with your doctor or a genetic counselor. They will ask you for information about your family's medical history and explain CVS to you in more detail. If you decide to have the procedure, they will help you in scheduling an appointment.

HOW CAN I FIND OUT MORE ABOUT CVS?

If you have further questions, please discuss them with your health care provider or genetic counselor.

WARNING

Following CVS, **IMMEDIATELY REPORT** any abnormal conditions to your doctor. Abnormal conditions should be properly monitored. Examples of abnormal conditions are:

- flu-like symptoms
- fever
- abdominal cramping
- pain
- severe bleeding, or
- vaginal discharge.

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Cook Medical Europe Ltd