EN 3	Olcott Torque Device Instruction for Use
CS 3	Rotátor Olcott Návod k použití
DA 4	Olcott momenthåndtag Brugsanvisning
DE 5	Olcott Steuerungshandgriff Gebrauchsanweisung
EL 5	Συσσκευή στρεψής Olcott Οδηγίες χρήσης
ES 6	Dispositivo de torque Olcott Instrucciones de uso
FR 7	Dispositif de contrôle Olcott Mode d'emploi
HU 7	Olcott forgatóeszköz Használati utasítás
IT 8	Dispositivo per il controllo direzionale Olcott Istruzioni per l'uso
NL 9	Olcott torsiehulpmiddel Gebruiksaanwijzing
NO 9	Olcott momentanordning Bruksanvisning
PL 10	Rotator Olcott Instrukcja użycia
PT 11	Dispositivo de torção Olcott Instruções de utilização
SV 11	Olcott vridningsanordning Bruksanvisning
ZH 12	Olcott扭控器 使用说明


Figure A


Figure B

OLCOTT TORQUE DEVICE

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DEVICE DESCRIPTION

- Unique single-handed operation permits rapid position changes during procedural use

INTENDED USE

The Olcott Torque Device accepts .015 to .045 wire guide. Used to grip and torque control the Roadrunner® PC Wire Guide, or similar hydrophilic wire guides.

CONTRAINDICATIONS

There are no known contraindications.

WARNINGS

None known

PRECAUTIONS

None known

INSTRUCTIONS FOR USE

1. With thumb-slide in a forward position, place Olcott Torque Device over proximal end of wire guide to desired location (**Figure A**).
2. Pull thumb-slide back until wire guide is held securely in place (**Figure B**).

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MANUFACTURER

COOK INCORPORATED
750 Daniels Way
Bloomington, IN 47404 U.S.A.

www.cookmedical.com
© COOK 2019

EC REPRESENTATIVE

Cook Medical Europe Ltd
O'Halloran Road
National Technology Park
Limerick, Ireland

2019-09
T_OCD_REV1

EN 3	Wire Guides Instruction for Use
CS 4	Vodící dráty Návod k použití
DA 5	Kateterledere Brugsanvisning
DE 6	Führungsdrähte Gebrauchsanweisung
EL 8	Συρμάτινοι οδηγοί Οδηγίες χρήσης
ES 9	Guías Instrucciones de uso
FR 11	Guides Mode d'emploi
HU 12	Vezetődrótok Használati utasítás
IT 13	Guide Istruzioni per l'uso
NL 15	Voerdraden Gebruiksaanwijzing
NO 16	Ledevaiere Bruksanvisning
PL 17	Prowadniki Instrukcja użycia
PT 19	Fios guia Instruções de utilização
SV 20	Ledare Bruksanvisning
ZH 21	导丝 使用说明

WIRE GUIDES

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

- Wire Guide (*size and configuration may vary*)

INTENDED USE

Used to gain ureteral access, to establish a tract, and to assist in the placement, replacement, and exchange of medical devices during urological and gynecologic procedures. The HiWire, BiWire and Roadrunner Wire Guides may also be used for catheter positioning and exchange in a tortuous or kinked ureter, traversing a large stone in route to the kidney, or in cases demanding enhanced control and high radiopacity.

CONTRAINDICATIONS

None known

WARNINGS

None known

PRECAUTIONS

- Manipulation of wire guide requires appropriate imaging control. Use caution not to force or over manipulate the wire guides when gaining access.
- When using wire guide through a metal cannula/needle, use caution as damage may occur to outer coating.
- When exchanging or withdrawing an instrument over the wire guide, secure and maintain the wire guide in place under fluoroscopy in order to avoid unexpected wire guide displacement.
- These wire guides are not intended for PTCA use.

INSTRUCTIONS FOR USING A WIRE GUIDE:

1. Using aseptic technique, remove the wire from its outer packaging and place in the sterile field.
2. Remove the wire guide from the holder.
3. If using movable core wire guide, adjust mandril position by advancing or withdrawing the mandril within the wire guide to achieve the required degree of flexibility.
4. Advance the flexible end of the wire guide into the device and position to the desired anatomical location.
 - Endhole size and length of the device must be taken into consideration to ensure proper fit between wire guide and device.
5. After wire guide use is complete, remove from patient and discard.

INSTRUCTIONS FOR ACTIVATING HYDROPHILIC COATING:

The hydrophilic coating on the wire guide is activated by immersion in sterile water or sterile saline solution.

1. Prior to using the wire guide, fill a syringe with sterile water or sterile saline solution and attach it to the flushing port on the wire guide.
2. Inject enough solution to wet the wire guide surface entirely. This will activate the hydrophilic coating.
 - For Roadrunner® PC Wire Guide Double Flexible, the pink end of wire guide is not hydrophilically coated or intended for insertion into body.

- Hydrophilic coated wires are very slippery when wet. Always maintain control of the wire guide when manipulating it through any device.
- For optimal performance, rehydrate the hydrophilic coated wire guide after exposure to ambient environment or after extended use, replace it with a new hydrophilic coated wire guide.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

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