- EN**
3 **Olcott Torque Device**
Instruction for Use
- CS**
3 **Rotátor Olcott**
Návod k použití
- DA**
4 **Olcott momenthåndtag**
Brugsanvisning
- DE**
5 **Olcott Steuerungshandgriff**
Gebrauchsanweisung
- EL**
5 **Συσσκευή στρεψής Olcott**
Οδηγίες χρήσης
- ES**
6 **Dispositivo de torque Olcott**
Instrucciones de uso
- FR**
7 **Dispositif de contrôle Olcott**
Mode d'emploi
- HU**
7 **Olcott forgatóeszköz**
Használati utasítás
- IT**
8 **Dispositivo per il controllo direzionale Olcott**
Istruzioni per l'uso
- NL**
9 **Olcott torsiehulpmiddel**
Gebruiksaanwijzing
- NO**
9 **Olcott momentanordning**
Bruksanvisning
- PL**
10 **Rotator Olcott**
Instrukcja użycia
- PT**
11 **Dispositivo de torção Olcott**
Instruções de utilização
- SV**
11 **Olcott vridningsanordning**
Bruksanvisning
- ZH**
12 **Olcott扭控器**
使用说明


Figure A


Figure B

OLCOTT TORQUE DEVICE

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DEVICE DESCRIPTION

- Unique single-handed operation permits rapid position changes during procedural use

INTENDED USE

The Olcott Torque Device accepts .015 to .045 wire guide. Used to grip and torque control the Roadrunner® PC Wire Guide, or similar hydrophilic wire guides.

CONTRAINDICATIONS

There are no known contraindications.

WARNINGS

None known

PRECAUTIONS

None known

INSTRUCTIONS FOR USE

1. With thumb-slide in a forward position, place Olcott Torque Device over proximal end of wire guide to desired location (**Figure A**).
2. Pull thumb-slide back until wire guide is held securely in place (**Figure B**).

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WIRE GUIDES

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

Wire Guide (*size and configuration may vary*)

INTENDED USE

Used to gain ureteral access, to establish a tract, and to assist in the placement, replacement, and exchange of medical devices during urological and gynecologic procedures. The Roadrunner wire guides may also be used for catheter positioning and exchange in a tortuous or kinked ureter, traversing a large stone in route to the kidney, or in cases demanding enhanced control and high radiopacity.

CONTRAINDICATIONS

None known

WARNINGS

- Do not pull the wire guide through a needle with a sharp edge when the wire is not aligned straight with the end of the needle. Extreme caution should be observed any time the wire guide is used within a needle. Manipulation or removal of the wire guide through a needle may result in destruction/separation of the polymer jacket of the wire requiring retrieval.
- Do not force the device during the procedure. Excessive force can lead to the wire guide fracturing, bending, or kinking. Additionally, excessive force can cause the wire guide to become entrapped. If any resistance is felt or if the wire behavior and/or location seem improper, stop manipulating the wire and determine the cause. Failure to exercise proper caution may result in bending, kinking, separation of the wire guide jacket, damage to ancillary devices or damage to the urinary system. If necessary, remove the wire guide and ancillary device or scope as a complete unit to avoid possible damage and/or complications.
- Do not pass the wire through a damaged scope or ancillary device. A compromised lumen may cause damage or separation to the polymer jacket of the wire guide.

PRECAUTIONS

- Manipulation of wire guide requires appropriate imaging control. Use caution not to force or over manipulate the wire guides when gaining access.
- Avoid retracting the wire guide through a metal cannula or scope when the devices are not aligned straight, as the metal device can damage the polymer jacket of the wire guide.
- The surface of the wire guide is not lubricious unless it is wet. Before taking it out of its holder, activate the wire guide hydrophilic coating with sterile water or sterile saline. Removing the wire guide before activating the coating may result in damage to the polymer jacket of the wire guide.
- End hole size and length of device must be taken into consideration to ensure a proper fit between the wire guide and the procedural device.
- When exchanging or withdrawing an instrument over the wire guide, secure and maintain the wire guide in place under fluoroscopy in order to avoid unexpected wire guide displacement.
- These wire guides are not intended for PTCA use.

DEVICE PREPARATION

1. Using aseptic technique, remove the wire guide dispenser with the wire guide from its outer packaging and place the wire guide dispenser with the wire guide in the sterile field.

2. Prior to removing the hydrophilic wire guide from its dispenser, fill a syringe with sterile water or sterile saline solution and inject the solution into the luer lock hub at the end of the dispenser.
3. Inject enough solution to fill the dispenser coil. The solution will completely cover the wire guide surface and activate the hydrophilic coating. (Fig 1.)
4. Remove the hydrophilic wire guide from its dispenser by gently withdrawing the wire guide from its tip.

NOTE: If the hydrophilic wire guide cannot be easily be removed from its dispenser, then inject more solution into the dispenser and attempt withdrawal again.

NOTE: For Roadrunner PC Wire Guide Double Flexible, the pink end of the wire guide is not hydrophilically coated or intended for insertion into the body.

INSTRUCTIONS FOR USE

1. Flush the procedural device with saline solution before and during use to ensure smooth movement of the hydrophilic wire guide within the device.
2. Advance the flexible end of the wire guide into the procedural device and position the wire guide in the desired anatomical location.

NOTE: For optimal performance, rehydrate the hydrophilic-coated wire guide after exposure to ambient environments or if the lubricity of the wire guide is diminished. Replace the wire guide with a new hydrophilic-coated wire guide if wire guide lubricity remains diminished.

3. After wire guide use is complete, remove the wire guide from the patient and discard the wire guide in an appropriate waste receptacle.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

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© COOK 2019
Cook Medical Europe Ltd