- EN**
3 **Advance® 35LP Low Profile**
PTA Balloon Dilatation Catheter
Instructions for Use
- CS**
4 **Balónkový dilatační katetr s nízkým**
profilem Advance® 35LP pro PTA
Návod k použití
- DA**
6 **Advance® 35LP lavprofils-PTA**
ballondilataionskateter
Brugsanvisning
- DE**
8 **Advance® 35LP PTA-Balldilataionskatheter**
mit flachem Profil
Gebrauchsanweisung
- EL**
10 **Καθετήρας διαστολής PTA Advance® με**
μπαλόνη, χαμηλού προφίλ, 35LP
Οδηγίες χρήσης
- ES**
12 **Catéter balón de dilatación para ATP de perfil**
bajo Advance® 35LP
Instrucciones de uso
- FR**
14 **Cathéter de dilatation à ballonnet et à profil**
bas pour ATP Advance® 35LP
Mode d'emploi
- HU**
15 **Advance® 35LP keskeny PTA ballonos**
tágítókatéter
Használati utasítás
- IT**
17 **Catetere a palloncino per dilatazione per PTA**
Advance® 35LP a basso profilo
Istruzioni per l'uso
- NL**
19 **Advance® 35LP laagprofiel-PTA-**
ballondilatatiekatheter
Gebruiksaanwijzing
- NO**
21 **Advance® 35LP lavprofils PTA-**
ballongdilatasjonskateter
Bruksanvisning
- PL**
23 **Niskoprofilowy balonowy cewnik**
rozszerzający Advance® 35LP do PTA
Instrukcja użycia
- PT**
25 **Cateter de dilatação com balão, de baixo**
perfil, para PTA Advance® 35LP
Instruções de utilização
- SV**
27 **Advance® 35LP PTA-ballongdilataionskateter**
med låg profil
Bruksanvisning

ADVANCE® 35LP LOW PROFILE PTA BALLOON DILATATION CATHETER

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Advance 35LP Low Profile PTA Balloon Dilatation Catheter is a double-lumen catheter with a balloon near its distal tip. The catheter consists of two independent lumens, which are labeled "DISTAL" and "BALLOON." The distal lumen extends the length of the catheter and is used for placement of wire guides. The balloon lumen is used to expand the balloon. Inscribed on the tip of the manifold are the balloon diameter (mm) and the balloon length (cm).

The balloon is manufactured from an extra-thinwall, high-strength, minimally-compliant material. **Particular care should be taken in handling the balloon to prevent damage. It will inflate to the indicated size parameters when utilizing proper pressure recommendations.** Adhere to balloon inflation pressure parameters indicated in the **Compliance Card insert**. Refer to label for further information. Use of a pressure gauge is recommended to monitor inflation pressures.

Platinum/iridium, radiopaque markers are positioned on the shaft within the balloon to enable visualization of the catheter/balloon under fluoroscopy. The catheter is compatible with .035 inch (0.89 mm) wire guides.

INTENDED USE

The Advance 35LP Low Profile PTA Balloon Dilatation Catheter is indicated for percutaneous transluminal angioplasty (PTA) of lesions in peripheral arteries including iliac, renal, popliteal, infrapopliteal, femoral and iliofemoral as well as obstructive lesions of native or synthetic arteriovenous dialysis fistulae.

CONTRAINDICATIONS

None known

WARNINGS

- Do not exceed rated burst pressure.¹ Rupture of balloon may occur. Adhere to balloon inflation pressure parameters in the **Compliance Card insert**. Over-inflation may cause rupture of the balloon, with resultant damage to the vessel wall. Use of a pressure gauge is recommended to monitor inflation pressures.
- Do not use a power injector for balloon inflation or injection of contrast medium through catheter lumen marked "DISTAL." Rupture may occur.

¹The burst pressure data was analyzed using factors for a one-sided tolerance to determine with 95% confidence that 99.9% of these balloons would not burst at or below the calculated rated burst pressure.

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of vascular access sheaths, angiographic catheters and wire guides should be employed.
- The Advance 35LP Low Profile PTA Balloon Dilatation Catheter has been designed for introduction into the vascular system utilizing percutaneous (Seldinger) technique. Alternatively, an introducer sheath may be utilized. Refer to label information for the appropriate introducer sheath size.
- The balloon is constructed of a heat-sensitive material. Do not heat or attempt to shape the catheter tip.
- In heavily scarred access sites, use of an introducer sheath is recommended.
- Manipulate the catheter using fluoroscopic control.
- Use only the recommended balloon inflation medium. Never use air or any gaseous medium to inflate the balloon.
- Use the catheter prior to the expiration date specified on the package.
- The catheter is not intended for the delivery of stents.
- All stents should be deployed in accordance with the manufacturer's indications and instructions for use.

POTENTIAL ADVERSE EVENTS

- Air embolism
- Aneurysm
- Arrhythmias
- Arteriovenous fistula
- Death
- Drug reactions, allergic reaction to contrast medium
- Hemorrhage or hematoma
- Hypo/hypertension
- Infection and pain at insertion site
- Systemic embolization
- Vascular thrombosis
- Vessel dissection, perforation, rupture or injury
- Vessel spasm

INSTRUCTIONS FOR USE

Balloon Preparation

1. Choose a balloon appropriate to lesion length and vessel diameter.

2. Upon removal from package, inspect the catheter to ensure no damage has occurred during shipping.
3. Remove protective balloon sleeve from balloon. Do not discard protective sleeve after initial introduction of catheter into the vascular system; it may be helpful in wrapping balloon prior to subsequent insertions.
4. Prepare balloon lumen with standard contrast-saline mixture as follows:
 - a. Fill a syringe of appropriate size with 1:1 mixture of contrast medium and normal saline.
 - b. Attach the syringe to the hub of the catheter marked "BALLOON" and pull back to apply negative pressure.
 - c. Release pressure, allowing negative pressure to draw mixture into balloon lumen.
 - d. Detach syringe, leaving a meniscus of mixture on the hub of the balloon lumen.
 - e. Prepare inflation device in standard manner and purge to remove all air from syringe and tubing.
 - f. Attach inflation device to the balloon lumen, ensuring no air bubbles remain at connection.
 - g. Pull negative pressure on the inflation device.

Balloon Introduction and Inflation

1. Flush the catheter lumen labeled "DISTAL" using heparinized saline solution.
2. Apply negative pressure to lumen labeled "BALLOON" prior to introduction. Advance the balloon dilatation catheter counter-clockwise over a pre-positioned .035 inch (0.89 mm) wire guide.
3. Under fluoroscopy, advance the balloon to the lesion site. Carefully position the balloon across the lesion using both the distal and proximal radiopaque balloon markers.

NOTE: If resistance is met while advancing the balloon dilatation catheter, determine the cause and proceed with caution.
4. Inflate balloon to desired pressure. **Adhere to recommended balloon inflation pressures. (See Compliance Card insert.)**
5. If balloon pressure is lost and/or balloon rupture occurs, deflate balloon and remove balloon and sheath as a unit.

Balloon Deflation and Withdrawal

1. Completely deflate the balloon using an inflation device or syringe. Allow adequate time for the balloon to deflate.

NOTE: Balloons with large diameters and/or longer lengths may require longer deflation times.
2. Deflate the balloon by pulling vacuum on the inflation syringe or inflation device. Maintain vacuum on the balloon and withdraw the catheter. Upon catheter withdrawal, a gentle counterclockwise rotation of the catheter will assist balloon rewrap, minimizing trauma to the percutaneous entry site.
3. If resistance is met during withdrawal, apply negative pressure with a larger syringe before proceeding. If resistance continues, remove balloon and sheath as a unit.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2015
Bloomington, IN 47404 U.S.A.