EN 3	Cook 810 Set Instructions for Use
CS 4	Cook 810 Set Návod k použití
DA 5	Cook 810 Set Brugsanvisning
DE 6	Cook 810 Set Gebrauchsanweisung
EL 7	Cook 810 Set Οδηγίες χρήσης
ES 8	Cook 810 Set Instrucciones de uso
FR 9	Cook 810 Set Mode d'emploi
HU 10	Cook 810 Set Használati utasítás
IT 11	Cook 810 Set Istruzioni per l'uso
NL 12	Cook 810 Set Gebruiksaanwijzing
NO 13	Cook 810 Set Bruksanvisning
PL 14	Cook 810 Set Instrukcja użycia
PT 15	Cook 810 Set Instruções de utilização
SV 16	Cook 810 Set Bruksanvisning

COOK 810 SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner). Please read all instructions before using this device.

DEVICE DESCRIPTION

The Cook 810 Set consists of a dilator and a sheath. The sheath has a hub on the proximal end to act as a handle and aide the introduction of the set dilator and other instruments. The distal tip of the sheath is tapered down to provide a smooth transition to the dilator. The dilator is able to pass through the lumen of the sheath. The two ends of the dilator are tapered down to provide a smooth transition to the wire guide. When the sheath and dilator are used coaxially together, they are known as an "Introducer".

NOTE: Placement of the dilator should be performed over a wire guide not to exceed 0.038" in diameter.

INTENDED USE

The Cook 810 Set is used for retrograde or antegrade ureteral access and to facilitate wire guide exchange.

CONTRAINDICATIONS

This product is contraindicated in the presence of conditions which create unacceptable risk during endourological procedures.

WARNINGS

Do not force the device during the procedure; proceed with caution if excessive resistance is encountered.

PRECAUTIONS

This product is intended for use by physicians trained and experienced in urological procedures. Standard techniques for placement of ureteral introducers should be employed.

POTENTIAL ADVERSE EVENTS

Complications associated with endourological procedures may include, but are not limited to: bleeding, urinary tract perforation, stricture, and urinary retention.

INSTRUCTIONS FOR USE

1. Using standard access techniques (retrograde or antegrade), establish wire guide position.
2. Backload the dilator over the wire guide. Introduce dilator into urinary tract over the wire guide to desired location. If there is difficulty passing a stricture with the dilator, remove the dilator and retry using the other end of the dilator.
3. Once dilator location is achieved, backload the sheath over the wire guide and dilator. Introduce sheath into urinary tract to desired location.
4. While maintaining sheath and wire guide position, completely withdraw the dilator. The sheath is now ready to accept an additional wire through the lumen of the sheath.
5. Once patient intervention with other devices is complete, remove all instruments and then remove the sheath.
6. Upon completion of the procedure, dispose of the device per institutional guidelines for biohazardous medical waste.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

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Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd