EN
3

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ATLAS WIRE™ STONE EXTRACTOR AND CAPTURA® HELICAL STONE EXTRACTOR

Read all instructions carefully. Failure to properly follow the information provided may lead to the device not performing as intended or injury to the patient.

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

The Atlas Wire™ Stone Extractor and The Captura® Helical Stone Extractor are 3- or 4-wire stainless steel baskets available in a range of French sizes, lengths and basket diameters. An endoscopic cap and inserter are also included.

Performance Characteristics

- Decreased material thickness at shaft end to increase flexibility and increase scope deflection and irrigation flow.
- Device shaft consisting of inner stainless steel coil to increase torque response.
- Length of product allows targets in the kidney to be manipulated and removed.
- Stainless steel braid incorporated into sheath material for resistance to kinking.
- Ergonomic handle design allows for use with one hand with open and close directions noted.
- Includes endoscopic cap to resist backflow of fluid from the working channel of the endoscope around the extractor.
- Includes inserter to assist insertion of the extractor into the working channel of the endoscope.

Device Compatibility

- The 3.2 Fr extractors are compatible with endoscope working channels with a minimum inner diameter of 3.6 Fr (1.2 mm) and a maximum length of 80 cm.
- 1.7 Fr and 2.8 Fr extractors are compatible with endoscope working channels with a minimum inner diameter of 3.0 Fr (1.0 mm) and a maximum length of 80 cm.
- The endoscopic cap can be connected to any standard female luer lock fitting or bulb-type endoscopic fitting.
- The inserter tip is compatible with the endoscopic cap and allows for insertion of the extractor.
- These devices are compatible with sterile saline and sterile water.

RPN	FR	Length (CM)	Basket Diameter / Configuration
AE-417115	1.7	115	1.7 cm / 4 Wire
AE-428115	2.8	115	1.7 cm / 4 Wire
CE-328115	2.8	115	1.0 cm / 3 Wire
CE-428115	2.8	115	1.0 cm / 4 Wire
AE-432115	3.2	115	1.4 cm / 4 Wire
CE-332115	3.2	115	1.0 cm / 3 Wire
CE-432115	3.2	115	1.0 cm / 4 Wire

Patient Population

These devices are not intended to be used in any specific target population unless outlined on the labelling, in the instructions, and/or in promotional materials. Patient population could include any person requiring treatment according to the device's intended use and as directed by a qualified physician (or properly licensed practitioner).

Intended User

These devices are intended for physicians or healthcare practitioners with proper authorization to perform medical procedures (e.g., licensing, qualifications, credentials) per their local administrative and regulatory guidelines.

Contact with Body Tissue

These devices interact with the structures and tissues comprising the urinary tract, which may include the urethra, bladder, ureters, and kidneys.

Operating Principle

The Atlas Wire Stone Extractor is manufactured with a layered material shaft consisting of a polytetrafluoroethylene core with stainless steel braid and polyimide underneath a layer of fluorinated ethylene propylene in 1.7, 2.8, and 3.2 Fr diameters and 115 cm in length. The extractors feature a four-wire construction in a flat wire design for maximum stone retention. The Atlas Wire Stone Extractor is comprised of the 4-wire extractor, an optional endoscopic cap and an optional inserter. The Atlas Wire Stone Extractor is intended to be introduced into the urinary tract through the working channel of an endoscope. An endoscopic cap is provided with the Atlas Wire Stone Extractor that can be connected to the endoscope working channel to resist fluid backflow around the extractor. An inserter is also provided that may be used to facilitate insertion of the extractor into the endoscopic cap. The distal end of the extractor features a large bulb-shaped basket with four widely spaced wires. The basket diameter and wire spacing help to maintain a clear working field which allows for the manipulation of larger stones or fragments. The basket is controlled by a handle that facilitates the opening and closing of the basket. Sliding the thumb lever in the "OPEN" direction opens the basket, while sliding the thumb lever in the "CLOSE" direction closes the basket. Once the extractor is in position alongside or beyond the stone, the basket is opened and manipulated to encapsulate the stone. The basket can then be closed around the stone and removed from the patient along with the endoscope.

The Captura Helical Stone Extractors are manufactured from stainless steel wires in 3 and 4 –wire configurations. In the open position the basket expands 10 mm wide and 22 and 25 mm in length, when in the closed position the device diameter measures are 2.8 and 3.2 Fr. The shaft is manufactured with layered material consisting of a polytetrafluoroethylene core with stainless steel braid and polyimide underneath a layer of fluorinated ethylene propylene. The Captura Helical Stone Extractors have a length of 115 cm with a Uni-Dex Handle, optional endoscopic cap and optional inserter. The Captura Helical Stone Extractor is intended to be introduced into the urinary tract through the working channel of an endoscope. An endoscopic cap is provided with the Captura Helical Stone Extractor that can be connected to the endoscope working channel to resist fluid backflow around the extractor. An inserter is also provided that may be used to facilitate insertion of the extractor into the endoscopic cap. The distal end of the extractor features a three- or four-wire helical basket. The basket is controlled by a handle that facilitates the opening and closing of the basket. Sliding the thumb lever in the "OPEN" direction opens the basket, while sliding the thumb lever in the "CLOSE" direction closes the basket. Once the extractor is in position alongside or beyond the stone, the basket is opened and manipulated to encapsulate the stone. The basket can then be closed around the stone and removed from the patient along with the endoscope.

INTENDED USE

These devices are used for stone manipulation and removal in the urinary tract through the working channel of an endoscope during cystoscopic or ureteroscopic procedures by actuating a stainless-steel wire basket.

DEVICE AND CLINICAL BENEFITS

To manipulate and remove stones from the urinary tract through the working channel of an endoscope.

INDICATIONS FOR USE

These devices are used to assist in extraction of stones (urolithiasis).

CONTRAINDICATIONS

No contraindications.

WARNINGS

- The single use device is not designed for re-use. Attempts to reprocess (re-sterilize) and/or to re-use may lead to chemical contamination, device failure, and/or transmission of disease.
- Do not use the device if the sterile packaging is damaged or unintentionally opened before use.
- These devices are designed for the removal of stones no larger than their maximum basket diameter according to the table located in the device compatibility section of this instruction.

PRECAUTIONS

- This device is conductive. Avoid contact with any electrified instrument.
- Enclose the device in the sheath before removing from the tray/holder.
- User should be familiar with and have experience in urological endoscopic surgery.
- Assess stone prior to instrument deployment to ensure that the stone is not too large to be removed through anatomy.
- If resistance is encountered while attempting to remove stone, release the stone. Do not exert excessive force on the device.
-  This symbol on the label indicates that the device contains cobalt (Co) at a level above 0.1% w/w, which is a substance that is toxic to reproduction and a carcinogen (Class IB). However, the device contains stainless-steel alloy containing cobalt, which does not cause an increased risk of cancer or adverse reproductive effect according to current scientific evidence.

POTENTIAL ADVERSE EVENTS

- Bleeding e.g. Hematoma, Hematuria, Hemorrhage
- Death
- Deep venous thrombosis
- Fever
- Fluid overload
- Hypotension
- Hypothermia
- Inflammation e.g. Systemic inflammatory response syndrome
- Infection e.g. Bacteremia, Bacteriuria, Peri-renal Abscess, Pneumonia, Pyelonephritis, Sepsis, Urinary Tract Infection
- Nausea
- Pain
- Pulmonary Embolism
- Tissue, organ or structure injury e.g. Arteriovenous fistula, False tract, Infundibular stenosis, Pleural injury, Pseudoaneurysm, Renal collecting system injury, Selective arterial embolization, Surrounding organs/tissues injury, Urinary fistula
- Ureteral injury e.g. Ureteral avulsion, Ureteral mucosal injury, Ureteral mucosal tear, Ureteral obstruction, Ureteral perforation, Ureteral stricture
- Urinary tract issues e.g. Dysuria, Extrarenal stone migration, Fluid collection, Hydronephrosis, Steinstrasse, Urinary frequency, Urinary leakage, Urinary retention, Urinary urgency

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Keep dry and away from sunlight. Upon removal from package, inspect the product to ensure no damage has occurred.

INSPECTION OF DEVICE

Visually inspect the device thoroughly including all levels of the packaging (as applicable) to verify that there is no damage prior to use. Visually inspect and confirm that the integrity of the sterile barrier has not been compromised in any way.

INSTRUCTIONS FOR USE

1. Move the thumb lever on the handle to the "Close" position to close the basket and remove the extractor from its packaging.
2. The inserter may be placed securely into the endoscopic cap on the scope, as desired.
3. Then carefully insert the retracted extractor into the scope.
4. Under direct visualization, locate the stone and advance the extractor to a location just proximal to it. Slide the thumb lever on the handle to the "Open" position to deploy the basket near the stone and maneuver the device until it is within the basket.
5. Once correctly positioned, draw the basket down around the stone, and slide the thumb lever into the "Close" position.
6. Maintaining position of the stone within the basket, retract the scope and basket.

DISPOSAL OF DEVICES

After the procedure, this device may be contaminated with potentially infectious substances of human origin and should be disposed of in accordance with institutional guidelines.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device, this should be reported to Cook Medical and this shall also be reported to the competent authority where the device was used.

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© COOK 2025
Cook Medical Europe Ltd
William A. Cook Australia Pty Ltd