COOK

MEDICAL

- EN**
3 **Amplatz Renal Dilator Sets**
Instructions for Use
- CS**
3 **Soupravy renálních dilatátorů Amplatz**
Návod k použití
- DA**
4 **Amplatz-nyredilatatorsæt**
Brugsanvisning
- DE**
5 **Amplatz Nierendilatator-Sets**
Gebrauchsanweisung
- EL**
6 **ΣΕΤ νεφρικών διαστολέων Amplatz**
Οδηγίες χρήσης
- ES**
7 **Equipos de dilatadores renales Amplatz**
Instrucciones de uso
- FR**
8 **Set de dilateurs rénaux d'Amplatz**
Mode d'emploi
- HU**
9 **Amplatz renalis dilatátorkészletek**
Használati utasítás
- IT**
10 **Set con dilatatori renali Amplatz**
Istruzioni per l'uso
- NL**
11 **Amplatz renale dilatatorsets**
Gebruiksaanwijzing
- NO**
12 **Amplatz nyredilatatorsett**
Bruksanvisning
- PL**
13 **Zestawy rozszerzacza nerkowego Amplatz**
Instrukcja użycia
- PT**
14 **Conjuntos de dilatador renal Amplatz**
Instruções de utilização
- SV**
15 **Amplatz njurdilatatorset**
Bruksanvisning

T _ A R D S 2 _ R E V 0

AMPLATZ RENAL DILATOR SETS

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner). Please read all instructions before using this device.

DEVICE DESCRIPTION

Amplatz Renal Dilator Sets contain:

- an 8.0 French radiopaque catheter, tapered on both ends
- various fascial dilators
- a 10.0 French introducer sheath
- 30 cm long tapered dilators, available in various French sizes
- Radiopaque PTFE sheaths, available in various French sizes and lengths

INTENDED USE

Amplatz Renal Dilator Sets are intended for nephrostomy tract dilation and the placement of a sheath.

CONTRAINDICATIONS

This device is contraindicated in the presence of conditions which create unacceptable risk during percutaneous nephrostomy.

WARNINGS

Sterilized by Ethylene Oxide gas. Do not use if sterile barrier is damaged. The device is intended for one-time use.

PRECAUTIONS

- Imaging guidance is required for use of this device.
- This product is intended for use by physicians trained and experienced in urologic techniques. Standard techniques for placement of wire guides, dilators, and catheters should be employed.

POTENTIAL ADVERSE EVENTS

- Tissue trauma
- Bleeding
- Perforation
- Infection
- Urine leakage
- Fistula
- Pain

INSTRUCTIONS FOR USE

NOTE: Prior to use, immerse AQ® dilator in sterile water or isotonic saline to allow the hydrophilic surface to absorb water and become lubricious. This will ease placement under standard conditions.

NOTE: Prior to use of this device, you must gain access to the renal pelvis via wire guide placement.

1. Under fluoroscopic guidance, progressing from the smallest to the largest size, sequentially pass the 6.0 French to 10.0 French dilators over the wire guide. Maintain the wire guide's position. **NOTE:** The 6.0 French to 10.0 French dilators are tapered to fit a .038 inch diameter wire guide.
2. Pass the 8.0 French radiopaque catheter over the wire guide.
3. Sequentially dilate from 12.0 French dilator to the largest appropriate size dilator/sheath to be used. **NOTE:** Dilation with 12.0 French and larger dilators should only be done over 8.0 French radiopaque catheter.
4. Load the Amplatz sheath over the dilator.
5. While maintaining the dilator's position, advance the sheath into the calyx.
6. Use fluoroscopy to confirm location of the sheath in the calyx.
7. While holding the sheath in place, remove the dilator.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2019
Bloomington, IN 47404 U.S.A.