
- EN**  
5 **Amplatz Ureteral Stent Set**  
Instructions for Use
- CS**  
7 **Souprava ureterálního stentu Amplatz**  
Návod k použití
- DA**  
9 **Amplatz ureterstentsæt**  
Brugsanvisning
- DE**  
11 **Amplatz Ureterstent-Set**  
Gebrauchsanweisung
- EL**  
13 **Σετ ενδοπρόσθεσης ουρητήρα Amplatz**  
Οδηγίες χρήσης
- ES**  
15 **Equipo de stent ureteral Amplatz**  
Instrucciones de uso
- FR**  
17 **Set d'endoprothèse urétérale d'Amplatz**  
Mode d'emploi
- HU**  
20 **Amplatz ureterális sztentkészlet**  
Használati utasítás
- IT**  
22 **Set con stent ureterale Amplatz**  
Istruzioni per l'uso
- NL**  
24 **Amplatz ureterale stent set**  
Gebruiksaanwijzing
- NO**  
26 **Amplatz ureterstentsett**  
Bruksanvisning
- PL**  
28 **Zestaw stentu moczowodowego Amplatz**  
Instrukcja użycia
- PT**  
30 **Conjunto de stent ureteral Amplatz**  
Instruções de utilização
- SV**  
32 **Amplatz uretärstentset**  
Bruksanvisning


**Fig. 3**


**Fig. 4**

## AMPLATZ URETERAL STENT SET

**CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician or a properly licensed practitioner. Please read all instructions for use before using this device.**

### DEVICE DESCRIPTION

The Amplatz Ureteral Stent Set consists of a hydrophilic coated ultrathane stent with positioning string and radiopaque markers, and a stent positioner with stent introducer catheter and loading stylet.

### INTENDED USE

The Amplatz Ureteral Stent Set is delivered percutaneously and is intended to establish drainage from the ureteropelvic junction to the bladder and stenting of the ureter for all patients in whom it is desirable to place a drain which does not extend externally. The Amplatz Ureteral Stents are not intended to remain indwelling more than 180 days.

### CONTRAINDICATIONS

This device is contraindicated in the presence of conditions which create unacceptable risk during catheterization.

### WARNINGS

This product contains **NATURAL RUBBER LATEX**, which may cause allergic reactions.

### PRECAUTIONS

- The product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of ureteral stents should be employed.
- This product is intended for use in adult patients only.
- Do not over-advance ureteral stent onto stent positioner. Crimping of stent could occur.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.
- Where long-term use is indicated, it is recommended that indwelling time not exceed 180 days and that the physician evaluate the catheter before this time has expired.

### POTENTIAL ADVERSE EVENTS

The complications that may result from the use of these devices include, but are not limited to:

- Extravasation
- Urinary tract infection
- Perforation or laceration of the kidney, renal pelvis, ureter, and/or bladder
- Abrasions
- Sepsis
- Hemorrhage
- Bowel injury
- Iatrogenic renal artery injury
- Pleural complications (including but not limited to Pneumothorax)
- Hematuria
- Fever
- Encrustation
- Stent dislodgement
- Stent obstruction
- Lower urinary tract symptoms
- Pain

## MRI SAFETY INFORMATION

Nonclinical testing has demonstrated that the Amplatz Ureteral Stent is **MR Conditional**. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 T and 3.0 T
- Maximum spatial field gradient of 1,900 gauss/cm (19 T/m)
- Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of 2 W/kg (Normal Operating Mode)

Under the scan conditions defined above, the Amplatz Ureteral Stent is expected to produce a maximum temperature rise of less than 2.5° C after 15 minutes of continuous scanning.

In nonclinical testing, the image artifact caused by the Amplatz Ureteral Stent extends approximately 2 mm from the device when imaged with a gradient echo pulse sequence and a 3.0 T MRI system.

## INSTRUCTIONS FOR USE

1. Using standard percutaneous access technique, establish wire guide position through the renal pelvis and ureter into the bladder using a support wire guide in a diameter appropriately sized to the endhole of the catheter.
2. Determine appropriate catheter size using imaging technique.
3. Load the non-tapered proximal end of the ureteral stent onto the distal end of the pre-assembled stent introducer.
4. Holding the loading stylet beyond the tapered end of the ureteral stent introducer, lightly grasp the stent itself and begin advancing the non-tapered proximal end of the stent onto the pre-assembled stent introducer. Advance the stent until its proximal end aligns with and covers the radiopaque marker on the stent introducer catheter. (**Fig. 1**)
5. **Remove the loading stylet from the pre-assembled stent introducer. The loading stylet is not to be used for introduction of ureteral stent into the patient.**
6. Loosen the Tuohy-Borst adapter and position the distal end of the stent positioner so that it abuts the proximal end of the ureteral stent at the radiopaque marker. (**Fig. 2**) Retighten the Tuohy-Borst adapter.
7. Activate the hydrophilic coating by wetting the surface of the device with sterile water or saline. For best results, maintain wetted condition of the device during placement. The ureteral stent is now ready for introduction.
8. Introduce the loaded, pre-assembled stent introducer over the wire guide through the ureter into the bladder. Final position of the ureteral stent, prior to actual placement, should be verified by fluoroscopy or radiography.
9. To release the stent for final placement:
  - a. Loosen the Tuohy-Borst adapter.
  - b. Holding the sidearm fitting of the stent positioner to prevent inadvertent advancement of the stent, slowly pull back and remove the stent introducer catheter from the stent positioner. (**Fig. 3**) The ureteral stent is now positioned at its predetermined site.

**NOTE:** Do not advance the stent positioning catheter. Doing so will advance the ureteral stent.
10. Carefully withdraw the wire guide from the ureteral stent, but not completely from the collecting system. This will allow the ureteral stent to resume its double "J" shape. (**Fig. 4**)
11. Cut the string and remove. The wire guide can now be advanced further into the renal pelvis. **NOTE:** If the string cannot be easily removed, tie an additional length of string to the cut string, then advance an appropriately sized dilator over the string. While the dilator holds the stent in place, slowly remove the string.
12. Remove the stent positioner and place a temporary external drainage catheter over the wire guide into the collecting system. The wire guide can now be removed.

## HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

## REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2021
Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd