COOK**MEDICAL****CE**
0123

EN 2	Amplatz Wire Guides Instructions for Use
BG 4	Телени водачи Amplatz Инструкции за употреба
CS 6	Vodicí dráty Amplatz Návod k použití
DA 8	Amplatz-kateterledere Brugsanvisning
DE 10	Amplatz Führungsdrähte Gebrauchsanweisung
EL 12	Συρμάτινοι οδηγοί Amplatz Οδηγίες χρήσης
ES 14	Guías Amplatz Instrucciones de uso
FI 16	Amplatz-johtimet Käyttöohjeet
FR 18	Guides d'Amplatz Mode d'emploi
HR 20	Žice vodilice Amplatz Upute za uporabu
HU 22	Amplatz vezetődrótok Használati utasítás
IT 24	Guide Amplatz Istruzioni per l'uso
LT 26	„Amplatz“ vieliniai kreipikliai Naudojimo instrukcija
LV 28	Amplatz vadītājstīgas Lietošanas instrukcija
NL 30	Amplatz voerdraden Gebruiksaanwijzing
NO 32	Amplatz ledevaiere Bruksanvisning
PL 34	Prowadniki Amplatz Instrukcja użycia
PT 36	Fios guia Amplatz Instruções de utilização
RO 38	Fire de ghidaj Amplatz Instrucțiuni de utilizare
SL 40	Žična vodila Amplatz Navodila za uporabo
SR 43	Amplatz žice vodiči Uputstvo za upotrebu
SV 45	Amplatz-ledare Bruksanvisning
TR 47	Amplatz Kılavuz Teller Kullanma Talimatı

T _ A W G - M _ R E V O

AMPLATZ WIRE GUIDES

Read all instructions carefully. Failure to properly follow the information provided may lead to the device not performing as intended or injury to the patient.

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

All Amplatz Wire Guides are manufactured as fixed core wire guides, designed with a gradual transition to a straight flexible tip. These devices are provided in stiff, extra stiff, and ultra stiff configurations. These devices are manufactured in 0.025, 0.035 and 0.038 inch diameters and lengths of 145 cm. These devices are manufactured with an 5-9 cm flexible tip length or a 8 cm double flexible tip length. These devices are designed with a polytetrafluoroethylene (PTFE) coating on the outer surface of the coils.

Performance Characteristics

- The devices identified by a part number ending with “-DF” are double flexible wire guides, which allow for increased flexibility.
- Stainless steel coiling material for resistance to kinking.

Device Compatibility

These devices are compatible with appropriately sized needles, sheaths, dilators, endoscopes, stents, and catheters.

Patient Population

These devices are used for adult patients with a urological condition who require a procedure using a wire guide.

Intended Users

These devices are intended for use by physicians or healthcare practitioners with proper authorization to perform medical procedures (e.g., licensing, qualifications, credentials) per their local administrative and regulatory guidelines.

Contact with Body Tissue

These devices interact with the structures and tissues of the upper and lower urinary tract.

Operating Principle

These devices are inserted into the body through the lumen of a percutaneous entry needle, guiding sheath, or other interventional device. Only the flexible tip is intended as the initially inserted end. After insertion, medical imaging is used to visualize the movement of the wire guide to the desired location within the body. Catheters or other devices can then be inserted over the wire guide and directed to the intended location.

INTENDED USE

These devices are stainless steel fixed core wire guides that assist the insertion, removal, and positioning of other devices within the urinary tract. These devices act as a guide to access a portion of the urinary tract during a urological procedure to allow placement, replacement, or exchange of medical devices.

INDICATIONS FOR USE

Multiple medical conditions may result in the need for a urological procedure. Urology wire guides are used as a guide to access a portion of the urinary tract to allow placement, replacement, or exchange of medical devices during urological procedures such as percutaneous nephrostomy (PCN), percutaneous nephrolithotomy (PNL), ureteroscopy (URS), and stent placement or exchange.

DEVICE AND CLINICAL BENEFITS

To assist with urinary tract access for minimally invasive placement, replacement, or exchange of medical devices (without the wire guide itself having a direct diagnostic or therapeutic function).

CONTRAINDICATIONS

None

WARNINGS

- The single use device is not designed for re-use. Attempts to reprocess (re-sterilize) and/or to re-use may lead to chemical contamination, device failure, and/or transmission of disease.
- Do not use the device if the sterile packaging is damaged or unintentionally opened before use.
- Do not force the device during the procedure; proceed with caution if excessive resistance is encountered. Excessive force can lead to the wire guide fracturing, bending, kinking, unraveling or uncoiling. Additionally, excessive force can cause the wire guide to become entrapped.

- Imaging and/or placement procedures may expose the user and patient to radiation. The user should consider potential risks associated with radiation and appropriately plan procedures to minimize exposure in patients and medical personnel. The occurrence and/or severity of undesirable side effects may be dose dependent or probabilistic.

PRECAUTIONS

-  This symbol on the label indicates that the device contains cobalt (Co) at a level above 0.1%w/w, which is a substance that is toxic to reproduction and a carcinogen (Class IB). However, the device contains cobalt as a stainless-steel alloy containing cobalt, which does not cause an increased risk of cancer or adverse reproductive effect according to current scientific evidence.
- When exchanging or withdrawing an instrument over the wire guide, secure and maintain the wire guide in place under fluoroscopy in order to avoid unexpected wire guide displacement or migration.
- Manipulation of wire guide requires appropriate imaging control. Use caution not to force or over manipulate the wire guide when gaining access. Over manipulation can cause the wire guide to become entrapped.
- When using wire guide through a metal cannula/needle, use caution as damage may occur to outer coating, such as the coating shearing off the wire guide.

POTENTIAL ADVERSE EVENTS

- Acute Bleeding / hematoma
- Acute pyelonephritis
- Injury to adjacent organs
- Hematuria
- Infection / UTI / urosepsis
- Mucosal tear / laceration / abrasion
- Pain
- Ureteral laceration / avulsion / tear / lesion
- Ureteral perforation / false passage / urinoma
- Ureteral stenosis
- Ureteral stricture
- Urinary tract injury
- Vascular injury

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Keep dry and away from sunlight. Upon removal from package, inspect the product to ensure no damage has occurred.

INSPECTION OF DEVICE

Visually inspect the device thoroughly including all levels of the packaging (as applicable) to verify that there is no damage prior to use. Visually inspect and confirm that the integrity of the sterile barrier has not been compromised in any way.

INSTRUCTIONS FOR USE

1. Using aseptic technique, remove the wire from its outer packaging and place in the sterile field.
2. Remove the wire guide from the holder.
3. Advance the flexible end of the wire guide into the device and position to the desired anatomical location.

NOTE: Endhole size and length of the device must be taken into consideration to ensure proper fit between wire guide and device.

4. After wire guide use is complete, remove from patient and discard.

DISPOSAL OF DEVICES

After the procedure, this device may be contaminated with potentially infectious substances of human origin and should be disposed of in accordance with institutional guidelines.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

PATIENT COUNSELING INFORMATION

Please inform the patient as necessary of the relevant warnings, precautions, contraindications, measures to be taken, and limitations of use.

SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device, this should be reported to Cook Medical and this shall also be reported to the competent authority where the device was used.

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Cook Medical Europe Ltd
William A. Cook Australia Pty Ltd