EN
3**Bander Ureteral Diversion Stent Set**

Instructions for Use

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Bruksanvisning

T _ B A N D S 2 _ R E V 2

BANDER URETERAL DIVERSION STENT SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

- Radiopaque Bander stent
- .038 inch (0.97 mm) diameter wire guide
- Catheter retainer
- Catheter adapter

NOTE: Set components may vary.

INTENDED USE

Used for intraoperative placement to stent the ureter and provide drainage from the kidney to the external stoma during ureteroileal conduit construction and continent urinary diversions.

CONTRAINDICATIONS

None known

WARNINGS

None known

PRECAUTIONS

- Manipulation of the wire guide requires appropriate imaging control. Use caution not to force or over manipulate the wire guide when gaining access.
- When using a wire guide through a metal cannula/needle, use caution as damage may occur to outer coating.
- When exchanging or withdrawing an instrument over the wire guide, secure and maintain the wire guide in place under fluoroscopy in order to avoid unexpected wire guide displacement.
- Hydrophilically coated wires are very slippery when wet. Always maintain control of the wire guide when manipulating it through any device.
- The maximum indwelling time for this ureteral stent is 90 days. All stents may be subject to varying degrees of encrustation when placed in the urinary tract. Periodic checks of the stent by cystoscopic and/or radiographic means are recommended.

MRI SAFETY INFORMATION

Nonclinical testing has demonstrated the Bander Ureteral Diversion Stent is **MR Conditional**. A patient with this device can be safely scanned in an MR system meeting the following conditions.

- Static magnetic field of 3.0 T
- Maximum spatial field gradient of 19 T/m

The Bander Ureteral Diversion Stent has been deemed to be MR Conditional under these scan conditions where the nonclinical tests of magnetically induced displacement angle/force and magnetically induced torque result in no risk to the patient. Furthermore, the results of the electrical conductivity test indicate that the device will not experience RF induced heating above the background heating in the MR environment.

In nonclinical testing, the image artifact caused by the device extends approximately 3.6 mm from the Bander Ureteral Diversion Stent when imaged with a gradient echo pulse sequence and a 3.0 T MRI system.

INSTRUCTIONS FOR USE

Wire Guide Preparation

1. Using aseptic technique, remove the wire from its outer packaging and place in the sterile field.
2. If using a hydrophilically coated wire guide, fill a syringe with sterile water or sterile saline solution and attach it to the flushing port on the wire guide holder. Inject enough solution to wet the wire guide surface entirely. This will activate the hydrophilic coating. **NOTE:** For optimal performance, rehydrate the hydrophilically coated wire guide after exposure to ambient environment or, after extended use, replace it with a new compatible hydrophilically coated wire guide.
3. Remove the wire guide from the holder.
4. If using a movable core wire guide, adjust mandril position by advancing or withdrawing the mandril within the wire guide to achieve the required degree of flexibility.

Bander Ureteral Diversion Stent

1. Pass the wire guide, flexible end first, through the open end of the ureteral stent to straighten the pigtail coil.
2. At the anastomosis site, place the ureteral stent, coil end first, through the ureter to the renal pelvis.
3. Place the opposite end of the stent through the ileal conduit and out through the stoma.
4. Remove the wire guide. The stent pigtail will form spontaneously.
5. Check to assure the stent has formed a full coil in the renal pelvis.
6. Slide the catheter retainer onto the externalized portion of the stent and tighten the retainer on the stent.
7. Suture the retainer to the skin near the stoma to provide added assurance against dislodgment during and following the procedure.
8. If necessary, use a connecting tube (not included in the set) to connect the stent to a drainage bag.
9. The correct position of the stent may be confirmed by radiographic means.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2020
Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd