EN
3

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Bruksanvisning

T _ B A N D S _ R E V 1

BANDER URETERAL DIVERSION STENT SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

- Single pigtail ureteral stent
- .038 inch (0.97 mm) diameter stainless steel stylet with 5 mm flexible tip
- Two catheter retainers
- Wire guide (*size and configuration may vary*)

NOTE: Set components may vary.

INTENDED USE

Sof-flex® and silicone bands are used for intraoperative placement to stent the ureter during ureteroileal conduit construction and continent urinary diversions.

Wire guides are used to gain ureteral access, to establish a tract, and to assist in the placement, replacement, and exchange of medical devices during urological and gynecologic procedures. The HiWire, BiWire and Roadrunner Wire Guides may also be used for catheter positioning and exchange in a tortuous or kinked ureter, traversing a large stone in route to the kidney, or in cases demanding enhanced control and high radiopacity.

CONTRAINDICATIONS

None known

WARNINGS

None known

PRECAUTIONS

- Manipulation of the wire guide requires appropriate imaging control. Use caution not to force or over manipulate the wire guide when gaining access.
- When using a wire guide through a metal cannula/needle, use caution as damage may occur to outer coating.
- When exchanging or withdrawing an instrument over the wire guide, secure and maintain the wire guide in place under fluoroscopy in order to avoid unexpected wire guide displacement.
- The included wire guide is not intended for PTCA use.
- Hydrophilically coated wires are very slippery when wet. Always maintain control of the wire guide when manipulating it through any device.
- The maximum indwelling time for this ureteral stent is six months.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

INSTRUCTIONS FOR USE

Wire Guide Preparation

1. Using aseptic technique, remove the wire from its outer packaging and place in the sterile field.
2. If using a hydrophilically coated wire guide, fill a syringe with sterile water or sterile saline solution and attach it to the flushing port on the wire guide holder. Inject enough solution to wet the wire guide surface entirely. This will activate the hydrophilic coating. **NOTE:** For optimal performance, rehydrate the hydrophilically coated wire guide after exposure to ambient environment or, after extended use, replace it with a new hydrophilically coated wire guide.

NOTE: For the Roadrunner® PC Wire Guide Double Flexible, the pink end of the wire guide is not hydrophilically coated or intended for insertion into the body.

3. Remove the wire guide from the holder.
4. If using a movable core wire guide, adjust mandril position by advancing or withdrawing the mandril within the wire guide to achieve the required degree of flexibility.

Bander Ureteral Diversion Stent

1. Pass the wire guide, flexible end first, through the open end of the ureteral stent to straighten the pigtail coil.
2. At the anastomosis site, place the ureteral stent, coil end first, through the ureter to the renal pelvis.
3. Place the opposite end of the stent through the ileal conduit and out through the stoma.
4. Remove the wire guide. The stent pigtail will form spontaneously.
5. Check to assure the stent has formed a full coil in the renal pelvis.
6. Slide one or both catheter retainers onto the externalized portion of the stent and tighten the retainers on the stent.
7. Suture the retainers to the skin near the stoma to provide added assurance against dislodgment during and following the procedure.
8. If necessary, use a connecting tube to connect the stent to a drainage bag.
9. The correct position of the stent may be confirmed by radiographic means.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2012
Bloomington, IN 47404 U.S.A.
COOK IRELAND LTD.