COOK

MEDICAL

EN
2

Advance® Biliary Balloon Catheter

Instructions for Use

ES
4

Catéter biliar con balón Advance®

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FR
6

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Advance®

Mode d'emploi

T _ B B D C _ R E V O

ADVANCE® BILIARY BALLOON CATHETER

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

The Advance Biliary Balloon Catheter is a double-lumen catheter with a balloon near its distal tip. The catheter consists of two independent lumens, which are labeled "DISTAL" and "BALLOON." The distal lumen extends the length of the catheter and is used for placement of wire guides. The balloon lumen is used to expand the balloon. Inscribed on the tip of the manifold are the balloon diameter (mm) and balloon length (cm).

The balloon is manufactured from an extra-thin wall, high-strength, minimally compliant material. Particular care should be taken in handling the balloon to prevent damage. It will inflate to the indicated size parameters when utilizing proper pressure recommendations. Refer to label for further information. Use of a pressure gauge is recommended to monitor inflation pressures.

Radiopaque markers are positioned on the shaft within the balloon to enable visualization of the catheter/balloon under fluoroscopy. The balloon catheter is compatible with 0.035 inch (0.89mm) wire guides.

INTENDED USE

The Advance Biliary Balloon Catheter is indicated for laparoscopic and general surgical procedures for the dilation of the cystic duct to facilitate common bile duct exploration.

CONTRAINDICATIONS

- Peritonitis
- Severe cholangitis
- Severe acute pancreatitis
- Previous abdominal operations which prevent safe abdominal access or progression of the procedure
- Inability to advance the balloon into the cystic duct
- Coagulopathy
- Known or suspected perforation of the cystic and/or bile ducts

WARNINGS

- This product is not intended for use in the vasculature.
- Do not exceed the recommended balloon inflation pressure as listed on the product label.
- Do not advance the balloon if resistance is encountered. Stop and assess the cause of the resistance to determine if dilation should be re-attempted.
- Do not use a power injector for balloon inflation or injection of contrast medium through catheter lumen marked "DISTAL." Rupture may occur.
- If package is opened or damaged when received, do not use. Visually inspect with particular attention to kinks, bends and breaks. If an abnormality is detected that would prohibit proper working condition, do not use. Please notify Cook for return authorization.
- Carefully read all instructions prior to use. Failure to observe all warnings and precautions may result in complications.

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for biliary access and procedures should be employed.
- The balloon catheter is constructed of heat sensitive material. Do not heat or attempt to shape the catheter tip.
- The balloon catheter should be inserted and manipulated under direct visualization and/or fluoroscopic guidance.
- The balloon catheter is used in conjunction with an inflation device or manometer and normal saline or a 1:1 mixture of contrast medium and normal saline. **DO NOT** use air or gaseous substances to inflate the balloon catheter, as this will result in reduced balloon effectiveness.
- This product is not recommended to inflate the balloon more than once during the procedure.

POTENTIAL ADVERSE EVENTS

- Perforation of cystic or bile ducts
- Cholangitis
- Hemorrhage
- Infection
- Allergic reaction to contrast medium
- Hypotension

PRODUCT RECOMMENDATIONS

**PTA Biliary Balloon Catheter
Compliance Table
Diameter vs. ATM**

		ATM										
Diameter (mm)	4	5	6	7	8	9	10	11	12	13	14	15
6.0	5.26	5.39	5.55	5.68	5.80	5.90	6.00	6.09	6.17	6.26	6.34	6.44
8.0	7.17	7.42	7.66	7.85	8.00	8.15	8.27	8.41				

Nominal Pressure and corresponding balloon diameter

Burst Pressure and corresponding balloon diameter

**ATB Biliary Balloon Catheter
Compliance Table
Diameter vs. ATM**

		ATM										
Diameter (mm)	4	5	6	7	8	9	10	11	12	13	14	15
8.0	7.79	8.00	8.11	8.25	8.38	8.49	8.58	8.66	8.74	8.81	8.89	8.97

Nominal Pressure and corresponding balloon diameter

Burst Pressure and corresponding balloon diameter

INSTRUCTIONS FOR USE

1. Remove the balloon catheter from the package using standard sterile technique.
2. Flush the balloon catheter lumen labeled "DISTAL" with normal saline solution.
3. Apply negative pressure to lumen labeled "BALLOON" prior to introduction.
4. Advance the balloon catheter over the prepositioned 0.035 inch (0.89mm) wire guide while visually monitoring. Continue advancing the balloon catheter into the cystic duct until the desired position is achieved.
NOTE: If resistance is met while advancing the balloon catheter, determine the cause and proceed with caution.
5. Inflate balloon to specified diameter without exceeding the rated burst pressure, as listed on the product label, using normal saline or a 1:1 mixture of contrast medium and normal saline.
NOTE: Position may be verified using the fluoroscopic markers located at either end of the balloon.
6. To deflate the balloon pull vacuum on the inflation syringe or inflation device.
7. Maintain vacuum on the balloon and observe for complete deflation. Carefully withdraw the balloon catheter after the balloon is completely deflated.
8. Continue common bile duct exploration procedures as standard.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.