COOK

MEDICAL

CE
0088

EN
2

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Bruksanvisning

T _ B B F S _ R E V 1

TRANSLUMINAL BILIARY BIOPSY FORCEPS SET

DESCRIPTION

The Transluminal Biliary Biopsy Forceps Sets consist of:

The BBFS-100 and BBFS-100-MS contain:

- A 5.2 French flexible biopsy forceps
- A 7.0 French Flexor sheath and dilator

The BBFS-100 also contains:

- A percutaneous access set (containing a 22 gage access needle and stylet, an 0.018 inch mandril wire guide and a dilator sheath assembly)
- An 0.035 inch fixed-core wire guide
- An 0.035 inch angled-tip hydrophilic-coated wire guide
- A 5.0 French hydrophilic-coated catheter

NOTE: Use of the BBFS-100-MS requires use of wire guides and catheters from institutional inventory. Once wire guide position is obtained, continue with the procedure using the Flexor sheath and forceps from the BBFS-100-MS set

INTENDED USE

The Transluminal Biliary Biopsy Forceps Set is intended for access to and biopsy of tissue within the biliary ductal system.

CONTRAINDICATIONS

None known

WARNINGS

None known

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in diagnostic and Interventional techniques.
- Imaging guidance is recommended when using this product.
- The mandril wire guide in the percutaneous access set should not be withdrawn through the needle once placed beyond the needle tip.
- Prior to use of the 0.035 hydrophilic-coated wire guide, hydrate the wire guide by injecting sterile heparinized saline solution into the Luer lock hub end of its dispenser.
- Prior to use of the hydrophilic-coated catheter, hydrate the catheter by wetting the distal end of the catheter with sterile water or saline.
- If resistance is encountered during wire guide or catheter manipulation, stop and determine the cause before proceeding any further.
- Manipulation of the biopsy forceps requires fluoroscopic visualization.
- Visually inspect the biopsy forceps before use to make sure the shaft is free of kinks or severe angulation.
- Verify the distal position of the sheath prior to each biopsy.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

INSTRUCTIONS FOR USE

Prep and drape the patient according to hospital protocols.

Initial Percutaneous Access

NOTE: If access to the biliary ducts is already present, steps 1-6 may be omitted.

1. Place the 22 gage needle from the access set percutaneously from skin puncture to the desired location within the biliary ductal system.
2. Once the needle tip position is confirmed, remove the stylet from the percutaneous puncture needle.
3. Advance the 0.018 inch mandril wire guide through the needle cannula into the biliary duct, confirming position with imaging guidance.
4. Maintaining position of the wire guide, remove the needle cannula.
5. Advance the dilator sheath assembly over the wire guide until the distal tip is within the biliary ductal system.
6. Maintaining position of the outer sheath, remove the wire guide, stiffening cannula and dilator. **NOTE:** If desired, the wire guide can remain in position as a safety wire.

Percutaneous Transluminal Biliary Biopsy Procedure

7. Insert the 0.035 inch angled-tip, hydrophilic-coated wire guide through the percutaneous access sheath or biliary drainage catheter.
8. Maintaining position of the wire guide in the biliary duct, remove the access sheath or catheter.
9. Insert the 5.0 French angiographic catheter over the wire guide and position it within the biliary ductal system.
10. Advance the wire guide and catheter as a unit, positioning the wire guide and catheter tips beyond the area to be biopsied.

11. Exchange the 0.035 inch hydrophilic-coated wire guide for the 0.035 inch fixed-core wire guide, keeping the catheter in position beyond the target biopsy anatomy.
12. Maintaining the wire guide in position, remove the catheter.
13. Insert the Flexor sheath assembly over the wire guide, positioning the sheath tip at the ductal area to be biopsied.
14. Remove the dilator from the Flexor sheath assembly, leaving the wire guide in position within the sheath and across the target biopsy anatomy.
15. Insert the biopsy forceps through the sheath, positioning the biopsy cup at the distal sheath tip.
16. Maintaining forceps position, withdraw the Flexor sheath to uncover the biopsy cup mechanism.
17. Position the sheath and forceps assembly at the intended biopsy site.
18. Move the thumb button on the handle of the forceps forward to open the biopsy cup.
19. Maintaining the biopsy cup in an open position, advance the sheath and forceps assembly as a unit to the area where the specimen for biopsy is to be obtained.
20. While maintaining forward pressure on the sheath and forceps assembly to keep the biopsy sample against the duct wall, retract the thumb button to close the jaws of the cup and secure the tissue specimen.
21. While maintaining the thumb button in the retracted position to ensure the cup remains closed and the specimen is kept intact, withdraw the biopsy forceps from the Flexor sheath and retrieve the specimen from the biopsy cup. **NOTE:** If more than one biopsy specimen is required, rinse the forceps with saline and repeat steps 13 – 20.
22. Once biopsy sampling is complete, the Flexor sheath may be exchanged over the wire guide for a biliary drainage catheter.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2014
Bloomington, IN 47404 U.S.A.