- EN**  
3  
**Brush Biopsy Set**  
Instructions For Use
- CS**  
5  
**Souprava kartáčku pro biopsii**  
Návod k použití
- DA**  
7  
**Biopsisæt med børste**  
Brugsanvisning
- DE**  
9  
**Bürstenbiopsie-Set**  
Gebrauchsanweisung
- EL**  
11  
**Σετ ψήκτρας βιοψίας**  
Οδηγίες χρήσης
- ES**  
13  
**Equipo de biopsia por cepillado**  
Instrucciones de uso
- FR**  
15  
**Set de biopsie à brosse**  
Mode d'emploi
- HU**  
17  
**Kefebiopsziás készlet**  
Használati utasítás
- IT**  
19  
**Set per biopsia con spazzolino**  
Istruzioni per l'uso
- NL**  
21  
**Borstel biopsieset**  
Gebruiksaanwijzing
- NO**  
23  
**Børstebiopsisett**  
Bruksanvisning
- PL**  
25  
**Zestaw biopsyjny ze szczoteczką**  
Instrukcja użycia
- PT**  
27  
**Conjunto de escova de biopsia**  
Instruções de utilização
- SV**  
29  
**Borstbiopsiset**  
Bruksanvisning

T \_ B B S \_ R E V 1

## BRUSH BIOPSY SET

**CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).**

### DEVICE DESCRIPTION

Brush Biopsy Set product configurations:

- 5 French and 67 cm catheter with wire brush only
- 3.2 French and 115 cm catheter with wire brush, inserter, and endoscopic cap

Deflectable Brush Biopsy Set product configurations:

- 5 French and 67 cm catheter with deflectable wire brush only
- 3.2 French and 115 cm catheter with deflectable wire brush, inserter, and endoscopic cap

### INTENDED USE

The Brush Biopsy Sets are intended to obtain pathology specimens from lesions in the ureter, renal pelvis, infundibula or calyces under direct vision.

### CONTRAINDICATIONS

None known

### WARNINGS

None known

### PRECAUTIONS

- This product is intended for use by physicians trained and experienced in endoscopic biopsy procedures. Standard techniques should be employed.
- The potential effects of phthalates on children or pregnant or nursing women have not been fully characterized and there may be concern for reproductive and developmental effects.

### POTENTIAL ADVERSE EVENTS

Complications associated with endoscopic biopsy procedures may include, but are not limited to:

- Bleeding
- Urinary tract perforation
- Stricture
- Urinary retention

Additional complications associated with the Brush Biopsy Set may include, but are not limited to:

- Inadvertent puncture with tissue damage or bleeding
- Infection
- Injury to the ureter or kidney

### INSTRUCTIONS FOR USE

1. If using the Deflectable Brush Biopsy Set, confirm deflection prior to use. Unsheath the brush and use the handle to deflect. If brush does not deflect, slide the two-rings to the distal end of handle assembly, and tighten set screw to secure.
2. If applicable, place the endoscopic cap on the scope port and place the inserter securely into the endoscopic cap. With the inserter in place, carefully insert the device into the scope.
3. Make sure to sheath the catheter over the brush before insertion. If necessary, loosen the catheter from the proximal assembly.
4. Under direct endoscopic vision, advance the device through the appropriate endoscopic working channel.

5. Position the endoscope near the intended biopsy site. **NOTE:** The biopsy sets should always be used under an operator's direct vision to reduce risks of knotting, perforation, and difficulty in device removal.
6. Advance the distal tip of the brush to the area of interest.
7. With the brush in position, deploy the brush and scrape tissue with multiple passes to obtain specimen. If using the Deflectable Brush Biopsy Set and if deflection is needed to reach the area, use the handle to deflect the brush.
8. Sheath the catheter over the brush prior to removing to preserve specimen. If necessary, loosen the catheter from the proximal assembly.
9. Remove brush from the endoscope.

#### **HOW SUPPLIED**

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

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Bloomington, IN 47404 U.S.A.