EN
4**Biliary Drainage Catheter**

Instructions for Use

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
Fig. 4

Biliary Drainage Catheter

Read all instructions carefully. Failure to properly follow the information provided may lead to the device not performing as intended or injury to the patient.

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

Biliary Drainage Catheters Ultrathane®, Ring Biliary Duct Drainage Catheter Polyethylene and Ring Biliary Duct Drainage Catheter Ultrathane® come in a range of French sizes, lengths and sideport quantities. The catheter may be supplied with a flexible catheter introduction stiffening cannula, a rigid stiffener, and a catheter securement device. Locking Loop Catheters are available with one of two types of locking mechanisms:

- Mac-Loc® Locking Loop (**Fig. 5**)
- Cook-Cope Type Loop (**Fig. 6**)

Performance Characteristics

These devices are manufactured in various French sizes and lengths; and 0.038 in end hole diameter, which allows the catheter to accept wire guides 0.038 in or smaller. Sideports allow fluid drainage facilitated by gravity. The tapered distal tip of the catheter facilitates over-the-wire trackability and is designed with a pigtail shape in order to maintain placement in the body. Products with part numbers including RH have a radiopaque band and products with part numbers including P8.3G have radiopaque tubing which aids in image guidance.

Device Compatibility

These devices are designed to be compatible with components such as stiffeners, fixation device, connecting tubes and drainage bags with standard Luer lock fittings, as well as standard size wire guides (0.038 inches or smaller).

Qualitative and Quantitative Information

Device Material		Weight (grams)
Polyurethane	Shaft Tubing (ULT devices only)	Up to 7.082
Polyethylene	Shaft Tubing (P8.3G devices only)	Up to 1.972
Hydrogel	Hydrophilic Coating (-HC and -RH devices only)	Up to 0.065
Silicone Lubricant	Silicone Coating (not included on -HC and -RH devices only)	Up to 0.010

Mac-Loc® Locking Loop Mechanism

Nylon	Suture	Up to 0.083
-------	--------	-------------

Cope-Loop Type Locking Loop Mechanism*

PTFE-Coated Polyester	Suture	Up to 0.068
-----------------------	--------	-------------

*The device includes a sleeve made of natural rubber latex. The sleeve is a skin-contacting component. **NATURAL RUBBER LATEX** may cause allergic reactions.

Patient Population

The target population of the device are adult patients who require placement of a catheter to allow for percutaneous biliary drainage according to the device's intended use.

Intended User

The Percutaneous Biliary Drainage Catheters and Sets are intended for biliary drainage procedures and should only be performed by licensed and trained clinical practitioners, such as physicians, and licensed and trained healthcare providers. Maintenance and care of the device may be performed by laypersons (e.g. patients, lay caregivers, or lay first responders), nurses, and technicians.

Contact with Body Tissue

These devices are externally communicating devices contacting tissue. They contact bile ducts, liver, duodenum, gallbladder, pancreatic ducts, and skin during use.

Operating Principle

These devices are placed with imaging guidance to identify and facilitate control of the distal tip of the catheter. Drainage occurs through the sideports and is facilitated by gravity. Placement of the catheter in the body is maintained through retention of the pigtail loop.

Expected Lifetime of the Device

The indwell time of the catheter is up to a maximum of 90 days.

INTENDED USE

The devices are catheters intended for percutaneous transhepatic drainage of bile by creating a conduit for internal or external bile drainage.

DEVICE AND CLINICAL BENEFITS

Provide a percutaneously placed pathway for biliary drainage.

INDICATIONS FOR USE

The most common indications for biliary drainage are to manage sequelae of biliary obstructions (malignant or benign tumors, internal or external to duct), biliary injury (leakage or stricture), biliary stones, or inflammatory processes such as pancreatitis.

CONTRAINDICATIONS

Uncorrectable coagulopathy

WARNINGS

- If a catheter has become malpositioned or if drainage ceases, the catheter should be promptly exchanged or removed.
- Minimize catheter manipulation during application and removal of the fixation device.
- See product label for latex indication. **NATURAL RUBBER LATEX** may cause allergic reactions.
- The single use device is not designed for re-use. Attempts to reprocess (re-sterilize) and or to re-use may lead to chemical contamination, device failure, and/or transmission of disease.
- Do not use if the sterile packaging is damaged or unintentionally opened before use.
- Do not use fixation device where loss of adherence could occur, such as with a confused patient or non-adherent skin, or if a patient has known tape or adhesive allergies.
- The Peel-Away® Pigtail Straightener, if present, is not to be used as a vascular introducer sheath.

PRECAUTIONS

-  This symbol on the label indicates that the device contains cobalt (Co) at a level above 0.1%w/w, which is a substance that is toxic to reproduction and a carcinogen (Class IB). However, the device contains cobalt as a stainless-steel alloy containing cobalt, which does not cause an increased risk of cancer or adverse reproductive effect according to current scientific evidence.
- Standard techniques for placement of percutaneous drainage catheters should be employed.
- Manipulation of products requires ultrasound, fluoroscopy, or other imaging guidance.
- When inserting a stiffening cannula into a catheter with retention suture, hold suture during cannula insertion to avoid bunching or tangling of suture.
- A TFE-coated wire guide must be used with Ultrathane catheters.
- Activate the hydrophilic coating, if present, by wetting the catheter with sterile water or saline. For best results, keep catheter surface wet during placement.
- Catheters should be irrigated on a routine basis to ensure function.
- Patients with indwelling drainage catheters should be evaluated routinely to ensure continuous function of the catheter and adherence of the fixation device, if applicable.
- Traction on the locking suture, if present, should be sufficient to ensure adequate retention of the tip, but should not be overly tight. Verify catheter tip configuration by fluoroscopy.
- It is recommended to use a wire guide when removing a Locking Loop catheter.
- Avoid fixation device contact with alcohol or acetone. Both can weaken bonding of components and the fixation device anchor pad adherence to the skin.

MRI SAFETY INFORMATION

Nonclinical testing has demonstrated that the Ultrathane catheter (ULT) is **MR Conditional** according to ISO 15223-1/2 ASTM F2503.

A patient with this device may be safely scanned after placement under the following conditions:

- Static magnetic field of 1.5 T or 3.0 T only
- Maximum magnetic field spatial gradient of 1,900 gauss/cm (19 T/m)
- Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of ≤ 2.0 W/kg (Normal Operating Mode) for 15 minutes of continuous scanning

Under the scan conditions provided above, this device is expected to produce a maximum temperature rise of less than 2.1°C after 15 minutes of continuous scanning in nonclinical testing.

The image artifact extends approximately 40.1 mm from the device as found during nonclinical testing when imaged with a gradient echo pulse sequence and a 3.0 T MRI system.

The Polyethylene catheter (P8.3G) is **MR Safe**.

POTENTIAL ADVERSE EVENTS

- Bile leak/biloma
- Bleeding
- Catheter-related events (e.g., blockage, dislodgement, fracture, pericatheter leakage, kinking)
- Death
- Fistula
- Fluid and electrolyte loss (associated with external drainage)
- Infection (e.g., abscess, cholangitis, fever)
- Inflammation (e.g., pancreatitis, cholecystitis)
- Injury/perforation of an adjacent organ or vessel (e.g., viscera perforation, pseudoaneurysm)
- Pain
- Tumor seeding

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Keep dry and away from sunlight. Upon removal from package, inspect the product to ensure no damage has occurred.

INSPECTION OF DEVICE

Visually inspect the device thoroughly including all levels of the packaging (as applicable) to verify that there is no damage prior to use. Visually inspect and confirm that the integrity of the sterile barrier has not been compromised in any way.

INSTRUCTIONS FOR USE

Catheter Placement

1. Under image guidance, perform standard techniques for placement of percutaneous drainage catheters.
NOTE: The provided rigid and flexible stiffeners are intended to be used in conjunction with the catheter for additional stability during catheter placement.
2. Once catheter is in desired location, remove any wire guides or stiffeners, allowing the catheter to form its configuration.

For Mac-Loc® Locking Loop Mechanism:

3. Lock the catheter in place using the Mac-Loc Locking Loop Mechanism, as described below.
 - a. Stabilize the Mac-Loc catheter hub assembly and pull back on the drawstring to form the distal catheter loop configuration. **(Fig. 1)**
 - b. While maintaining traction on the drawstring, push the locking cam lever down until a distinct 'snap' is felt. **(Fig. 2)**
 - c. Trim the excess drawstring. **(Fig. 3)**

For Cook-Cope Type Locking Loop Mechanism:

- a. Pull suture end tight to form loop configuration in catheter, and tie securely.
- b. Trim off excess suture, and slide latex sleeve over suture to prevent leakage.

Unlocking Catheter Loop

For Mac-Loc® Locking Loop Mechanism:

- a. While stabilizing the Mac-Loc catheter hub assembly with one hand, position a small, blunt object (approximately the shape and size of a ball point pen or small forceps) into the Mac-Loc release notch.
- b. Pry upward until the locking cam lever is free. **(Fig. 4)**
NOTE: For catheter exchange, advance the distal end of a wire guide into the locked loop configuration of the catheter before unlocking the Mac-Loc assembly. Release the Mac-Loc as described above. Advance the wire guide through the catheter end hole. Catheter exchange can now be performed.

For Cook-Cope Type Locking Loop Mechanism:

- a. Advance wire guide into catheter.
- b. Uncover suture by folding back the latex sleeve, and release suture.
- c. Withdraw catheter over wire guide.
NOTE: For catheter exchange, insert new catheter over wire guide. For removal, withdraw wire guide following removal of catheter.

Applying the Catheter Fixation Device

1. Bring the fixation device anchor pad into position. Stabilize catheter between fingers. Place catheter into fixation device retainer.
2. Close lid, peel off paper backing, and place on skin.
3. Suture the fixation device to the skin if desired or deemed necessary.

Removing the Catheter Fixation Device

1. To open retainer, press tab gently with thumb and lift lid.
2. Remove the fixation device anchor pad adhesive with alcohol.
NOTE: The fixation device should be replaced at least every 7 days.

DISPOSAL OF DEVICES

This device may be contaminated with potentially infectious substances of human origin and should be disposed of in accordance with institutional guidelines.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

PATIENT COUNSELING INFORMATION

Please inform the patient as necessary of the relevant warnings, precautions, contraindications, measures to be taken, and limitations of use that the patient should be aware of. You can use the basic unique device identifier for this device seen in **Table 1** below to find the Summary of Safety and Clinical Performance (SSCP) on the EUDAMED website. When the EUDAMED is available, use the following link: ec.europa.eu/tools/eudamed

Part Number Prefix	Part Number Suffix	Basic Unique Device Identifier (BUDI)
ULT...	...BCL, ...CLB or ...CLBS	0827002CINC050-051_A-1JP
ULT...	...MEH	0827002CINC050-051_A-4JV
ULT...	...RING	0827002CINC050-051_A-2JR
P8.3G...	...RING	0827002CINC050-051_A-3JT

SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device, this should be reported to Cook Medical and this shall also be reported to the competent authority where the device was used.

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© COOK 2025
Cook Medical Europe Ltd.
William A. Cook Australia Pty. Ltd.