EN 3 Black Silicone Filiform Double Pigtail Ureteral Stent

Instructions for Use

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
Fig 1


Fig 2

BLACK SILICONE FILIFORM DOUBLE PIGTAIL URETERAL STENT

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner). Please read all instructions before using this device.

DEVICE DESCRIPTION

- Radiopaque silicone double pigtail stent
- Stent positioner with lockable fitting
- T inserter
- Dimensions include 6 French, 7 French, and 8.5 French with fixed lengths from 20 cm to 30 cm

NOTE: 6, 7, and 8.5 Fr stents are compatible with an 0.038" wire guide.

INTENDED USE

This product is used for temporary internal drainage from the ureteropelvic junction to the bladder. The device is not intended to remain indwelling more than twelve months.

CONTRAINDICATIONS

This device is contraindicated in the presence of conditions which create unacceptable risk during catheterization.

WARNINGS

Sterilized by Ethylene Oxide gas. Do not use if sterile barrier is damaged. The device is intended for one-time use.

PRECAUTIONS

- Complications of ureteral stent placement are documented. Use of this device should be based upon consideration of risk-benefit factors as they apply to your patient. Informed consent should be obtained to maximize patient compliance with follow-up procedures.
- Do not force set components during placement, replacement, or removal. Carefully remove the set components if any resistance is encountered.
- The tether should be removed if the stent is to remain indwelling longer than 14 days.
- The stent must not remain indwelling more than twelve months. If the patient's status permits, the stent may be replaced with a new stent.
- The included stent is not intended as a permanent indwelling device.
- A pregnant patient must be more closely monitored for possible stent encrustation due to calcium supplements.
- Improper handling can seriously weaken the stent. Acute bending or overstressing during placement may result in subsequent separation of the stent at the point of stress after a prolonged indwelling period. Angulation of the stent should be avoided.
- Individual variations of interaction between stents and the urinary system are unpredictable.
- Periodic evaluation via cystoscopic, radiographic, or ultrasonic means is suggested. The stent must be replaced if encrustation hampers drainage.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.
- If problems occur using this device, please call your Cook Urological sales representative or contact our Customer Service department at the address/phone listed at www.cookmedical.com.

POTENTIAL ADVERSE EVENTS

- Migration and dislodgement
- Pain and discomfort
- Urinary frequency and urgency
- Perforation and fistula formation

- Stent obstruction by stone or tissue
- Stent fragmentation

MRI SAFETY INFORMATION

Nonclinical testing has demonstrated the Black Silicone Filiform Double Pigtail Ureteral Stent is **MR conditional**. A patient with these devices can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 3.0 tesla or 1.5 tesla only
- Maximum magnetic field spatial gradient of 1900 gauss/cm (19.0 T/m) or less

Nonclinical testing demonstrated that the Black Silicone Filiform Double Pigtail Ureteral Stent is not expected to induce RF heating greater than that of background heating.

In nonclinical testing, the image artifact caused by the device extends approximately 2.4 mm from the Black Silicone Filiform Double Pigtail Ureteral Stent when imaged with a gradient echo pulse sequence and a 3.0 tesla MR system.

INSTRUCTIONS FOR USE

One-Pass Placement

1. Introduce the wire guide into the fitting end of the stent positioner and advance it until the wire guide extends half its length beyond the blunt end of the positioner.
NOTE: The T inserter can be utilized to help introduce the wire guide into the stent positioner.
NOTE: The lockable fitting on the distal end of the stent positioner must be loosened to allow the wire guide to be introduced.
2. Carefully introduce the wire guide into the blunt end of the stent and advance it until the stent is straightened and the tip of the wire guide extends beyond the tapered end of the stent.
3. While maintaining the position of the stent and the wire guide, advance the stent positioner until its blunt end rests against the blunt end of the stent.
4. Check to ensure the tip of the wire guide extends beyond the tip of the straightened stent and then tighten the lockable fitting to secure the position of the stent positioner on the wire guide.
5. Advance the assembly to the desired location using fluoroscopy and direct vision. Stop the advancement of the assembly once the last ink mark on the stent is at the ureteral orifice.
NOTE: If the stent is inadvertently advanced too far, the attached tether may be used to gently pull the stent back down to the proper location.
6. Confirm fluoroscopically that the tip of the wire guide and the proximal end of the stent are in the renal pelvis, stabilize the positioner, loosen the lockable fitting, and gently withdraw the wire guide to allow the retentive coil to form in the renal pelvis. (**Fig 1**)
7. Maintaining the placement of the positioner, continue to gently withdraw the wire guide until the intravesicle coil forms.
8. Completely remove the positioner and the wire guide. (**Fig 2**)
9. Fluoroscopically confirm the position of the stent coils in the bladder and renal pelvis.
NOTE: If necessary, final adjustments can be made with endoscopic forceps. The stent may be removed by gentle withdrawal traction using endoscopic forceps.

Traditional Endoscopic Placement

1. With the wire guide of choice in the proper position within the renal pelvis, pass the stent over the wire guide through the cystoscope. Under direct vision, advance the stent into the ureter with the stent positioner. Have an assistant hold the wire guide in position to prevent advancement of the wire guide into the renal parenchyma.
2. Stop the advancement of the stent once the last ink mark on the stent is at the ureteral orifice.

3. Confirm fluoroscopically that the tip of the wire guide and the proximal end of the stent are in the renal pelvis, stabilize the positioner, loosen the lockable fitting, and gently withdraw the wire guide to allow the retentive coil to form in the renal pelvis. (Fig 1)
4. Maintaining the placement of the positioner, continue to gently withdraw the wire guide until the intravesicle coil forms.
5. Completely remove the positioner and the wire guide. (Fig 2)
6. Fluoroscopically confirm the position of the stent coils in the bladder and renal pelvis.

NOTE: If necessary, final adjustments can be made with endoscopic forceps. The stent may be removed by gentle withdrawal traction using endoscopic forceps.

Tether Options:

A. Removal after stent positioning:

After stent placement, hold the knot of the tether and cut one strand. While holding the knot, gently pull on the tether to remove, maintaining stent position with positioner and wire guide.

B. Stent removal:

The tether may be used for subsequent stent removal. To avoid possible tether encrustation, the tether should be removed if the stent is to remain indwelling longer than 14 days.

C. Removed prior to stent placement:

Remove tether by holding the knot, cutting one strand, and while holding the knot, gently pull on the tether.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2022
Bloomington, IN 47404 U.S.A.