en
3

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BALLOON URETERAL DILATOR AND BALLOON URETERAL DILATOR SET

Read all instructions carefully. Failure to properly follow the information provided may lead to the device not performing as intended or injury to the patient.

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

The Balloon Ureteral Dilator is manufactured with outer diameters of 5.0 or 7.0 Fr with a length of 65 cm. The 5.0 Fr balloons accept a 0.028-inch wire guide and the 7.0 Fr balloons accept a 0.038-inch wire guide. The balloon dilator is available with an inflated balloon diameter of 5.0 or 6.0 mm and a balloon length of 4.0 cm.

Balloon Ureteral Dilators consist of a balloon catheter, endoscopic cap, and syringe. Balloon Ureteral Dilator Sets also include a pin vise handle assembly and wire guide.

Performance Characteristics

- The balloon on the Balloon Ureteral Dilator has radiopaque markers at the proximal and distal ends of the balloon to provide visibility.
- The endoscopic cap can be attached to the accessory channel of the cystoscope to seal around the balloon dilator.
- The pin vise provided in the Balloon Ureteral Dilator Set stabilizes the wire guide during introduction through the catheter.

Device Compatibility

- Syringe
- Sterile saline
- Contrast media
- Water-soluble lubricant
- Cystoscope with appropriately sized working channel
- 0.028" wire guide (5.0 Fr devices)
- 0.038" wire guide (7.0 Fr devices)

Patient Population

These devices are not intended to be used in any specific target populations. Patient populations could include any adult who requires treatment according to the device's intended use and as directed by a qualified physician (or properly licensed practitioner).

Intended Users

These devices are used by physicians or healthcare practitioners with proper authorization to perform medical procedures (e.g., licensing, qualifications, credentials) per their local administrative and regulatory guidelines.

Contact with Body Tissue

These devices interact with the structures and tissues comprising the trans-urethral pathway to the ureter, which includes the urethra, urinary bladder, ureteral orifice, and ureter.

Operating Principle

The Balloon Ureteral Dilator is constructed of two coaxially assembled polyethylene tubes with a polyethylene balloon at the distal end. Of the two coaxial tubes, the inner tube is comprised of yellow radiopaque polyethylene tubing and the outer tube is comprised of a clear polyethylene tubing. The Balloon Ureteral Dilator is shipped with a protective sleeve made of fluorinated ethylene propylene over the balloon portion of the dilator. The protective sheath is to be discarded prior to the use of the device.

The balloon on the Balloon Ureteral Dilator has radiopaque markers at the proximal and distal ends of the balloon to provide visibility.

The balloon dilator is manufactured with outer diameters of 5.0 or 7.0 Fr with a length of 65 cm. The balloon dilator is available with an inflated balloon diameter of 5.0 or 6.0 mm and a balloon length of 4.0 cm.

The endoscopic cap can be attached to the accessory channel of the cystoscope to seal around the balloon dilator once it is introduced into the cystoscope. The balloons are also provided with a 10 mL Luer-lock syringe for balloon inflation.

Depending on the Fr size of the balloon dilator in the Balloon Ureteral Dilator Sets, the wire guide provided is either a 0.028-inch (with 5.0 Fr Balloon Dilator) or a 0.038-inch (with 7.0 Fr Balloon Dilator) outer diameter stainless steel straight safety wire guide with a length of 145 cm.

Also provided in the Balloon Ureteral Dilator Sets is a pin-vise handle to stabilize the wire guide during introduction through the dilator.

INTENDED USE

These balloon dilators and sets facilitate ureteral dilation to permit completion of subsequent urologic procedures.

DEVICE AND CLINICAL BENEFITS

These devices facilitate minimally invasive dilation of the ureter.

INDICATIONS FOR USE

These devices are indicated for use in adult patients with urinary tract conditions that require ureteral dilation.

CONTRAINDICATIONS

None.

WARNINGS

- Always inflate the balloon with a sterile liquid. Never inflate with air, carbon dioxide, or any other gas.
- Removal of the catheter is facilitated by rotating the shaft counterclockwise during withdrawal. Using excessive force to withdraw the balloon can inflict trauma to tissue and/or damage the cystoscope.
- Imaging and/or placement procedures may expose the user and patient to radiation. The user should consider potential risks associated with radiation and appropriately plan procedures to minimize exposure in patients and medical personnel. The occurrence and/or severity of undesirable side effects may be dose dependent or probabilistic.
- The single use device is not designed for re-use. Attempts to reprocess (re-sterilize) and/or to re-use may lead to chemical contamination, device failure and/or transmission of disease.
- Do not use the device if the sterile packaging is damaged or unintentionally opened before use.

PRECAUTIONS

- Do not pre-inflate the balloon prior to insertion.
- To ensure proper regulation of balloon pressure, use of a balloon inflation device and pressure gauge is recommended
- Refer to product label for maximum inflation pressure. Do not overinflate the balloon.
-  This symbol on the label indicates that the device contains cobalt (Co) at a level above 0.1%w/w, which is a substance that is toxic to reproduction and a carcinogen (Class IB). However, the device contains cobalt as a stainless-steel alloy containing cobalt, which does not cause an increased risk of cancer or adverse reproductive effect according to current scientific evidence.

POTENTIAL ADVERSE EVENTS

- Bleeding
- Infection
- Injury to (or inflammation of) surrounding structures
- Pain

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Keep dry and away from sunlight. Upon removal from package, inspect the product to ensure no damage has occurred.

INSPECTION OF DEVICE

Visually inspect the device thoroughly including all levels of the packaging (as applicable) to verify that there is no damage prior to use. Visually inspect and confirm that the integrity of the sterile barrier has not been compromised in any way.

INSTRUCTIONS FOR USE

1. Remove the protective sleeve from the balloon and discard.

CAUTION: Do not pre-inflate the balloon.

2. Lock the syringe onto the side-arm stopcock of the balloon catheter and aspirate until all air is removed. Close the stopcock.
3. Lubricate the deflated balloon with a water-soluble lubricant.
4. Attach the endoscopic cap to the accessory channel of the cystoscope.
5. Under X-ray or fluoroscopic control, cystoscopically pass the wire guide the desired distance into the ureter.
6. Pass the balloon catheter over the wire guide. Position the tip of the balloon within the distal ureter, being sure to leave a small portion of the balloon in the ureteral orifice and extending into the bladder.

NOTE: Alternatively, the balloon and catheter can be passed through the cystoscope as a unit, following the instructions below:

- a. Attach the pin vise handle to the female Luer lock fitting of the balloon catheter.
 - b. Place the wire guide, flexible tip leading, into the central lumen and pass it the length of the catheter until only the 3-cm flexible tip extends beyond the catheter tip.
 - c. Maintaining position of the wire guide, twist the pin vise handle clockwise to lock the wire guide in place.
 - d. Cystoscopically pass the wire guide/catheter assembly so that the tip of the balloon is within the distal ureter, being sure to leave a small portion of the balloon in the ureteral orifice and extending into the bladder.
7. Fill the syringe with dilute (30% solution) contrast medium and eliminate air from the syringe.

CAUTION: Always inflate the balloon with a sterile liquid. NEVER inflate with air, carbon dioxide or any other gas.

8. Open the stopcock on the inflation device and inflate the balloon.

NOTE: To ensure proper regulation of balloon pressure, use of a balloon inflation device and pressure gauge is recommended.

CAUTION: Do not exceed 60 psi on the 7.0 Fr balloon ureteral dilators and 40 psi on the 5.0 Fr balloon ureteral dilators.

9. Close the stopcock for at least 30 seconds to maintain pressure and allow adequate dilation.

NOTE: Time, rather than excessive pressure, is the key factor.

10. Use the syringe to aspirate balloon contents and deflate the balloon. Reposition the balloon at a higher level and repeat the process until the desired location is reached.
11. Aspirate the balloon completely, then gently remove the device.

CAUTION: Removal of the catheter is facilitated by rotating the shaft counterclockwise during withdrawal. Using excessive force to withdraw the balloon can inflict trauma to tissue and/or damage the cystoscope.

DISPOSAL OF DEVICES

After the procedure, this device may be contaminated with potentially infectious substances of human origin and should be disposed of in accordance with institutional guidelines.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

SERIOUS INCIDENT REPORTING

If any serious incident has occurred in relation to the device, this should be reported to Cook Medical, and this shall also be reported to the competent authority of the member state in which the device was used.

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© COOK 2025
Cook Medical Europe Ltd
William A. Cook Australia Pty Ltd