**COOK**

**MEDICAL**

**EN**  
**3**

**Amplatz Sheaths**

Instructions for Use

**CS**  
**4**

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**DA**  
**5**

**Amplatz-sheaths**

Brugsanvisning

**DE**  
**6**

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**EL**  
**8**

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T - C A S - R E V O

## Amplatz Sheaths

**CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner). Please read all instructions before using this device.**

### DEVICE DESCRIPTION

The Amplatz sheaths are constructed of either radiopaque PTFE material (VS- prefix) or vinyl with a radiopaque stripe that runs the length of the device (CAS- prefix). The distal end of the sheath is cut at a 45-degree angle, and the proximal end is cut at a 90-degree angle. The sheaths are available in a variety of lengths and French sizes.

### INTENDED USE/INDICATION FOR USE

This device is intended to be used as a working channel to maintain a previously established nephrostomy tract.

### CONTRAINDICATIONS

This device is contraindicated in the presence of conditions which create unacceptable risk during percutaneous nephrostomy.

### WARNINGS

Sterilized by Ethylene Oxide gas. Do not use if sterile barrier is damaged. The device is intended for one-time use.

### PRECAUTIONS

- This product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of devices in angiographic procedures should be employed.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

### POTENTIAL ADVERSE EVENTS

- Tissue trauma
- Bleeding
- Perforation
- Infection
- Urine leakage
- Fistula
- Pain

### INSTRUCTIONS FOR USE

1. Remove the support tubing from the sheath before use and discard.
2. With a nephrostomy tract established by an Amplatz Dilator of the same French size as the sheath, load the sheath over the dilator.
3. While maintaining dilator's position, advance the sheath into the calyx.
4. Use fluoroscopy to confirm the location of the sheath in the calyx.
5. While holding the sheath in place, remove the dilator.
6. If a wire guide anchor is needed and the sheath is equipped with a notch, slide a wire guide into the notch on the proximal end of the sheath.

## HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

## REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2019
Bloomington, IN 47404 U.S.A.