COOK

MEDICAL

CE
0088

EN
3

Osteo-Site® Coaxial Bone Biopsy Needle

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Bruksanvisning

T _ C B B N _ R E V O

OSTEO-SITE® COAXIAL BONE BIOPSY NEEDLE

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Osteo-Site Coaxial Bone Biopsy Needle contains a stainless steel biopsy needle cannula with an accompanying stainless steel obturator. The device is to be used in conjunction with either 11 or 13 gage Osteo-Site Vertebroplasty Needles.

INTENDED USE

The Osteo-Site Coaxial Bone Biopsy Needle is intended for use in bone access, biopsy and sampling procedures.

CONTRAINDICATIONS

None known

WARNINGS

None known

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in bone access and biopsy techniques. Standard techniques should be employed.
- Ensure that the lumen of the external cannula being used will accept the outer diameter of the cutting cannula.

INSTRUCTIONS FOR USE

1. Under fluoroscopic guidance, using standard technique, position the Osteo-Site Vertebroplasty Needle just proximal to the area to be biopsied.
2. Remove the obturator from the cannula of the Osteo-Site Vertebroplasty Needle.
3. Advance the Osteo-Site Coaxial Bone Biopsy Needle coaxially down the outer cannula.
4. The Osteo-Site Coaxial Bone Biopsy Needle may be advanced past the tip of the vertebroplasty needle by rotating the coaxial needle cannula around its long axis while imparting forward pressure on the coaxial needle cannula.
5. Once the coaxial needle cannula is in the desired biopsy location, remove the coaxial needle. **NOTE:** A syringe may be attached to the cannula to perform aspiration during withdrawal if desired.
6. The biopsy specimen may be cleared from the coaxial needle by passing the obturator through the lumen of the coaxial needle.
7. If multiple specimens are desired, steps 3-6 may be repeated through the cannula.
8. If no further diagnostic or therapeutic procedures are desired, the cannula may be removed and appropriate dressing applied to the skin puncture site.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

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Bloomington, IN 47404 U.S.A.