- EN**
3 **Cook-Cope Loop Nephrostomy Set**
Instructions for Use
- CS**
5 **Souprava smyčkového nefrostomického katetru Cook-Cope**
Návod k použití
- DA**
7 **Cook-Cope nefrostomisæt med løkke**
Brugsanvisning
- DE**
9 **Cook-Cope Sperrschlaufen-Nephrostomie-Set**
Gebrauchsanweisung
- EL**
12 **Σετ νεφροστομίας με βρόχο Cook-Cope**
Οδηγίες χρήσης
- ES**
15 **Equipo de nefrostomía de lazo Cook-Cope**
Instrucciones de uso
- FR**
17 **Set à boucle de Cook-Cope pour néphrostomie**
Mode d'emploi
- HU**
20 **Cook-Cope nefrostomiás hurokkatéter-készlet**
Használati utasítás
- IT**
22 **Set con catetere nefrostomico a mono J Cook-Cope**
Istruzioni per l'uso
- NL**
24 **Cook-Cope lusnefrostomieset**
Gebruiksaanwijzing
- NO**
27 **Cook-Cope-slyngenefrostomisett**
Bruksanvisning
- PL**
29 **Zestaw cewnika nefrostomijnego Cook-Cope z pętlą**
Instrukcja użycia
- PT**
31 **Conjunto de cateter com ansa Cook-Cope para nefrostomia**
Instruções de utilização
- SV**
34 **Cook-Cope nefrostomiset med ögla**
Bruksanvisning

COOK-COPE LOOP NEPHROSTOMY SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner). Please read all instructions before using this device.

DEVICE DESCRIPTION

- 21 gage stainless steel needle, 15 cm long
- .018 inch (0.46 mm) diameter stainless steel wire guide, 60 cm long, with flexible coudé tip
- 6.3 French radiopaque introducing catheter with 20 gage stainless steel stiffening cannula, 22 cm long
- .038 inch (0.97 mm) diameter TFE-coated stainless steel wire guide, 100 cm long, with 3 mm curved Safe-T-J[®] tip
- Radiopaque polyethylene dilators, 20 cm long
- Radiopaque polyurethane loop nephrostomy catheter, 19.5 cm long, with latex cuff
- 18 gage stainless steel catheter introduction stiffening cannula, 35 cm long
- Silicone retention disc with pull tie
- 14.0 French clear polyvinylchloride connecting tube, 30 cm long, with one-way stopcock and drainage bag connector

NOTE: Set components may vary.

INTENDED USE

This device is intended for percutaneous placement of a loop catheter in the renal pelvis for nephrostomy drainage.

CONTRAINDICATIONS

This device should not be used in the presence of:

- Transitional cell carcinoma of the upper urinary tract
- Coagulopathy

WARNINGS

None known

PRECAUTIONS

- This product contains **natural rubber latex** which may cause allergic reactions in some individuals.
- Periodic evaluation is advised; the catheter should not remain indwelling more than four weeks.
- This device is intended for use by physicians trained and experienced in percutaneous access and placement of drainage catheters. Standard techniques should be employed.
- Do not forcefully advance any components during removal or replacement of the catheter. Carefully remove the components if any resistance is encountered.

MRI SAFETY INFORMATION

Nonclinical testing has demonstrated that the Cook-Cope Loop Nephrostomy Catheter is **MR Conditional** according to ASTM F2503. A patient with this device may be safely scanned after placement under the following conditions.

- Static magnetic field of 3.0 tesla or 1.5 tesla only
- Maximum magnetic field spatial gradient of 1900 gauss/cm (19T/m) or less

- Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of ≤ 2.0 W/kg (Normal Operating Mode)

Under the scan conditions provided above, the Cook-Cope Loop Nephrostomy Catheter is not expected to experience RF induced heating above background heating.

The image artifact extends approximately 2.4 mm from the Cook-Cope Loop Nephrostomy Catheter as found during nonclinical testing when imaged with a gradient echo pulse sequence and a 3.0 tesla MRI system.

INSTRUCTIONS FOR USE

Catheter Placement

1. Under fluoroscopic control, identify and anesthetize the skin site overlying the collecting system.
2. Introduce the 21 gage needle below the twelfth rib, about five finger-breadths from the midline. Advance the needle in short steps toward the lower pole until the needle is seen moving with respiration. While the patient holds his or her breath, advance the needle no further than 3 cm. If urine is not aspirated, repeat the process using a slightly different angle until a calyx is successfully punctured.
3. Thread the .018 inch (0.46 mm) diameter wire guide through the needle; and maneuver the wire guide into the pelvis or proximal ureter.
4. Remove the 21 gage needle and advance the 6.3 French introducing catheter with 20 gage cannula over the wire guide into the calyx. Holding the cannula stationary, further advance the catheter to follow the wire guide into the pelvis and ureter.
5. Remove the wire guide and cannula. Decompress and irrigate the collecting system.
6. Align the curved tip of the Safe-T-J wire guide with the curve of the introducing catheter. Introduce the flexible tip of the wire guide into the catheter and advance it until the wire guide exits the sideport of the introducing catheter. If the advancing wire guide encounters resistance, withdraw it slightly and then advance it with a slight rotary motion.
7. When the tip of the Safe-T-J wire guide is within the renal pelvis, carefully withdraw the catheter.
NOTE: Care must be exercised to avoid withdrawal of the wire guide.
8. To facilitate passage of the Cope nephrostomy catheter, dilate the musculofascial tract by progressing from the smallest to the largest dilator provided.
9. **IMPORTANT: To avoid breaking the catheter tether, unwind several centimeters of tether from the attached spool and gently straighten the catheter loop.** Pass the Cope nephrostomy catheter over the external end of the wire guide; gradually advance the catheter into the collecting system. If it becomes difficult to advance the catheter, insert the catheter introduction stiffening cannula over the wire guide and into the catheter. Some lubricant may be helpful.
10. Confirm catheter position fluoroscopically. While holding the shaft of the catheter securely with one hand, withdraw the wire guide (and stiffening cannula, if used) with the other.
11. While pulling the drawstring, maneuver the catheter back and forth under fluoroscopic monitoring until there is a looping of the distal tip. Wrap the drawstring around the catheter shaft; and cut the loose ends. Cover the knot by unrolling the rubber sleeve.
12. Attach the retention disc to the catheter shaft 1 to 2 cm from the skin surface. To prevent the catheter from being accidentally displaced forward, fix the retention disc in position by gluing or suturing.
13. The connecting tube provides a convenient, secure method of attachment to standard urine drainage pouches. The optional stopcock may be used for intermittent patency check.

Catheter Removal

1. Transect the catheter and tether immediately below the latex cuff to free the catheter loop.
2. Gently withdraw the catheter. The soft, retentive coil will straighten during removal.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2023
Bloomington, IN 47404 U.S.A.