EN 2	Angiography Catheters Instructions for Use
BR 3	Cateteres angiográficos Instruções de uso
CS 4	Angiografické katetry Návod k použití
DA 6	Angiografikatetre Brugsanvisning
DE 7	Angiographiekatheter Gebrauchsanweisung
EL 9	Αγγειογραφικοί καθετήρες Οδηγίες χρήσης
ES 11	Catéteres angiográficos Instrucciones de uso
FR 12	Cathéters d'angiographie Mode d'emploi
HU 14	Angiográfias katéterek Használati utasítás
IT 15	Cateteri angiografici Istruzioni per l'uso
NL 17	Angiografiekatheters Gebruiksaanwijzing
NO 18	Angiografikatetre Bruksanvisning
PL 20	Cewniki angiograficzne Instrukcja użycia
PT 21	Cateteres angiográficos Instruções de utilização
SV 23	Angiografikatetrar Bruksanvisning

ANGIOGRAPHY CATHETERS (CEREBRAL)

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

Catheters for use in angiographic procedures are available in a variety of French sizes, endhole sizes, lengths, materials and designs (e.g., polyethylene or nylon, non-braided or braided with 1:1 torque).

The intended purpose of the cervicocerebral catheter is to gain access to vessels that can be appropriately accessed in the head, neck, brain, and spine, allowing for angiographic procedures.

The choice of catheter material for different angiographic procedures should be based on the physician's experience.

INTENDED USE / INDICATIONS FOR USE

The catheters are intended for use in angiographic procedures by physicians trained and experienced in angiographic techniques. Standard techniques for placement of vascular access sheaths, angiographic catheters and wire guides should be employed.

CONTRAINDICATIONS

There are no known absolute contraindications.

WARNINGS

None known.

POTENTIAL ADVERSE EVENTS

- Access site injury
- Thrombosis/embolism/transient ischemic attack/stroke
- Death
- Infection
- Bleeding/hematoma/hemorrhage
- Vessel Injury (e.g., perforation, rupture, pseudoaneurysm, dissection, arteriovenous fistula)
- Vessel spasm
- Vessel occlusion (e.g., air embolism, thromboembolism, thrombosis, vessel dissection, disruption of atherosclerotic plaque)
- Renal failure (e.g., contrast-induced nephropathy)
- Allergic reaction to contrast media

PRECAUTIONS

- Manipulation of catheter requires fluoroscopic control.
- Due to thinwall construction, extreme care must be exercised during manipulation and withdrawal. Catheter insertion through a synthetic vascular graft should be avoided whenever possible.
- Do not attempt to heat or reshape the catheter curve; the catheter tip is made from a heat-sensitive material.
- The possible whiplash effect of the long, soft catheter tip must be considered during selective angiography.
- Activate hydrophilic coating, if present, by wetting the distal end of the catheter with sterile water or saline. For best results, keep the catheter surface wet during placement.
- If an over-the-wire catheter design is used, where the wire is left in place during injections, the proximal end of the wire should be clamped to protect against accidental injection of the wire guide.

INSTRUCTIONS FOR USE

The catheter has been designed for percutaneous introduction into the vascular system over an appropriately sized wire guide or through an appropriately sized sheath introducer.

Due to the catheter construction, the use of a wire guide is recommended during advancement.

If resistance is encountered during manipulation, stop and determine the cause before proceeding any further.

Prior to advancement of the cerebral catheter, confirm that the vessel can sufficiently accommodate the 1.6764 mm diameter catheter.

Peel-Away® Straightener

Catheters with pigtail shaped or similar configurations are supplied with Peel-Away straighteners to be used as follows.

1. Slide Peel-Away straightener forward to open catheter tip. (Fig. 1)

Fig. 1

2. With catheter tip straight, advance catheter over wire guide. (**Fig. 2**)


Fig. 2

3. After initial advancement of catheter over the wire guide, pull back and peel away the straightener for catheter introduction. (**Fig. 3**)


Fig. 3

Peel-Away straightener must not be used as a vascular introducer sheath.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd