COOK

MEDICAL

CE
0088

EN Hilal Embolization Microcoils™

4 Instructions for Use

DA Hilal Microcoils™ til embolisering

5 Brugsanvisning

DE Hilal-Embolisationsspiralen Microcoils™

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Fig. 3

HILAL EMBOLIZATION MICROCOILS™

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

Hilal Embolization Microcoils are made of platinum with spaced synthetic fibers and supplied preloaded in a loading cartridge.

INTENDED USE

Hilal Embolization Microcoils are intended for embolization of selective vessel supply to arteriovenous malformations and other vascular lesions of the brain, spinal cord and spine.

CONTRAINDICATIONS

None known

WARNINGS

- Embolization Microcoils should be positioned with particular care. Microcoils should not be left too close to the inlets of arteries and should be intermeshed with previously placed Microcoils if possible. A minimal, but sufficient, arterial blood flow should remain to hold the Microcoils against the previously placed Microcoils until a solid clot ensures permanent fixation. The purpose of these suggestions is to minimize the possibility of loose Microcoils becoming dislodged and obstructing a normal and essential arterial channel.
- Embolization Microcoils are not recommended for use with polyurethane catheters or catheters with sideports. If a catheter with sideports is used, the embolus may lodge in the sideport or pass inadvertently through it. Use of a polyurethane catheter may also result in lodging of the embolus within the catheter.

PRECAUTIONS

- The product is intended for use by physicians trained and experienced in embolization techniques. Standard techniques for placement of vascular access sheaths, angiographic catheters, and wire guides should be employed.
- To determine correct catheter position, perform an angiogram prior to embolization.
- Prior to introduction of the Microcoil, flush the angiographic catheter with saline.
- Manipulation of products requires fluoroscopic control.

PRODUCT RECOMMENDATIONS

- Microcoils are recommended for use through catheters designed for use with .018 inch (0.46 mm) diameter wire guides and whose inner diameter (ID) does not exceed .027 inch (0.69 mm) diameter.
- Wire guides recommended for loading and positioning Microcoils are TFE-coated .018 inch (0.46 mm) diameter with flexible tapered tips.
- Microcoils may be used in conjunction with particulate or liquid embolization materials. Final positioning of Microcoils creates a “platinum cast” effect within the vessel lumen.

INSTRUCTIONS FOR USE

1. Perform an angiogram prior to embolization to determine optimal catheter position.
2. Firmly grasp Microcoil loading cartridge between thumb and forefinger at point where right-angle shipping stylet exits. (**Fig. 1**)
3. While maintaining firm finger grip, remove shipping stylet. This will prevent Microcoil from exiting cartridge. Verify its position inside cartridge by direct vision.
4. Position loading cartridge into base of catheter hub. (**Fig. 2**)
5. Using .018 inch (0.46 mm) diameter wire guide, push Microcoil out of loading cartridge and into catheter lumen. (**Fig. 3**)
6. Remove loading cartridge.
7. Once positioned into delivery catheter, coil design permits delivery into the target vessel by saline flush or by push technique using an appropriately sized wire guide or pusher.
8. Perform final angiogram to confirm coil position within target vessel.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2012
Bloomington, IN 47404 U.S.A.