## **MReye® Embolization Coil**

Instructions for Use

**EN  
4**

## **MReye® emboliseringscoil**

Brugsanvisning

**DA  
6**

## **MReye® Embolisationsspirale**

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**DE  
8**

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

**Fig. 2**

## **MREYE® EMBOLIZATION COIL**

**MR Conditional with an MR scanner not exceeding 3.0 Tesla.**

**CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).**

### **DEVICE DESCRIPTION**

MReye Embolization Coils are constructed of an MR conditional super alloy, which has similar physical characteristics to stainless steel, but has no ferromagnetic properties. The coils are easily detected radiographically and are considered to be MR conditional, with the potential for minor artifacts to be observed in the immediate vicinity of the device. The coils feature spaced synthetic fibers to maximize thrombogenicity. Embolization coils are delivered to the target vessel using a soft, straight wire guide through a standard angiographic catheter.

### **INTENDED USE**

MReye Embolization Coils are intended for use in peripheral arterial and venous vessel embolization procedures.

### **MR STATEMENTS**

Through non-clinical testing, embolization coils manufactured of this super alloy have been shown to be MR conditional at static magnetic field strengths of 3.0 Tesla or less, a maximum spatial gradient of 330 gauss/cm, and a maximum whole body averaged specific absorption rate (SAR) of 2.0 W/kg for 20 minutes of MRI. Embolization coils manufactured of this MR conditional super alloy will not migrate in this MRI environment. Non-clinical testing has not been performed to rule out the possibility of device migration at static magnetic field strengths higher than 3.0 Tesla and a maximum spatial gradient higher than 330 gauss/cm.

In this testing, embolization coils manufactured of this MR conditional super alloy produced a temperature rise of  $\leq 0.6$  degrees C at a maximum whole body averaged specific absorption rate (SAR) of 2.0 W/kg for 20 minutes of MRI. The effect of heating in the MRI environment for other conditions, multiple coils or overlapping coils is unknown.

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the embolization coils manufactured of this MR conditional super alloy.

### **CONTRAINDICATIONS**

None known

### **WARNINGS**

- Positioning of Embolization Coils should be done with particular care. Coils should not be left too close to the inlets of arteries and should be intermeshed with previously placed coils if possible. A minimal but sufficient arterial blood flow should remain to hold the coils against the previously placed coils until a solid clot ensures permanent fixation. The purpose of these suggestions is to minimize the possibility of loose coils becoming dislodged and obstructing a normal and essential arterial channel.
- MReye Embolization Coils are not intended for neurovascular use.
- MReye Embolization Coils are not recommended for use with polyurethane catheters or catheters with sideports. If a catheter with sideports is used, the embolus may lodge in the sideport or pass inadvertently through it. Use of a polyurethane catheter may also result in lodging of the embolus within the catheter.
- If difficulties occur when deploying the embolization coil, withdraw the wire guide, coil and angiographic catheter simultaneously as a unit.

### **PRECAUTIONS**

- This product is intended for use by physicians trained and experienced in embolization techniques. Standard techniques for placement of vascular access sheaths, angiographic catheters and wire guides should be employed.
- Perform an angiogram prior to embolization to determine correct catheter position.
- Prior to introduction of the embolization coil, flush the angiographic catheter with saline.

PRODUCT RECOMMENDATIONS

Delivery Catheter and Wire Guide Recommendations for Embolization Coils

| Coil Size Diameter | Wire Guide Type & Size | Catheter Type & Size |
|--------------------|------------------------|----------------------|
| .025 inch          | TSFB-25                | SCR3.0B              |
|                    | TSF-25                 |                      |
| .035 inch          | TSFNB-35               | HNB(R)4.1-35         |
|                    | TSFNA-35               | HNB(R)5.0-35         |
|                    | TSFB-35                | HNB6.0-35            |
|                    | TSFBP-35               | SCBR4.0-35           |
|                    | TSF-35                 | SCBR5.0-35           |
| .038 inch          | TSFNB-38               | HNB(R)5.0-38         |
|                    | TSFNA-38               | HNB6.0-38            |
|                    | TSFB-38                | SCBR4.0-38           |
|                    | TSF-38                 | SCBR5.0-38           |

INSTRUCTIONS FOR USE

**NOTE:** The embolization coil is supplied preloaded in the loading cartridge.

1. Perform an angiogram prior to embolization to determine optimal catheter position.
2. Firmly grasp the loading cartridge between thumb and forefinger. Introduce the metal end of the loading cartridge into the base of the catheter hub. Lock loading cartridge onto catheter hub by rotating Luer lock adapter clockwise. (**Fig. 1**)
3. Maintaining position of the cartridge, advance the stiff portion of the wire guide into the loading cannula. Push the coil into the first 20 to 30 centimeters of the angiographic catheter. (**Fig. 2**) Remove the wire guide and loading cartridge.
4. With the flexible tip of the wire guide, advance the embolization coil to the tip of the catheter. Verify position of the angiographic catheter prior to deployment.
5. Deploy the coil by advancing the wire guide past the tip of the catheter.
6. Perform final angiogram to confirm coil position within target vessel.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

V. P. Chuang, S. Wallace, C. Gianturco: "A New Improved Coil for Tapered Tip Catheter for Arterial Occlusion," *Radiology*, 135 (1980), 507-509.

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Bloomington, IN 47404 U.S.A.