EN
8

MReye® Embolization Coil

Instructions for Use

CS
10

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DA
11

MReye® emboliseringscoil

Brugsanvisning

DE
14

MReye® Embolisationsspirale

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EL
16

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- 1 Vävning
- 2 Hylsa
- 3 Kateter


Fig. 5

- 1 Weaving
- 2 Sheath
- 3 Catheter
- 4 Inconel Coil (High Radial Force)
- 5 Inconel Coil Deployed (Scaffold)
- 6 Platinum Coil (Soft Coil)
- 7 Platinum Coil Deployed, Cross-Sectional Occlusion Completed

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MREYE® EMBOLIZATION COIL

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

MReye Embolization Coils are made of Inconel, an MR conditional super alloy with spaced synthetic fibers and are supplied preloaded in a loading cartridge. They are designed to be delivered to the target vessel using a soft, straight wire guide through a standard angiographic catheter.

INTENDED USE/INDICATIONS FOR USE

MReye Embolization Coils are intended for use in peripheral arterial and venous vessel embolization procedures.

MRI INFORMATION

Nonclinical testing has demonstrated that single and multiple MReye® Embolization Coils are **MR Conditional** according to ASTM F2503. A patient with this device may be safely scanned after placement under the following conditions.

- Static magnetic field of 3.0 tesla or 1.5 tesla only
- Maximum magnetic field spatial gradient of 1600 gauss/cm (16 T/m) or less
- Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of ≤ 2.0 W/kg (Normal Operating Mode)

Under the scan conditions provided above, the MReye® Embolization Coil is not expected to result in a temperature rise of more than **3.1°C** after 15 minutes of continuous scanning.

The image artifact extends approximately 12 mm from the MReye® Embolization Coil as found during nonclinical testing when imaged with a gradient echo pulse sequence and a 3.0 Tesla MRI system.

For US Patients Only

Cook recommends that the patient register the MR conditions disclosed in this IFU with the MedicAlert Foundation. The MedicAlert Foundation can be contacted in the following manners:

Mail: MedicAlert Foundation International
2323 Colorado Avenue
Turlock, CA 95382, USA

Phone: 888-633-4298 (toll free)
209-668-3333 from outside the US

Fax: 209-669-2450

Web: www.medicalert.org

CONTRAINDICATIONS

None known

WARNINGS

- Positioning of embolization coils should be done with particular care. Coils should not be left too close to the inlets of arteries and should be intermeshed with previously placed coils if possible. A minimal but sufficient arterial blood flow should remain to hold the coils against the previously placed coils until a solid clot ensures permanent fixation. The purpose of these suggestions is to minimize the possibility of loose coils becoming dislodged and obstructing a normal and essential arterial channel.
- MReye Embolization Coils are not intended for neurovascular use.
- MReye Embolization Coils are not recommended for use with polyurethane catheters or catheters with sideports. If a catheter with sideports is used, the embolization coil may lodge in the sideport or pass inadvertently through it. Use of a polyurethane catheter may also result in lodging of the embolization coil within the catheter.
- If difficulties occur when deploying the embolization coil, withdraw the wire guide, coil and angiographic catheter simultaneously as a unit.

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in embolization techniques. Standard techniques for placement of vascular access sheaths, angiographic catheters and wire guides should be employed.
- Perform an angiogram prior to embolization to determine correct catheter position.
- Prior to introduction of the embolization coil, flush the angiographic catheter with saline.

PRODUCT RECOMMENDATIONS

Delivery Catheter and Wire Guide Recommendations for Embolization Coils

Coil Size Diameter	Catheter Type & Size	Wire Guide Type & Size
.025 inch	SCR3.0B	TSF-25 TSFB-25
.035 inch	HNB4.0-35 HNB4.1-35 HNB[R]5.0-35 SCBR4.0-35 SCBR5.0-35	TSF-35 TSFB-35 TSFBP-35 TSFNA-35 TSFNB-35
.038 inch	HNB[R]5.0-38 HNB[R]6.0-38 SCBR4.0-38 SCBR5.0-38	TSF-38 TSFB-38 TSFNA-38 TSFNB-38

Coil Delivery Technique and Coil Size Selection

Long-term occlusion depends on achieving cross-sectional occlusion of the blood vessel, and coaxial catheters provide the ability to control placement of coils and permanent occlusion. The combination of the coaxial technique and either the anchor or scaffold technique significantly enhances stability of coil deployment.

- **Coaxial technique:** The use of an outer guiding sheath/catheter is the most important step in preventing coil elongation and uncertain long-term occlusion. The outer guiding sheath/catheter provides support, and the inner catheter provides finer selective maneuvers. (**Fig. 3**)
- **Anchor technique:** The anchor technique provides safe and distal occlusion when there is a question about instability of coils. At least 2 cm of a coil is advanced into the side branch, which is normally sacrificed. The rest of the coil is then deployed just proximal to that side branch, and additional coils are packed. (**Fig. 4**) **Note :** Use a coil longer than 10cm for anchor technique.
- **Scaffold technique:** The scaffold technique is used for high-flow vessels when there is concern about migration of a softer coil. An Inconel coil is placed initially. Then, several Inconel coils or platinum coils (soft coils) may be packed within the scaffold. (**Fig. 5**)

In general, the first coil selected should have a diameter that is 20% larger, or at least 2 mm oversized, than the vessel that is being occluded.

INSTRUCTIONS FOR USE

NOTE: The embolization coil is supplied preloaded in the loading cartridge.

1. Perform an angiogram prior to embolization to determine optimal catheter position.
2. Firmly grasp the loading cartridge between thumb and forefinger. Introduce the metal end of the loading cartridge into the base of the catheter hub. Lock loading cartridge onto catheter hub by rotating Luer lock adapter clockwise. (**Fig. 1**)
3. Maintaining position of the cartridge, advance the stiff portion of the wire guide into the loading cannula. Push the coil into the first 20 to 30 centimeters of the angiographic catheter. (**Fig. 2**) Remove the wire guide and loading cartridge.
4. With the flexible tip of the wire guide, advance the embolization coil to the tip of the catheter. Verify position of the angiographic catheter prior to deployment.
5. Deploy the coil by advancing the wire guide past the tip of the catheter.
6. Perform final angiogram to confirm coil position within target vessel.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

V. P. Chuang, S. Wallace, C. Gianturco: "A New Improved Coil for Tapered Tip Catheter for Arterial Occlusion," *Radiology*, 135 (1980), 507-509.

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Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd