Nester® Embolization Coils	Instructions for Use	EN 4
Embolizační spirálky Nester®	Návod k použití	CS 5
Nester® emboliseringscoils	Brugsanvisning	DA 6
Nester®-Embolisationsspiralen	Gebrauchsanweisung	DE 7
Σπειράματα εμβολής Nester®	Οδηγίες χρήσης	EL 9
Espiraless de embolización Nester®	Instrucciones de uso	ES 10
Spirales d'embolisation Nester®	Mode d'emploi	FR 12
Nester® embolizációs spirálok	Használati utasítás	HU 13
Spirali per embolizzazione Nester®	Istruzioni per l'uso	IT 14
Nester® embolisatiecoils	Gebruiksaanwijzing	NL 16
Nester® emboliseringsspiraler	Bruksanvisning	NO 17
Spirale embolizacyjne Nester®	Instrukcja użycia	PL 18
Espiraiss de embolização Nester®	Instruções de utilização	PT 20
Nester® emboliseringsspiraler	Bruksanvisning	SV 21
Nester® 栓塞弹簧圈	使用说明	ZH 22

CE
0123

T _ C E _ N E C _ R E V 4

Fig. 1

Fig. 2

NESTER® EMBOLIZATION COILS

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

Nester Embolization Coils are made of platinum with spaced synthetic fibers, and are supplied preloaded in a loading cartridge. They are designed to be delivered to the target vessel using a soft, straight wire guide through a standard angiographic catheter.

INTENDED USE

Nester Embolization Coils are intended for arterial and venous vessel embolization procedures. The product is intended for use by physicians trained and experienced in embolization techniques. Standard techniques for placement of vascular access sheaths, angiographic catheters and wire guides should be employed.

CONTRAINDICATIONS

None known.

WARNINGS

- Positioning of Embolization Coils should be done with particular care. Coils should not be left too close to the inlets of arteries and should be intermeshed with previously placed coils if possible. A minimal but sufficient arterial blood flow should remain to hold the coils against the previously placed coils until a solid clot ensures permanent fixation. The purpose of these suggestions is to minimize the possibility of loose coils becoming dislodged and obstructing a normal and essential arterial channel.
- Nester Embolization Coils are not recommended for use with polyurethane catheters or catheters with sideports. If a catheter with sideports is used, the embolus may lodge in the sideport or pass inadvertently through it. Use of a polyurethane catheter may also result in lodging of the embolus within the catheter.
- If difficulties occur when deploying the embolization coil, withdraw the wire guide, coil and angiographic catheter simultaneously as a unit.

PRECAUTIONS

- Perform an angiogram prior to embolization to determine correct catheter position.
- Prior to introduction of the embolization coil, flush the angiographic catheter with saline.
- **If using a .018 inch MicroNester™ embolization coil, ensure that the delivery catheter has an internal diameter (ID) between .018 and .025 inch.**

PRODUCT RECOMMENDATIONS

Appropriately sized diameter TFE-coated wire guides with flexible tapered tips are recommended for positioning Nester Embolization Coils. The following table offers specific recommendations.

Coil Size Diameter	Wire Guide Type & Size	Catheter Type & Size
.018 inch	TSF-18 STF-18	MF-2.5-18 MF-2.8-NT MF-3.0-18*
.035 inch	TSFNB-35 TSFNA-35 TSFB-35 TSFBP-35 TSF-35	HNB(R)4.1-35 HNB(R)5.0-35 HNB6.0-35 SCBR5.0-35
.038 inch	TSFNB-38 TSFNA-38 TSFB-38 TSF-38	HNB(R)5.0-38 HNB6.0-38 SCBR5.0-38

*Available outside U.S.

INSTRUCTIONS FOR USE

1. Perform an angiogram prior to embolization to determine optimal catheter position.
2. Firmly grasp the loading cartridge between thumb and forefinger. Introduce the metal end of the loading cartridge into the base of the catheter hub. Lock loading cartridge onto catheter hub by turning Luer lock adapter clockwise. (**Fig. 1**)

3. Maintaining position of the cartridge, advance the stiff portion of the wire guide into the loading cannula. Push the coil into the first 20 to 30 cm of the angiographic catheter. **(Fig. 2)** Remove the wire guide and loading cartridge.
- NOTE:** If you are placing a .018 inch embolization coil, you will advance the coil into the first few centimeters of the angiographic catheter with the loading stylet that is provided with the coils.
4. With the flexible tip of the wire guide, advance the embolization coil to the tip of the catheter. Verify position of the angiographic catheter prior to deployment.
5. Deploy the coil by advancing the wire guide past the tip of the catheter.
6. Perform final angiogram to confirm coil position within target vessel.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2019
Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd