COOK

MEDICAL

EN
3

Check-Flo Performer® Introducers and Guiding Sheaths

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NL
20

Check-Flo Performer® introducers en geleidesheaths

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Bruksanvisning

T _ C F P G S N O C E _ R E V O

CHECK-FLO PERFORMER® INTRODUCERS AND GUIDING SHEATHS

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

Check-Flo Performer Introducers and Guiding Sheaths have a Check-Flo® valve and a radiopaque tip incorporated within the sheath material to identify the location of the sheath's distal tip. The intended purpose of these devices is to serve as a pathway for other devices to be inserted into the cardiovascular system, either through the vasculature or via a transseptal approach following transseptal puncture.

INTENDED USE / INDICATIONS FOR USE

Check-Flo Performer Introducers and Guiding Sheaths are indicated for the percutaneous introduction of therapeutic or diagnostic devices to all heart chambers, including the left atrium via transseptal puncture.

CONTRAINDICATIONS

None known

WARNINGS

- When the sheath is left in the vessel, a continuous heparinized infusion under pressure is strongly recommended through the sheath side port.
- Infusion through the side port should only be done after all air is removed from the unit.
- Aspiration of the side port is recommended when withdrawing the catheter, probe, or dilator to remove any fibrin deposition which may have accumulated in or on the tip of the sheath.
- Careful sheath manipulation must be performed in the presence of an implantable cardiac device of any kind to minimize the potential to displace or dislodge lead placement.
- Fluoroscopic monitoring of the location of the distal tip of the sheath using the radiopaque marker, especially when used in a transseptal approach, is recommended.
- Before withdrawing the sheath through tortuous anatomy, insert the performer introducer dilator to avoid possible breakage.

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of vascular access sheaths should be employed.
- Determine the maximum diameter of the instrument or catheter to be introduced in order to ensure that it will pass through the introducer and/or dilator.
- All instruments or catheters used with this product should move freely through the valve and introducer or dilator. Damage may result when the fit is tight.
- When inserting, manipulating, or withdrawing a device through the introducer assembly, always maintain introducer assembly position.
- Before removing or inserting devices through the introducer, aspirate through the side-arm of the valve to clear the introducer, then flush with heparinized saline.
- When inflating a balloon at, or close to, the introducer tip, ensure the balloon is not inside the distal tip of the introducer.
- When puncturing, suturing, or incising the tissue near the introducer, use caution to avoid damaging the introducer.
- If resistance is encountered during manipulation, stop and determine the cause before proceeding any further.
- The following conditions require that special care be taken when using this product involving the transseptal approach:
 - Enlarged aortic root
 - Marked right atrial enlargement
 - Small left atrium
 - Marked distortion of the thoracic configuration (e.g. kyphosis or scoliosis)
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

POTENTIAL ADVERSE EVENTS

Adverse reactions to cannulation of the peripheral vasculature and intracardiac placement of the sheath and dilator may include, but are not limited to:

- Infection
- Local nerve damage
- Perforation
- Dissection
- AV fistula formation

- Pseudoaneurysm formation
- Arrhythmias
- Hematoma
- Hemorrhage
- Thromboembolic events
- Catheter entrapment
- Valve damage
- Pacemaker/defibrillator lead displacement
- Air embolus
- Vasovagal reaction
- Vessel trauma
- Vessel spasm
- Atrial septal defect
- Aortic puncture
- Perforation and/or tamponade
- Coronary artery spasm and/or damage
- Stroke
- Myocardial infarction
- Pericardial/pleural effusion
- Pulmonary edema
- Death

INSTRUCTIONS FOR USE

Sheath Introduction

1. Upon removal from package, ensure the inner diameter (ID) of the introducer and/or dilator is appropriate for the maximum diameter of the instrument or catheter to be introduced.
2. Using the side-arm of the valve, flush the introducer by filling the introducer assembly completely with heparinized saline.
3. Flush the dilator with heparinized solution.
4. Insert the dilator completely into the introducer.
5. Using standard Seldinger technique, access the target vessel with the appropriate needle.
6. Insert an appropriate wire guide into the vessel through the needle, then remove the needle, leaving the wire guide in place.
7. Insert the dilator/introducer combination over the wire guide under fluoroscopic guidance.
8. Use standard technique to position the introducer/dilator assembly into the desired location within the vasculature.
9. Remove the wire guide and dilator; aspirate and flush the introducer side-arm.
10. Insert an appropriately sized device as needed.
11. If transseptal puncture is required, refer to the Instructions for Use of the transseptal puncture device.

Sheath Removal

1. Insert a wire guide at least 10 cm past the tip of the introducer.
2. Insert the dilator over the wire guide into the introducer.
3. Withdraw the introducer and dilator as a unit.
4. After removal of the introducer, use standard technique to achieve hemostasis.
5. Remove the wire guide.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.