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
Fig. 4

ULTRATHANE® NEPHROSTOMY SET WITH MAC-LOC®

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner). Please read all instructions before using this device.

DEVICE DESCRIPTION

- Ultrathane® nephrostomy catheter with Mac-Loc® locking loop mechanism, AQ® hydrophilic coating on distal 5 cm, and radiopaque band at base of pigtail
- .038 inch (0.97 mm) diameter TFE-coated stainless steel Safe-T-J® wire guide, 100 cm long
- Non-radiopaque polyethylene flexible stylet
- Stainless steel rigid stylet
- .018 inch (0.46 mm) diameter stainless steel Cope mandril wire guide, 60 cm long
- Radiopaque polyethylene dilators, 20 cm long
- 21 gage stainless steel thinwall needle, 15 cm long
- Soft, low-profile catheter support with pull tie
- 6.3 French radiopaque polyethylene introducing catheter, 22 cm long
- 20 gage stainless steel stiffening cannula, 22 cm long
- 14.0 French polyvinylchloride connecting tube, 30 cm long, with one-way stopcock

INTENDED USE

This device is intended for percutaneous placement of a loop catheter in the renal pelvis for nephrostomy drainage.

CONTRAINDICATIONS

This device should not be used in the presence of:

- Transitional cell carcinoma of the upper urinary tract
- Coagulopathy

WARNINGS

None known

PRECAUTIONS

- Periodic evaluation is advised; the catheter should not remain indwelling more than four weeks.
- This device is intended for use by physicians trained and experienced in percutaneous access and placement of drainage catheters. Standard techniques should be employed.
- Do not forcefully advance any components during removal or replacement of the catheter. Carefully remove the components if any resistance is encountered.

MRI SAFETY INFORMATION

Nonclinical testing has demonstrated that the Ultrathane Nephrostomy Catheter with Mac-Loc is **MR Conditional** according to ASTM F2503. A patient with this device may be safely scanned after placement under the following conditions.

- Static magnetic field of 3.0 tesla or 1.5 tesla only

- Maximum magnetic field spatial gradient of 1900 gauss/cm (19 T/m) or less
- Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of ≤ 2.0 W/kg (Normal Operating Mode)

Under the scan conditions provided above, the Ultrathane Nephrostomy Catheter with Mac-Loc is not expected to result in a temperature rise of more than 2.5°C after 15 minutes of continuous scanning.

The image artifact extends approximately 40.1 mm from the Ultrathane Nephrostomy Catheter with Mac-Loc as found during nonclinical testing when imaged with a gradient echo pulse sequence and a 3.0 tesla MRI system.

INSTRUCTIONS FOR USE

NOTE: A TFE-coated wire guide must be used with this catheter.

1. Under fluoroscopic control, identify and anesthetize the skin site overlying the collecting system.
2. Introduce the 21 gage needle below the twelfth rib. Advance the needle in short steps toward the lower pole until the needle is seen moving with respiration. While the patient holds his or her breath, advance the needle no further than 3 cm. If urine is not aspirated, repeat the process using a slightly different angle until a calyx is successfully punctured.
3. Thread the .018 inch (0.46 mm) diameter wire guide through the needle; maneuver the wire guide into the pelvis or proximal ureter. Remove the 21 gage needle and advance the 6.3 French introducing catheter with 20 gage cannula over the wire guide into the calyx. Holding the cannula stationary, further advance the catheter to follow the wire guide into the pelvis and ureter.
4. Remove the wire guide and cannula. Decompress and irrigate the collecting system.
5. Align the curved tip of the Safe-T-J wire guide with the curve of the introducing catheter then introduce the flexible tip of the wire guide into the catheter and advance it until the wire guide exits the sideport of the introducing catheter. If the advancing wire guide encounters resistance, withdraw it slightly and then advance it with a slight rotary motion. When the tip of the Safe-T-J wire guide is within the renal pelvis, carefully withdraw the catheter. **NOTE:** Take care to avoid withdrawing the wire guide.
6. To facilitate passage of the nephrostomy catheter, dilate the musculofascial tract by progressing from the smallest to the largest dilator provided.

NOTE: Prior to use, immerse AQ catheter in sterile water or isotonic saline to allow the hydrophilic surface to absorb water and become lubricious. This will ease placement under standard conditions.

7. Pass the Ultrathane nephrostomy catheter over the external end of the wire guide; gradually advance the catheter into the collecting system. If it becomes difficult to advance the catheter, insert the 22 cm long stiffening cannula over the wire guide and into the catheter. Some lubricant may be helpful.
8. Confirm the position fluoroscopically, then hold the shaft of the catheter securely in position with one hand, withdraw the wire guide (and stiffening cannula, if used) with the other.
9. Lock the distal loop in place. **NOTE:** See the Mac-Loc locking loop mechanism instructions at the end of these instructions.
10. Attach the retention disc to the catheter shaft 1 to 2 cm from the skin surface and fix the retention disc in position by gluing or suturing. This will prevent the catheter from being accidentally displaced forward.
11. Use the connecting tube for attachment to standard urine drainage pouches. The optional stopcock may be used for intermittent patency check.

Mac-Loc® Locking Loop Mechanism

To Lock

1. Stabilize the Mac-Loc catheter hub assembly with one hand and pull back on the monofilament to form the distal catheter loop configuration (**Fig. 1**).
2. While maintaining traction on the monofilament, push the locking cam lever down until a distinct “snap” is felt. The catheter’s distal loop is now locked into position (**Fig. 2**).

3. Trim off the excess monofilament (**Fig. 3**).

To Unlock

1. While stabilizing the Mac-Loc catheter hub assembly with one hand, position the end of a small, sterile straight or curved forceps into the Mac-Loc release notch.
2. Pry upward until the locking cam lever is free (**Fig. 4**).
3. It is now possible to remove the catheter.

Catheter Exchange

1. Before unlocking the Mac-Loc catheter hub assembly, advance the distal end of a wire guide into the locked loop configuration of the catheter.
2. Release the Mac-Loc as directed above.
3. Advance the wire guide through the catheter end hole. Catheter exchange may now be performed in a standard fashion.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2020
Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd