COOK

MEDICAL

CE
0088

EN
3

Approach® PRO LT, Approach® CTO, Approach® PRO ST, and Approach® Hydro ST Micro Wire Guides

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PT
14

Microfios guia Approach® PRO LT, Approach® CTO, Approach® PRO ST e Approach® Hydro ST

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SV
15

Approach® PRO LT, Approach® CTO, Approach® PRO ST och Approach® Hydro ST mikroledare

Bruksanvisning

ZH
16

Approach® PRO LT、Approach® CTO、Approach® PRO ST 和 Approach® Hydro ST 微细导丝

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T _ C M W _ R E V 6

APPROACH® PRO LT, APPROACH® CTO, APPROACH® PRO ST AND APPROACH® HYDRO ST MICRO WIRE GUIDES

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

Approach PRO LT, Approach CTO, Approach PRO ST, and Approach Hydro ST Micro Wire Guides are 0.014 inches (0.36mm) in diameter and are available in a variety of lengths and configurations. The Approach PRO LT is also available in 0.018 inches (.46mm) in diameter and is available in a variety of lengths and configurations. Refer to product label for product specifications (wire guide length, wire guide diameter, shaft configuration and tip length).

Approach PRO LT incorporates a unique nitinol shaft resulting in exceptional torque control and steerability. The tip, made of stainless steel and palladium, provides excellent visibility. The nitinol shaft is coated with polytetrafluoroethylene (PTFE).

Approach CTO incorporates a stainless steel shaft along with distal stainless steel and platinum coil tips of variable stiffness, used for general peripheral vascular interventions.

Approach PRO ST incorporates a stainless steel, polytetrafluoroethylene (PTFE) coated shaft, along with distal stainless steel and palladium coils, used for general peripheral vascular interventions.

Approach Hydro ST incorporates a stainless steel shaft with distal stainless steel and platinum coils and a hydrophilic floppy tip, used for general peripheral vascular interventions.

INTENDED USE

Approach PRO LT, Approach CTO, Approach PRO ST, and Approach Hydro ST Wire Guides are intended for use in facilitating delivery of percutaneous catheters into the peripheral vasculature.

CONTRAINDICATIONS

None known

WARNINGS

- This wire guide is a delicate instrument and should be handled carefully; forceful angulation should be avoided.
- The wire guide should be advanced only when visualizing the tip of the wire guide fluoroscopically. The wire guide should not be torqued without evidence of corresponding movement of the distal tip.
- Refer to the labeling of the percutaneous catheter or other therapeutic devices for contraindications and potential complications associated with use.
- Approach CTO series wire guides have stiff distal ends. Operate these wire guides carefully to avoid injury to the blood vessel, observing information in these instructions.
- Approach Hydro ST wire guides should never be used in conjunction with or through a metal needle cannula.

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of vascular access sheaths, angiographic catheters and wire guides should be employed.
- Avoid manipulating or withdrawing the wire guide back through a metal needle or cannula. A sharp edge may scrape the coating or shear the wire guide. A catheter, introducer sheath or vessel dilator should replace the needle as soon as the wire guide has been inserted in the vessel.
- Excessive tightening of the torque device (if provided) onto the wire guide may result in abrasion of the coating on the wire.
- Always flush the Hydro ST before use.
- Always maintain a wet surface on the hydrophilic portion of the Hydro ST for optimal results.

INSTRUCTIONS FOR USE

1. Flush the wire guide holder by attaching or pressing a syringe with heparinized saline or sterile water to the fitting of the wire guide holder and injecting enough solution to wet the wire guide surface entirely.
2. Carefully remove the wire guide from the holder.
3. If desired, wire guide tip may be carefully shaped using standard tip shaping techniques. Do not use a shaping instrument with a sharp edge.
4. Insert a wire guide insertion tool through the valve assembly of the guiding sheath or guiding catheter.
5. Carefully insert the distal tip of the wire guide through the insertion tool.
6. Attach the torque device (if provided) to the wire guide.
7. Under fluoroscopy and while maintaining position of the peripheral vascular access catheter, advance the wire guide to the targeted site.

NOTE: Under fluoroscopy, observe all wire guide movement in the vessels. The wire guide should not be torqued without evidence of corresponding movement of the distal tip. Further torquing against resistance may cause vessel trauma or wire guide damage which may lead to device fracture. The wire guide should be advanced only when visualizing the tip of the wire guide fluoroscopically.

NOTE: If resistance is noted tactilely or visually under fluoroscopy, determine the cause and take action necessary to relieve the resistance. Advancement and withdrawal of the wire guide should be performed slowly and carefully.

8. With the wire guide secured in place, advance the percutaneous catheter or other therapeutic device to the target site.

NOTE: Always maintain a continuous flush while exchanging wire guides to prevent air entry into the system.

NOTE: Prior to re-insertion of the exchanged wire guide, verify that the distal tip of the cardiovascular device is free within the vessel. If the tip is against the wall, vessel trauma may occur upon re-introduction of the wire guide.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.