**COOK**


MEDICAL

**CE**  
0123

<b>EN</b> 3	<b>Coda® Balloon Catheter</b> Instructions for Use
<b>BG</b> 6	<b>Балонен катетър Coda®</b> Инструкции за употреба
<b>CS</b> 10	<b>Balónkový katetr Coda®</b> Návod k použití
<b>DA</b> 14	<b>Coda® ballonkateter</b> Brugsanvisning
<b>DE</b> 17	<b>Coda® Ballonkatheter</b> Gebrauchsanweisung
<b>EL</b> 21	<b>Καθετήρας με μπαλόνι Coda®</b> Οδηγίες χρήσης
<b>ES</b> 25	<b>Catéter balón Coda®</b> Instrucciones de uso
<b>FR</b> 29	<b>Cathéter à ballonnet Coda®</b> Mode d'emploi
<b>HU</b> 33	<b>Coda® ballonkatéter</b> Használati utasítás
<b>IT</b> 37	<b>Catetere a palloncino Coda®</b> Istruzioni per l'uso
<b>LT</b> 40	<b>„Coda®“ balioninis kateteris</b> Naudojimo instrukcija
<b>NL</b> 44	<b>Coda®-ballonkatheter</b> Gebruiksaanwijzing
<b>NO</b> 48	<b>Coda® ballongkateter</b> Bruksanvisning
<b>PL</b> 52	<b>Cewnik balonowy Coda®</b> Instrukcja użycia
<b>PT</b> 56	<b>Cateter de balão Coda®</b> Instruções de utilização
<b>SK</b> 59	<b>Balónikový katéter Coda®</b> Návod na použitie
<b>SL</b> 63	<b>Balonski kateter Coda®</b> Navodila za uporabo
<b>SV</b> 67	<b>Coda® ballongkateter</b> Bruksanvisning

T \_ C O D A 4 6 - M \_ R E V 3

## CODA® BALLOON CATHETERS

**Read all instructions carefully. Failure to properly follow the information provided may lead to the device not performing as intended or injury to the patient.**

**CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).**

### DEVICE DESCRIPTION

The Coda Balloon Catheters each consist of two independent lumens. The "Distal" lumen extends the length of the catheter and is used for placement over wire guides. The "Balloon" lumen is used to inflate and deflate the balloon. The balloon is manufactured from a compliant polyurethane material. Particular care should be taken in handling the balloon to prevent damage. The balloon will inflate to the indicated size parameters when utilizing proper volume recommendations. Radiopaque bands are placed on the balloon catheter to assist with positioning of the device under fluoroscopy.

### Performance Characteristics

- Marker Bands – Placed on the balloon catheter to assist with positioning of the device under fluoroscopy.
- Hub – A connection to facilitate the attachment of a standard Luer Syringe.
- Tapered Distal Tip – minimize vessel trauma and to facilitate over-the-wire trackability.
- 10 Fr Balloon – maximum inflation diameter is 46 mm, with a corresponding maximum balloon inflation volume of 60 cc, respectively.

### Device Compatibility

- The catheter is compatible with .035 inch wire guides.
- For introduction, it is recommended to use a minimum 14 Fr introducer sheath for the 10 Fr Coda Balloon.

### Patient Population

The target population for the Coda Balloon Catheter is based on the need for temporary occlusion of large vessels, or expansion of vascular prostheses, and on the physician's discretion on whether this occlusion balloon catheter is appropriate for a given adult patient population.

### Intended User

Coda Balloon Catheter should only be used by licensed and trained physicians with experience in diagnostic and interventional techniques.

### Contact with Body Tissue

The device is tissue contacting in line with intended use.

### Operating Principle

The Coda Balloon Catheter is a 10 Fr device measuring 140 cm in usable length and contains two independent lumens within a single extrusion. The distal lumen extends the length of the catheter and is used for placement over a 0.035 inch diameter wire guide. The balloon (inflation and deflation) lumen extends from the proximal hub to the two inflation media exit ports within the balloon. Two radiopaque bands are located within the balloon bonds to assist with positioning of the balloon under fluoroscopy. The catheter is designed with a tapered distal tip to minimize vessel trauma and to facilitate over the wire trackability. The Coda Balloon Catheter is recommended for use with a 14 Fr introducer sheath. A one-way low pressure stopcock is provided with the Coda Balloon Catheter.

### INTENDED USE

The intended use is to temporarily occlude large blood vessels to control hemorrhage from vessels distal to the balloon's location or to expand and mold the sealing stents of endovascular stent-grafts to the vessel wall.

### INDICATIONS FOR USE

- Severe hemorrhage or prevention of severe hemorrhage from vessel(s) that may be minimized by large vessel balloon occlusion proximal to the hemorrhage.
- Deployed endovascular stent-grafts that require balloon expansion and molding of the sealing stents to the large vessel wall.

### CLINICAL BENEFITS

- To maintain or prevent worsening of hemodynamic status
- To reduce the risk of endoleak and rupture of an aneurysmal sac

### CONTRAINDICATIONS

For Resuscitative Endovascular Balloon Occlusion of the Aorta (REBOA), when:

- Suspected thoracic aortic injury
- Cardiac injury

- Other vascular injury above the diaphragm
- Injuries above the balloon site

## WARNINGS

- Do not exceed maximum inflation volume. Adhere to balloon inflation parameters as shown below. Over-inflation of balloon may result in:
  - Damage to vessel wall and/or vessel rupture.
  - Rupture of balloon.
- Do not use a pressure inflation device for balloon inflation.
- Do not use a power injector for injection of contrast medium through distal lumen. Rupture may occur.
- The Coda 46 mm balloon catheter should not be used for dilation of vascular prostheses in iliac or other non-aortic vessels. Injury to vessel wall and/or rupture may occur.
- The Coda 46 mm balloon catheter should not be used in vessels less than 24 mm in diameter.
- When used to expand a vascular prosthesis, the balloon radiopaque markers should remain within the prosthesis.
- Not for use as a valvuloplasty balloon catheter.
- The single use device is not designed for re-use. Attempts to reprocess (re-sterilize) and/or to re-use may lead to chemical contamination, device failure, and/or transmission of disease.
- Do not use the device if the sterile packaging is damaged or unintentionally opened before use.

### Maximum Inflation Volumes

Catheter Size	Max. Volume
CODA-2-10.0-35-140-46	60 cc

## PRECAUTIONS

- This product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of vascular access sheaths, angiographic catheters and wire guides should be employed.
- Prolonged duration of occlusion may cause paralysis or damage to vital organs.
- The balloon is constructed of heat-sensitive material. Do not heat or attempt to shape the catheter tip.
- Always manipulate catheter using fluoroscopic control.
- Use only the recommended balloon inflation medium. Never use air or any gaseous medium to inflate balloon.
- Do not use after labeled expiration date.
- Always monitor balloon inflation using fluoroscopic control.

## POTENTIAL ADVERSE EVENTS

- Access site events (arteriovenous fistula, bleeding, dissection, embolism, extremity ischemia, hematoma, infection, lymphatic injury/leak, nerve injury, pain, perforation, pseudoaneurysm, spasm, thrombosis)
- Amputation, limb
- Arrhythmia
- Balloon catheter events (insertion failure, migration, malposition, rupture)
- Bleeding/hemorrhage
- Cardiac adverse events (electrocardiogram changes, elevated pulmonary artery pressure, heart failure, increased afterload, infarction, ischemia, ventricular dilation)
- Cerebrovascular (i.e., stroke) and neurological events and/or injury
- Compartment syndrome (abdominal, limb)
- Conversion to open repair or open occlusion
- Death
- Embolism (air, atheromatous plaque, debris, thrombus), arterial or venous
- Endoleak(s)
- Hypertension (arterial or venous) and/or hyperperfusion
- Hypotension and/or shock
- Infection, sepsis, septic shock
- Inflammation/Inflammatory response
- Organ and/or tissue ischemia leading to death, hepatic failure, intestinal ischemia, limb loss, multiple organ failure, myonecrosis, renal failure, reperfusion injury, or spinal cord injury
- Pain
- Pulmonary adverse events (acute respiratory distress syndrome, edema, infection, respiratory failure)
- Stent-graft events (damage, erosion, fracture, migration, thrombosis)
- Thrombosis, and/or thromboembolism (arterial or venous)
- Vessel injury (arteriovenous fistula, artery spasm, dissection, intimal damage, perforation, pseudoaneurysm, rupture)

## HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Keep dry and out of sunlight. Upon removal from package, inspect the product to ensure no damage has occurred.

## INSPECTION OF DEVICE

Visually inspect the device thoroughly including all levels of the packaging to verify that there is no damage prior to use. Visually inspect and confirm that the integrity of the sterile barrier has not been compromised in any way.

## INSTRUCTIONS FOR USE

### Balloon Preparation

**NOTE:** "Balloon" lumen is for inflation and deflation of balloon.  
"Distal" lumen is for wire guide insertion and withdrawal.

The balloon and balloon lumen of the Coda Balloon Catheter contains air that must be removed before insertion. The required preparation is as follows:

1. Remove protective balloon sleeve.
2. Prepare balloon lumen with standard 3:1 saline and contrast mixture as follows:
  - a. Attach syringe, with appropriate amount of 3:1 saline and contrast mixture, to stopcock on balloon lumen.
  - b. Purge all air from balloon in standard fashion.
  - c. Completely deflate balloon and close stopcock.


### Balloon Introduction and Inflation

1. Flush the Distal lumen using heparinized saline solution.
2. Advance the balloon catheter over a pre-positioned .035 inch wire guide, using a minimum 14 Fr introducer sheath for a 10 Fr Coda Balloon Catheter.

**NOTE:** If resistance is met while advancing the wire guide or balloon catheter, determine the cause and proceed with caution.

**CAUTION:** Prior to introduction, determine the amount of standard 3:1 saline and contrast mixture needed to inflate the balloon to the desired inflation diameter. Refer to the Balloon Inflation Parameters chart. Over-inflation of the balloon may result in damage to vessel wall and/or vessel rupture.

**Coda  
Balloon Catheter Inflation Parameters**

◆ = CODA-2-10.0-35-[140]-46

3. Under fluoroscopy, advance the balloon to the desired position using radiopaque markers.

**CAUTION:** If the Coda Balloon Catheter is being utilized to expand a vascular prosthesis, use the radiopaque markers to ensure that the entire balloon is positioned within the prosthesis.

4. Inflate balloon with standard 3:1 saline and contrast mixture using a 20 cc or larger syringe. Adhere to recommended balloon inflation volumes.
5. If balloon pressure is lost and/or balloon rupture occurs, deflate the balloon and remove balloon and sheath as a unit.

**NOTE:** Care should be taken to monitor balloon manipulations and inflation using fluoroscopy at all times.

6. Prepare balloon lumen with standard 3:1 saline and contrast mixture, as follows:
  - a. Attach syringe, with appropriate amount of 3:1 saline and contrast mixture, to stopcock on balloon lumen.
  - b. Purge all air from balloon in standard fashion.

- c. Completely deflate balloon and close stopcock.

### **Balloon Deflation and Withdrawal**

1. Completely deflate the balloon using an appropriately sized syringe. Allow adequate time for the balloon to deflate.
2. Deflate the balloon by pulling vacuum on the inflation syringe.
3. Maintain vacuum on the balloon and withdraw the catheter. If resistance is met during the withdrawal, apply negative pressure with a larger syringe before proceeding. If resistance continues, remove balloon and sheath as a unit.

### **DISPOSAL OF DEVICES**

After the procedure, this device may be contaminated with potentially infectious substances of human origin and should be disposed of in accordance with institutional guidelines.

### **REFERENCES**

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

### **PATIENT COUNSELING INFORMATION**

Please inform the patient as necessary of the relevant warnings, precautions, contraindications, measures to be taken, and limitations of use that the patient should be aware of. You can use the basic unique device identifier for this device – (0827002CINC012-262 BHV) – to find the Summary of Safety and Clinical Performance (SSCP) on the Eudamed website. When the Eudamed is available, use the following link: ([ec.europa.eu/tools/eudamed](http://ec.europa.eu/tools/eudamed)).

### **SERIOUS INCIDENT REPORTING**

If any serious incident has occurred in relation to the device, this should be reported to Cook Medical and this shall also be reported to the competent authority the device was used in.

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© COOK 2023
**Cook Medical Europe Ltd**
**William A. Cook Australia Pty Ltd**