COOK

MEDICAL

CE
0088

EN 2	Coda® and Coda® LP Balloon Catheters Instructions for Use
DA 4	Coda® og Coda® LP ballonkatetre Brugsanvisning
DE 7	Coda® und Coda® LP Ballonkatheter Gebrauchsanweisung
EL 10	Καθετήρες με μπαλόνι Coda® και Coda® LP Οδηγίες χρήσης
ES 13	Catéteres balón Coda® y Coda® LP Instrucciones de uso
FR 16	Cathéters à ballonnet Coda® et Coda® LP Mode d'emploi
HU 19	Coda® és Coda® LP ballonkatéterek Használati utasítás
IT 22	Cateteri a palloncino Coda® e Coda® LP Istruzioni per l'uso
NL 24	Coda® en Coda® LP ballonkatheters Gebruiksaanwijzing
NO 27	Coda® og Coda® LP ballongkatetere Bruksanvisning
PT 30	Cateteres de balão Coda® e Coda® LP Instruções de utilização
SV 33	Coda® och Coda® LP ballongkatetrar Bruksanvisning

T _ C O D A 4 6 _ R E V O

CODA® AND CODA® LP BALLOON CATHETERS

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Coda and Coda LP Balloon Catheters each consist of two independent lumens. The "Distal" lumen extends the length of the catheter and is used for placement over wire guides. The "Balloon" lumen is used to inflate and deflate the balloon.

The balloon is manufactured from a compliant polyurethane material. Particular care should be taken in handling the balloon to prevent damage. The balloon will inflate to the indicated size parameters when utilizing proper volume recommendations.

Radiopaque bands are placed on the balloon catheter to assist with positioning of the device under fluoroscopy.

INTENDED USE

The Coda and Coda LP Balloon Catheters are intended for temporary occlusion of large vessels, or to expand vascular prostheses.

CONTRAINDICATIONS

None known

WARNINGS

- Do not exceed maximum inflation volume. Adhere to balloon inflation parameters as shown in **Figs. 1 and 2**. Over-inflation of balloon may result in:
 - Damage to vessel wall and/or vessel rupture.
 - Rupture of balloon.
- Do not use a pressure inflation device for balloon inflation.
- Do not use a power injector for injection of contrast medium through Distal lumen. Rupture may occur.
- **The Coda 46 mm balloon catheter should not be used for dilation of vascular prostheses in iliac or other non-aortic vessels. Injury to vessel wall and/or rupture may occur.**
- **The Coda 46 mm balloon catheter should not be used in vessels less than 24 mm in diameter.**
- When used to expand a vascular prosthesis, the balloon radiopaque markers should remain within the prosthesis.
- Not for use as a valvuloplasty balloon catheter.

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of vascular access sheaths, angiographic catheters and wire guides should be employed.
- Prolonged duration of occlusion may cause paralysis or damage to vital organs.
- The balloon is constructed of heat-sensitive material. Do not heat or attempt to shape the catheter tip.
- Always manipulate catheter using fluoroscopic control.
- Use only the recommended balloon inflation medium. Never use air or any gaseous medium to inflate balloon.
- Do not use after labeled expiration date.
- Always monitor balloon inflation using fluoroscopic control.

POTENTIAL ADVERSE EVENTS

- Vessel dissection, perforation, rupture or injury
- Occlusion at some locations may cause arrhythmia
- Drug reactions
- Infection and pain at insertion site
- Arterial thrombosis and/or embolism
- Death

PRODUCT RECOMMENDATIONS

Wire Guide Selection

The catheter is compatible with .035 inch wire guides.

Introducer Sheath Selection

For introduction, it is recommended to use a minimum 14.0 French introducer sheath for the 10.0 French Coda Balloon Catheter and a minimum 12.0 French introducer sheath for the 9.0 French Coda LP Balloon Catheter.

Balloon Inflation Volume

Do not exceed maximum inflation volume. Adhere to balloon inflation parameters as shown in **Figs. 1 and 2**. Over-inflation of balloon may result in:

- Damage to vessel wall and/or vessel rupture.
- Rupture of balloon.

Maximum Inflation Volumes

Catheter Size	Max. Volume
CODA-2-10.0-35-120-46	60 cc
CODA-2-10.0-35-140-46	60 cc
CODA-2-9.0-35-100-32	30 cc
CODA-2-9.0-35-120-32	30 cc


Fig. 1


Fig. 2

INSTRUCTIONS FOR USE


Fig. 3

NOTE: "Balloon" lumen is for inflation and deflation of balloon.

NOTE: "Distal" lumen is for wire guide insertion and withdrawal.

Balloon Preparation

NOTE: Balloon and balloon lumen of the Coda and Coda LP Balloon Catheters contain air. The air must be removed from balloon and balloon catheter prior to insertion using standard technique.

1. Remove protective balloon sleeve.
2. Prepare balloon lumen with standard 3:1 saline and contrast mixture as follows:
 - a. Attach syringe, with appropriate amount of 3:1 saline and contrast mixture, to stopcock on balloon lumen.
 - b. Purge all air from balloon in standard fashion.
 - c. Completely deflate balloon and close stopcock.
3. To increase ease of insertion, balloon may be lubricated with a thin layer of sterile, biocompatible lubricant.

Balloon Introduction and Inflation

1. Flush the Distal lumen using heparinized saline solution.
2. Advance the balloon catheter over a pre-positioned .035 inch wire guide, **utilizing a minimum 14.0 French introducer sheath for the 10.0 French Coda Balloon Catheter or a minimum 12.0 French introducer sheath for the 9.0 French Coda LP Balloon Catheter.**

NOTE: If resistance is met while advancing the wire guide or balloon catheter, determine the cause and proceed with caution.

CAUTION: Prior to introduction, determine the amount of standard 3:1 saline and contrast mixture needed to inflate the balloon to the desired inflation diameter. Refer to the Balloon Inflation Parameters chart in Figs. 1 and 2. Over-inflation of the balloon may result in damage to vessel wall and/or vessel rupture.

3. Under fluoroscopy, advance the balloon to the desired position using radiopaque markers.

CAUTION: If the Coda or Coda LP Balloon Catheter is being utilized to expand a vascular prosthesis, use the radiopaque markers to ensure that the entire balloon is positioned within the prosthesis.
4. Inflate balloon with standard 3:1 saline and contrast mixture using a 20 cc or larger syringe. **Adhere to recommended balloon inflation volumes.**
5. If balloon pressure is lost and/or balloon rupture occurs, deflate the balloon and remove balloon and sheath as a unit.

NOTE: Care should be taken to monitor balloon manipulations and inflation using fluoroscopy at all times.

Balloon Deflation and Withdrawal

1. Completely deflate the balloon using an appropriately sized syringe. **Allow adequate time for the balloon to deflate.**
2. Deflate the balloon by pulling vacuum on the inflation syringe.
3. Maintain vacuum on the balloon and withdraw the catheter. If resistance is met during withdrawal, apply negative pressure with a larger syringe before proceeding. If resistance continues, remove balloon and sheath as a unit.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2012
Bloomington, IN 47404 U.S.A.