COOK

MEDICAL

CE
0088

EN
3

TriForce™ Peripheral Crossing Set

Instructions for Use

DA
4

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T _ C T F _ R E V 1

TRIFORCE™ PERIPHERAL CROSSING SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The TriForce Peripheral Crossing Set combines the CXI Support Catheter and Flexor® technology in one device. The set includes three components: An outer 5.0 French Flexor sheath with radiopaque tip and hydrophilic coating, an inner 4.0 French CXI Support Catheter, and a Peel-Away® sheath. The outer Flexor sheath is tapered down to a 4.0 French inner diameter at the tip to ensure smooth transition to the CXI Support Catheter. This set is compatible with a 0.035" wire guide.

CXI Support Catheter Flow Rates (ml/sec)*

The infusion dynamic pressure should not exceed 1200 psi (82.7 bar). The static pressure should not exceed 700 psi (48.3 bar).

CXI Model	Sterile Saline			Contrast**		
	100 psi	300 psi	600 psi	100 psi	300 psi	600 psi
65 cm, straight/ angled	5.9	11.3	16.2	3.4	8.3	13.2
100 cm, straight/ angled	4.9	9.3	13.1	2.6	6.3	10.3

*Catheter flow rates are approximate values.

**100% Omnipaque. Omnipaque is a registered trademark of Amersham Health, Princeton, New Jersey.

INTENDED USE

The TriForce Peripheral Crossing Set is intended to be percutaneously introduced into blood vessels and support a wire guide while performing percutaneous peripheral interventions. This device is also intended for injection of radiopaque contrast media for the purpose of angiography.

CONTRAINDICATIONS

None known

WARNINGS

- The maximum recommended infusion pressure is **1200 psi**.
- Avoid introducing air or any other gas through the catheter into the vascular system.

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in percutaneous peripheral interventions. Standard techniques for placement of vascular access sheaths, angiographic catheters, and wire guides should be employed.
- The device is designed and intended for one time use only. **Do not re-sterilize and/or reuse.**
- Device manipulation should only occur under fluoroscopy.
- If the device is used for infusion, reference the table of flow rates in the Device Description section of these Instructions for Use and ensure infusion pressure does not exceed the recommendations.
- The device should not be advanced into a vessel having a reference vessel diameter smaller than the sheath outer diameter.
- The device should not be forced through an area of resistance unless the source of resistance is identified by fluoroscopy and appropriate steps are taken to reduce or remove the obstruction.
- Due to its tapered design, the outer sheath of this set can only accommodate an inner device with an outer diameter of 4.0 French or smaller.

POTENTIAL ADVERSE EVENTS

Vascular catheterization and/or vascular interventions may result in complications including, but not limited to:

- Vessel dissection
- Perforation
- Total occlusion
- Embolism
- Vessel spasm
- Death

INSTRUCTIONS FOR USE

1. Flush the Flexor sheath through its side arm with heparinized saline solution or sterile water.

2. Flush the CXI Support Catheter with heparinized saline solution or sterile water.
3. Wet the entire surface of both the sheath and the support catheter with heparinized saline solution or sterile water. This will activate the hydrophilic coating, making the surface of both devices lubricious.
4. Preassemble the set as follows:
 - a. Insert the Peel-Away sheath into the Flexor sheath's Check-Flo valve.
 - b. Gently advance the CXI Support Catheter through the Peel-Away sheath into the Flexor sheath. **CAUTION:** Advancing the support catheter into a sheath without a wire guide can easily cause kinking of the catheter. Use extra caution during this process.
 - c. When the CXI Support Catheter is advanced completely into the Flexor sheath, peel away and discard the Peel-Away sheath.
5. Under fluoroscopic guidance, introduce the preassembled device into the vascular system over a previously placed wire guide.
6. Under fluoroscopic guidance, carefully advance the device to the desired vascular location. Final positioning is accomplished by stepwise, short advances of the wire guide, CXI Support Catheter, and Flexor sheath.

CAUTION: Be sure the wire guide tip is extended beyond the CXI Support Catheter tip at all times, and the CXI Support Catheter tip is extended beyond the Flexor sheath tip at all times.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2011
Bloomington, IN 47404 U.S.A.