EN
2

Approach® CTO Micro Wire Guide
Instructions for Use

CS
3

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Approach® CTO-mikroledare
Bruksanvisning

APPROACH® CTO MICRO WIRE GUIDE

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

The Approach CTO Micro Wire Guide is 0.014 inches (0.36 mm) in diameter and is available in a variety of lengths and tip configurations. The product has a TFE-coated stainless steel shaft and distal stainless steel and platinum coils on the tip. The tips' configurations differ in stiffness.

INTENDED USE

The Approach CTO Microwire Guide is intended for use in facilitating delivery of percutaneous catheters into the peripheral vasculature and is also indicated for the intra-luminal placement of percutaneous catheters or other therapeutic devices beyond stenotic lesions (including chronic total occlusions) in the peripheral vasculature prior to further percutaneous intervention.

CONTRAINDICATIONS

This device is not intended for coronary and neurovascular use.

WARNINGS

- This product is a delicate instrument. Avoid forceful angulation.
- Avoid manipulating or withdrawing the wire guide back through a metal needle or cannula. A sharp edge may scrape or shear material from the wire guide.
- Altering the tip's configuration or curve manually may damage the wire guide.
- Excessive tightening of a torque device may abrade the coating on the wire guide.
- The Approach CTO Micro Wire Guide has stiff distal tip configurations. Operate the wire guide carefully, according to the information in these instructions, to avoid injuring the blood vessel.

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of devices in angiographic procedures should be employed.
- Use medical imaging when you manipulate the wire guide. Do not advance or torque the wire guide without visual evidence of the corresponding movement of the distal tip.
- When you use the wire guide with another device, consider the end-hole size and the length of the device in order to ensure a proper fit between the wire guide and the device.
- Gram tip loads associated with the products in the Approach CTO series are measured with the distal tip of the wire guide exposed 1 cm beyond the distal tip of a catheter.

INSTRUCTIONS FOR USE

1. Flush the wire guide holder by attaching or pressing a syringe with heparinized saline or sterile water to the fitting or opening of the wire guide holder. Inject enough solution to wet the entire wire guide surface.
2. Carefully remove the wire guide from the holder.
3. Insert a wire guide insertion tool through the valve assembly of the guiding sheath or guiding catheter. Insert the distal tip of the wire guide through the insertion tool.
4. Attach a torque device (if provided) to the wire guide.
5. Advance the wire guide while maintaining the position of the percutaneous catheter or other therapeutic device.
NOTE: Torqueing or advancing against resistance may cause vessel trauma or device damage that can lead to device fracture. If resistance is noted tactilely or visually under fluoroscopy, determine the cause and relieve the resistance. Advance and withdraw the wire guide slowly and carefully.
6. Once the wire guide has crossed the target lesion, advance the percutaneous catheter or other therapeutic device to the target site.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd