
COOK**MEDICAL****CE**
0088


EN 4	Turbo-Ject® Peripherally Inserted Central Venous Catheters With Micropuncture® Peel-Away Introducers Instructions for Use
DA 11	Turbo-Ject® perifert indførte centralvenekatetre med Micropuncture® Peel-Away indførere Brugsanvisning
DE 19	Turbo-Ject® peripher eingeführte zentrale Venenkatheter mit Micropuncture® Peel-Away-Einführungsschleusen Gebrauchsanweisung
EL 27	Περιφερικά εισαγόμενοι κεντρικοί φλεβικοί καθετήρες Turbo-Ject® με εισαγωγείς Micropuncture® Peel-Away Οδηγίες χρήσης
ES 36	Catéteres venosos centrales de inserción periférica Turbo-Ject® con introductores Micropuncture® Peel-Away Instrucciones de uso
FR 44	Cathéters veineux centraux insérés en périphérie Turbo-Ject® avec introducteurs Micropuncture® Peel-Away Mode d'emploi
HU 53	Turbo-Ject® perifériásan bevezetett centrális vénás katéterek Micropuncture® Peel-Away bevezetőeszközökkel Használati utasítás
IT 61	Cateteri venosi centrali ad inserimento periferico Turbo-Ject® con introduttori Micropuncture® Peel-Away Istruzioni per l'uso
NL 69	Turbo-Ject® perifeer ingebrachte centraal-veneuze katheters met Micropuncture® Peel-Away introducers Gebruiksaanwijzing
NO 78	Turbo-Ject® sentralvenekatetre til perifer innføring med Micropuncture® Peel-Away innføringsenheter Bruksanvisning
PT 86	Cateteres venosos centrais de introdução periférica Turbo-Ject® com introdutores Peel-Away Micropuncture® Instruções de utilização
SV 94	Turbo-Ject® perifert inlagda centrala venkatetrar med Micropuncture® Peel-Away-införare Bruksanvisning


T _ C T P I C C T T _ R E V 1

1

2

3

TURBO-JECT® PERIPHERALLY INSERTED CENTRAL VENOUS CATHETERS WITH MICROPUNCTURE® PEEL-AWAY INTRODUCERS

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Turbo-Ject Peripherally Inserted Central Venous Catheter (PICC) is a polyurethane catheter with a polyurethane suture wing, reinforced polyurethane extensions, and plastic clamps. Catheters can be trimmed to fit the patient's anatomy. PICC sets and trays also contain a Peel-Away® introducer, appropriately sized access needle, wire guides, injection caps, and other accessories for percutaneous vascular placement.

Single Lumen Catheters

Catheter Size (Fr)	Outer Diameter		Inner Diameter		
	Inch	Gage	Inch	Gage	Lumen Vol. (mL)
4	0.055	17	0.037	17	0.76
5	0.066	15	0.048	17	1.00

Double Lumen Catheters

Catheter		Outer Diameter		Inner Diameter		
Size (Fr)	Hub	Inch	Gage	Inch	Gage	Lumen Vol. (mL)
4	Purple (CT)	0.055	17	0.022/0.038*	22	0.61
	White			0.014/0.028*	24	0.46
5	Purple (CT)	0.066	15	0.024/0.048*	18	0.72
	White (CT)					

*height/width

Catheter		Maximum Flow Rate*	Injection Pressure Limit Setting (psi)
Size (Fr)/ Lumen	Hub		
4 / Single		5 mL/sec	325
5 / Single		7 mL/sec	325
4 / Double	Purple (CT)	3 mL/sec	325
	White	Do Not Power Inject This Lumen	
5 / Double	Purple (CT)	5 mL/sec	325
	White (CT)		

*Flow rates achieved using room temperature Omnipaque 300® contrast and verified using a Medrad Stellant® CT injector system. Omnipaque 300 has a viscosity of 11.8 centipoise at room temperature (20 degrees C). A change in temperature or viscosity of the contrast medium used will result in a change in achievable flow rates.

Omnipaque 300® is a registered trademark of GE Healthcare.

INTENDED USE

Turbo-Ject Peripherally Inserted Central Venous Catheter (PICC) Sets and Trays are intended for short- or long-term use for venous pressure monitoring, blood sampling, administration of drugs and fluids, and for use with power injectors for delivery of contrast in CT studies. The Turbo-Ject PICC is indicated for multiple injections of contrast media through a power injector. The maximum pressure limit setting for Power Injectors used with the Turbo-Ject PICC may not exceed 325 psi and the flow rate may not exceed the maximum flow rate indicated.

CONTRAINDICATIONS

None known

WARNINGS

- **Peripherally Inserted Central Venous Catheters play an important role in treatment of critically ill patients. However, catheter tips can erode or perforate vascular walls. Extreme caution must be used in placement and monitoring of catheters.**
- Catheter tip position should be monitored by X-ray on a routine basis. Periodic lateral view X-ray is suggested to assess tip location in relation to vessel wall. Tip position should appear to be parallel to vessel wall. (**Reference 1**)
- The safe and effective use of Turbo-Ject PICC lines with power injector pressures set above 325 psi has not been established.
- **Do not power inject if maximum injection rate cannot be verified to meet limit printed on catheter hub or extension tube.**
- To safely use Turbo-Ject PICC lines with a power injector, the technician/ health care professional must verify:

- The catheter lumen has “CT” on the hub to indicate the lumen is power injectable.
- Prior to use that the maximum pressure limit is set at or below 325 psi and that the maximum flow rate is at or below that which is listed on the catheter. Dynamic and static pressure test results are shown in the following table.

Turbo-Ject Power Injectable PICC Dynamic and Static Pressure Results

Catheter		Priming Volume (mL)	Maximum Labeled Flow Rate (mL/sec)	Average Max Catheter Pressure During Maximum Flow (psi)*	Average Static Burst Pressure in 37°C Water (psi)**	Range of Static Burst Pressure in 37°C Water (psi)**
Size (Fr)/ Lumen	Hub					
4 / Single		0.76	5	204	296	292-310
5 / Single		1.00	7	125	258	250-266
4 / Double	Purple (CT)	0.61	3	178	261	258-264
	White	0.46	Do Not Power Inject This Lumen			
5 / Double	Purple (CT)	0.72	5	191	286	277-300
	White (CT)					

*Maximum flow rate pressures are determined with pump safety cut-off set at 325 psi, using contrast media with a viscosity of 11.8 cP.

**Static Burst pressure is the failure point of the catheter when totally occluded. Warning: Power injector machine pressure limiting feature may not prevent over pressurization of an occluded catheter.

PRECAUTIONS

- This product is intended for use by health care practitioners trained and experienced in proper positioning of catheters in the central venous system using percutaneous entry (Seldinger) technique. Standard techniques for placement of vascular access sheaths, catheters and wire guides should be employed.
- Select puncture site and length of catheter needed by assessing patient anatomy and condition.
- If lumen flow is impeded, do not force injection or withdrawal of fluids. Notify attending physician immediately.
- Patient movement can cause catheter tip displacement. Catheters placed

via an antecubital vein have shown tip movement of up to 10 cm with motion of the extremity.

- Catheter size should be as small as the use will allow.

PRODUCT RECOMMENDATIONS

Catheter Size and Puncture Site


Preliminary reports indicate that catheter size can influence clotting. Larger diameter catheters have more tendency to promote clots. As reported by Amplatz, Gianturco and others, (Reference 2) clot formation has less relation to type of catheter material than to size of catheter.

The angle of the catheter tip to the vessel wall should be checked carefully. Blackshear reviewed the medical literature of catheter perforations, which have confirming X-rays, and found that an incident angle of the catheter to vessel wall greater than 40 degrees was more likely to perforate. (Reference 3)

The following variables must also be considered in selecting appropriate catheter and length:

1. Patient history
2. Patient age and size
3. Access site available
4. Unusual anatomical variables
5. Proposed use and duration of treatment plan

Catheter Tip Positioning

Verify catheter tip position using radiography or appropriate technology. In order to guarantee extrapericardial location, catheter tip should be located in the lower 1/3 of the SVC and not beyond the extrapericardial cava. **Every effort must be made to ascertain proper tip position in order to prevent erosion or perforation of the central venous system and to ensure proper delivery of infusates.**

Catheter Maintenance

Catheter entry site must be prepared and maintained in a manner consistent with standard procedure for central venous catheterization. After catheter placement and prior to use, tip position and lumen patency should be confirmed by free aspiration of venous blood. **If blood is not freely aspirated, catheter tip position should be immediately reevaluated by physician.** If catheter is not to be used immediately, its lumen should be maintained by continuous saline or heparinized saline drip or locked with a catheter locking solution. **NOTE:** If MicroClave or other needleless adapters approved for saline only lock are used, saline only catheter lock may be used. Catheter heparinization should be determined by institutional protocol and clinical judgement. Heparin concentrations of 10 Units/mL to 100 Units/mL have been reported adequate to maintain lumen patency. **(Reference 4)** Catheter lock should be reestablished after every use or at least every 24 hours if unused. Before using catheter lumen already locked with heparin, lumen should be flushed with twice the indicated lumen volume using normal saline. Lumen should be flushed with normal saline between administration of different infusates. After use, lumen should again be flushed with twice the indicated lumen volume using normal saline before reestablishing catheter lock. Strict aseptic technique must be adhered to while using and maintaining catheter.

INSTRUCTIONS FOR USE

Catheter Obturator Preparation

1. Flush catheter with heparinized saline solution or sterile water. **NOTE:** The catheter may be trimmed if a shorter length is required.
2. Attach a syringe with heparinized saline solution or sterile water to the Luer lock fitting of the catheter obturator holder.
3. Inject enough solution to wet the obturator surface entirely. This will activate the AQ® hydrophilic coating, making the obturator surface very lubricious. **NOTE:** If the surface of the obturator becomes dry after removal from the holder, wetting with additional heparinized saline or sterile water will renew the hydrophilic effect.
4. The catheter/obturator assembly can now be introduced as described in the following section.

Catheter Placement (Fluoroscopic Method)

5. After prepping the access site, introduce the access needle into the vessel. **NOTE:** The use of ultrasound is helpful to determine suitability for vessel access and patency. The EchoTip® marking on the needle is used to help visualize the tip of the needle during vessel access.
6. Using fluoroscopic guidance, introduce the wire guide through the needle and advance it 15-20 cm into the vessel.
7. Withdraw the needle, leaving wire guide in place. If necessary, enlarge the puncture site with scalpel blade.

8. Introduce the Peel-Away introducer assembly (sheath and dilator) over the wire guide. With a twisting motion, advance the assembly into the vessel. (**Fig. 1**)
9. Using fluoroscopic control, determine the correct catheter length by advancing the wire guide to the desired catheter tip location. Once the wire guide tip is in proper position, mark the length by clamping forceps onto the wire guide at the skin site.
10. Withdraw the wire guide and measure it from the forceps to the distal tip to determine correct catheter length. Trim catheter if necessary. **NOTE:** Withdraw obturator prior to trimming the catheter; reinsert for catheter introduction.
11. Leaving the sheath in place, remove the dilator. (**Fig. 2**) **NOTE:** To prevent inadvertent air aspiration after removing wire guide and dilator, place thumb or finger over the cuffed proximal end of the sheath.
12. Introduce the catheter/obturator assembly into the sheath as far as possible. **NOTE:** Resistance may be felt approximately 1-2 cm distal to the catheter suture wing when introducing the assembly into the sheath due to an increase in outer diameter. (**Fig. 3**)
13. Peel the sheath away from the catheter by grasping the two wings of the sheath and pulling outward and upward. **NOTE:** Final catheter position is accomplished by alternately advancing the catheter into the sheath and then further pulling on the two wings. Once the sheath is removed, a slight advancement of the catheter may be needed for final positioning.
14. After the catheter is in final position, remove the obturator, secure the catheter to skin and dress in standard fashion.
15. Verify catheter tip position using radiography or appropriate technology. In order to guarantee extrapericardial location, catheter tip should be located above the SVC-RA junction, within the lower 1/3 of the SVC.

Catheter Placement (Non-Fluoroscopic Method)

1. After prepping the access site, introduce the access needle into the vessel. **NOTE:** The use of ultrasound is helpful to determine suitability for vessel access and patency. The EchoTip® marking on the needle is used to help visualize the tip of the needle during vessel access.
2. Introduce the wire guide through the needle and advance it 15-20 cm into the vessel.
3. Withdraw the needle, leaving wire guide in place. If necessary, enlarge the puncture site with scalpel blade.
4. Introduce the Peel-Away introducer assembly (sheath and dilator) over the wire guide. With a twisting motion, advance the assembly into the vessel. (**Fig. 1**)
5. Using a tape measure, clinical judgment, or other institutional protocol, determine the correct catheter length and trim catheter as needed. **NOTE:** Withdraw obturator prior to trimming catheter; reinsert for catheter introduction.
6. Leaving the sheath in place, remove the dilator. (**Fig. 2**) **NOTE:** To prevent inadvertent air aspiration after removing wire guide and dilator, place thumb or finger over the proximal end of the sheath.
7. Introduce the catheter/obturator assembly into the sheath as far as possible. **NOTE:** Resistance may be felt approximately 1-2 cm distal to the

catheter suture wing when introducing the assembly into the sheath due to an increase in outer diameter. (**Fig. 3**)

8. Peel the sheath away from the catheter by grasping the two wings of the sheath and pulling outward and upward. **NOTE:** Final catheter position is accomplished by alternately advancing the catheter into the sheath and then further pulling on the two wings. Once the sheath is removed, a slight advancement of the catheter may be needed for final positioning.
9. After catheter is in final position, remove obturator, secure catheter to skin and dress in standard fashion.
10. Verify catheter tip position using radiography or appropriate technology. In order to guarantee extrapericardial location, catheter tip should be located above the SVC-RA junction, within the lower 1/3 of the SVC.

Power Injection Procedure

1. Use only lumens marked "CT" for power injection of contrast media.
WARNING: Use of lumens not marked "CT" for power injection may result in catheter failure.
2. Confirm proper catheter tip position radiographically prior to injection.
3. Remove any injection/needleless caps from the Turbo-Ject PICC.
4. Attach a 10 mL (or larger) syringe filled with sterile normal saline to the hub of the extension tube to be used for power injection.
5. Ensure adequate blood return and flush catheter thoroughly with the entire 10 mL of sterile normal saline to ensure lumen patency.
WARNING: Failure to ensure patency of the catheter lumen prior to injection may result in catheter failure.
6. Remove syringe and attach power injection device to the catheter using the manufacturer's recommendations.
7. Conduct study using the power injector, making sure not to exceed the maximum flow rate or pressure limit for the catheter.
8. Disconnect the power injection device and flush the catheter again with 10 mL of sterile normal saline.
9. Place a new injection/needleless cap on the Turbo-Ject PICC, flush and lock the catheter with saline or heparinized saline per institutional protocol.
10. Confirm proper catheter tip position radiographically following power injection.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

1. Tocino IM, Watanabe A: "Impending Catheter Perforation of Superior Vena Cava: Radiographic Recognition," *American Journal of Roentgenology*, 1986; 146:487-490.

2. Amplatz K. "A Simple Non-thrombogenic Coating." *Investigative Radiology*, 1971; 6:280-289.
3. Welch RH, Gravenstein N, Blackshear RH: "Multilumen Central Venous Catheters in Children: Relative Potential to Perforate Vessels. An In Vitro Study," *Journal of Clinical Monitoring*, 1997; 13:75-59.
4. Kelly C, Dumenko L, McGregor SE, McHutchion ME: "Change in Flushing Protocols of Central Venous Catheters," *Oncology Nursing Forum*, 1992; 19:599-605.

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1. Tocino IM, Watanabe A: "Impending Catheter Perforation of Superior Vena Cava: Radiographic Recognition," *American Journal of Roentgenology*, 1986; 146:487-490.
2. Amplatz K. "A Simple Non-thrombogenic Coating." *Investigative Radiology*, 1971; 6:280-289.
3. Welch RH, Gravenstein N, Blackshear RH: "Multilumen Central Venous Catheters in Children: Relative Potential to Perforate Vessels. An In Vitro Study," *Journal of Clinical Monitoring*, 1997; 13:75-59.

4. Kelly C, Dumenko L, McGregor SE, McHutchion ME: "Change in Flushing Protocols of Central Venous Catheters," *Oncology Nursing Forum*, 1992; 19:599-605.

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1. Tocino IM, Watanabe A: Impending Catheter Perforation of Superior Vena Cava: Radiographic Recognition, *American Journal of Roentgenology*, 1986; 146:487-490.
2. Amplatz K. A Simple Non-thrombogenic Coating. *Investigative Radiology*, 1971; 6:280-289.
3. Welch RH, Gravenstein N, Blackshear RH: Multilumen Central Venous Catheters in Children: Relative Potential to Perforate Vessels. An In Vitro Study, *Journal of Clinical Monitoring*, 1997; 13:75-59.

4. Kelly C, Dumenko L, McGregor SE, McHutchion ME: Change in Flushing Protocols of Central Venous Catheters, *Oncology Nursing Forum*, 1992; 19:599-605.

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1. Tocino IM, Watanabe A: "Impending Catheter Perforation of Superior Vena Cava: Radiographic Recognition," *American Journal of Roentgenology*, 1986; 146:487-490.
2. Amplatz K. "A Simple Non-thrombogenic Coating." *Investigative Radiology*, 1971; 6:280-289.
3. Welch RH, Gravenstein N, Blackshear RH: "Multilumen Central Venous Catheters in Children: Relative Potential to Perforate Vessels. An In Vitro Study," *Journal of Clinical Monitoring*, 1997; 13:75-59.
4. Kelly C, Dumenko L, McGregor SE, McHutchion ME: "Change in Flushing Protocols of Central Venous Catheters," *Oncology Nursing Forum*, 1992; 19:599-605.

MANUFACTURER

COOK INCORPORATED
750 Daniels Way
Bloomington, IN 47404 U.S.A.

EC REPRESENTATIVE

WILLIAM COOK EUROPE ApS
Sandet 6, DK-4632
Bjaeverskov, DENMARK

www.cookmedical.com

© COOK 2013

2013-10

T_CTPICCTT_REV1