COOK

MEDICAL

CE
0123

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T _ C T R U W C _ R E V 1

CONE TIP URETERAL CATHETER & RUTNER UNIVERSAL WEDGE CATHETER

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner). Please read all instructions before using this device.

DEVICE DESCRIPTION

The Cone Tip Ureteral Catheter is available in 4.8 and 6.0 French. The Rutner Universal Wedge Catheter is available in 4.0 and 5.0 French. The distal tips of both catheters are radiopaque and designed to occlude the ureteral orifice.

INTENDED USE

The Cone Tip Ureteral Catheter and Rutner Universal Wedge Catheter are intended for retrograde pyelogram. The 4.0 French Rutner Universal Wedge Catheter is indicated for pediatric population 2 years and older.

CONTRAINDICATIONS

These devices are contraindicated for use in patients unsuitable for catheterization.

WARNINGS

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use.

PRECAUTIONS

- This device is intended for use by physicians trained and experienced in urological techniques. Standard techniques should be employed.
- Avoid bending or kinking the catheter prior to placement. Doing so could damage the integrity of the catheter and result in patient injury.
- If you encounter resistance while advancing or withdrawing the catheter, stop. Determine the cause of the resistance before proceeding.
- The potential effects of phthalates on children or pregnant or nursing women have not been fully characterized and there may be concern for reproductive and developmental effects.

POTENTIAL ADVERSE EVENTS

Adverse events that may result from ureteral catheter placement include, but are not limited to:

- Extravasation
- Hemorrhage
- Sepsis
- Edema
- Loss of renal function
- Perforation or laceration of the kidney, renal pelvis, ureter, and/or bladder
- Peritonitis
- Urinary tract infection
- Abrasion

INSTRUCTIONS FOR USE

1. Place the cystoscope within the bladder to locate the ureteral orifice.
2. Advance the catheter through the cystoscope and into the bladder.
3. Under direct vision, advance the catheter tip to the ureteral orifice. The catheter tip will occlude the ureteral orifice.

NOTE: Depending upon the French size of the catheter, the catheter may need to be backloaded into the cystoscope before advancing.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

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Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd