EN
3

CrossCath™ Support Catheter

Instructions for Use

DA
4

CrossCath™ støttekateter

Brugsanvisning

DE
6

CrossCath™ Support-Katheter

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EL
8

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T _ C X C _ R E V 2

CROSSCATH™ SUPPORT CATHETER

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DESCRIPTION

CrossCath Support Catheters with hydrophilic coating are single-lumen intravascular catheters, designed to support a wire guide during access of vasculature, allow for exchange of wire guides, and provide a conduit for the delivery of saline solutions or diagnostic contrast agents.

CrossCath Support Catheters are available in nine different configurations as shown in the table below. All models have 3 radiopaque markers spaced equally along the distal shaft to aid in estimating lengths within the vascular system. The distal radiopaque marker is positioned approximately 3 mm from the catheter tip. Each model has a smaller distal portion for passing through smaller vasculature and a larger proximal portion for additional support.

Catheter Flow Rates (ml/sec)*

The infusion dynamic pressure should not exceed 300 psi. The static pressure should not exceed 150 psi.

Model	Size (in)	Length (cm)	Coating Length (cm)	Sterile Saline		Contrast**	
				150 psi	300 psi	150 psi	300 psi
CXC3.0-1.9-14-65-P-NS-0	0.014	65	40	1.5	2.2	0.6	0.9
CXC3.0-1.9-14-90-P-NS-0	0.014	90	40	1.3	2.1	0.3	0.7
CXC3.0-1.9-14-135-P-NS-0	0.014	135	40	1.2	1.9	0.3	0.6
CXC3.0-1.9-14-150-P-NS-0	0.014	150	40	1.2	1.9	0.2	0.6
CXC3.4-2.2-18-65-P-NS-0	0.018	65	40	2.5	3.5	1.0	1.7
CXC3.4-2.2-18-90-P-NS-0	0.018	90	40	1.9	3.2	0.5	1.4
CXC3.4-2.2-18-135-P-NS-0	0.018	135	40	1.6	2.6	0.3	1.0
CXC3.4-2.2-18-150-P-NS-0	0.018	150	40	1.6	2.5	0.3	1.0
CXC4.8-3.7-35-65-P-NS-0	0.035	65	40	9.4	13.8	5.1	8.6
CXC4.8-3.7-35-90-P-NS-0	0.035	90	60	8.5	13.1	4.7	8.7
CXC4.8-3.7-35-135-P-NS-0	0.035	135	60	6.9	10.9	3.4	7.3
CXC4.8-3.7-35-150-P-NS-0	0.035	150	60	6.2	10.3	3.0	7.3

* Catheter flow rates are approximate values.

**100% Omnipaque. Omnipaque is a registered trademark of Amersham Health, Princeton, New Jersey.

INTENDED USE

The CrossCath Support Catheter is intended for use in small vessel or superselective anatomy for diagnostic and interventional procedures, including peripheral use.

CONTRAINDICATIONS

None known

WARNINGS

- The maximum recommended infusion pressure is **300 psi**.
- Avoid introducing air or any other gas through the catheter into the vascular system.

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in small vessel access and interventional procedures. Standard techniques for placement of percutaneous catheters should be employed.
- If the catheter is used for infusion, reference the Catheter Flow Rates table and ensure infusion pressure does not exceed the recommendations.
- The catheter is designed and intended for one time use only. **Do not re-sterilize and/or reuse.**

- Catheter manipulation should only occur under fluoroscopy.
- The catheter should not be advanced into a vessel with a diameter smaller than the catheter outer diameter.
- The catheter should not be advanced through an area of resistance unless the source of resistance is identified by fluoroscopy and appropriate steps are taken to reduce or remove the obstruction.
- Do not advance catheter without the wire extending out of the distal tip.
- Only use wire guides of the recommended diameter and length.
- This catheter should only be used by physicians qualified to perform percutaneous vascular interventions.

POTENTIAL ADVERSE EVENTS

Vascular catheterization and/or vascular interventions may result in complications including, but not limited to:

- Death
- Embolism
- Perforation
- Total occlusion
- Vessel dissection
- Vessel spasm

INSTRUCTIONS FOR USE

1. Remove the CrossCath spiral holder from its sealed packaging.
2. Activate hydrophilic coating by wetting the catheter with heparinized saline solution or sterile water. To ensure hydrophilic activation, wet the entire surface of the catheter.
NOTE: The surface of the catheter may become dry if not used immediately after activation. Additional wetting with heparinized saline solution or sterile water will reactivate the hydrophilic coating.
3. Flush the catheter with heparinized saline solution.
4. Through a previously introduced, appropriately sized guiding catheter or sheath, insert the CrossCath over an appropriately sized wire guide using standard technique.
5. Under fluoroscopic guidance, carefully advance the catheter to the desired location, making sure that the wire guide tip extends beyond the catheter tip at all times. Final positioning is accomplished by short advances of the wire guide and catheter until the desired position is achieved and then confirmed by fluoroscopic visualization.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2011
Bloomington, IN 47404 U.S.A.