<b>EN</b> 3	<b>CXI® Support Catheter</b> Instructions for Use
<b>CS</b> 4	<b>Podpůrný katetr CXI®</b> Návod k použití
<b>DA</b> 6	<b>CXI®-støttekateter</b> Brugsanvisning
<b>DE</b> 8	<b>CXI® Support-Katheter</b> Gebrauchsanweisung
<b>EL</b> 10	<b>Καθετήρας στήριξης CXI®</b> Οδηγίες χρήσης
<b>ES</b> 12	<b>Catéter de soporte CXI®</b> Instrucciones de uso
<b>FR</b> 14	<b>Cathéter support CXI®</b> Mode d'emploi
<b>HU</b> 16	<b>CXI® támaszkatéter</b> Használati utasítás
<b>IT</b> 18	<b>Catetere di supporto CXI®</b> Istruzioni per l'uso
<b>NL</b> 20	<b>CXI® ondersteuningskatheter</b> Gebruiksaanwijzing
<b>NO</b> 22	<b>CXI® støttekateter</b> Bruksanvisning
<b>PL</b> 24	<b>Cewnik podpierający CXI®</b> Instrukcja użycia
<b>PT</b> 26	<b>Cateter de suporte CXI®</b> Instruções de utilização
<b>SV</b> 28	<b>CXI®-stödkateter</b> Bruksanvisning

## CXI® SUPPORT CATHETER

**Caution: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).**

### DEVICE DESCRIPTION

The CXI Support Catheter with hydrophilic coating is a braided, kink-resistant catheter designed to facilitate wire guide exchange, infusion, and wire guide support.

The catheter has 4 radiopaque markers spaced 5 cm apart, starting at the distal tip.

### Catheter Flow Rates (mL/sec)\*

2.3 French Model	Sterile Saline			Contrast**		
	100 psi	300 psi	600 psi	100 psi	300 psi	600 psi
65 cm, straight/angled	0.6	1.1	1.6	0.1	0.4	0.9
90 cm, straight/angled	0.5	0.9	1.4	0.1	0.4	0.7
135 cm, straight/angled	0.5	0.9	1.3	0.1	0.2	0.5
150 cm, straight/angled	0.4	0.7	1.1	0.1	0.2	0.4

2.6 French Model	Sterile Saline			Contrast**		
	100 psi	300 psi	600 psi	100 psi	300 psi	600 psi
65 cm, straight/angled	1.3	2.4	3.3	0.4	1.1	1.7
90 cm, straight/angled	1.1	2.1	3.0	0.3	0.9	1.7
135 cm, straight/angled	1.0	1.6	2.4	0.2	0.6	1.1
150 cm, straight/angled	0.9	1.6	2.3	0.2	0.5	1.0

4.0 French Model	Sterile Saline			Contrast**		
	100 psi	300 psi	600 psi	100 psi	300 psi	600 psi
65 cm, † straight/angled	6.7	12.2	25.0	4.7	10.9	14.2
90 cm, straight/angled	6.0	11.2	23.0	3.7	8.4	12.1
135 cm, straight/angled	4.7	8.6	17.5	2.6	6.2	9.7
150 cm, straight/angled	4.6	8.5	17.5	2.1	4.9	8.7

\*Catheter flow rates are approximate values

\*\*100% Omnipaque. Omnipaque is a registered trademark of Amersham Health, Princeton, New Jersey

† Values for 65 cm device are extrapolated from other device data

### INTENDED USE

The CXI Support Catheter is intended for use in small vessel or superselective anatomy for diagnostic and interventional procedures, including peripheral use.

## CONTRAINDICATIONS

None known

## WARNINGS

- Avoid introducing air or any other gas through the catheter into the vascular system.

## PRECAUTIONS

- This product is intended for use by physicians trained and experienced in small vessel access and interventional procedures. Standard techniques for placement of percutaneous catheters should be employed.
- Catheter manipulation should only occur under fluoroscopy.
- If the catheter is used for infusion, reference the table of flow rates in the Device Description section of these Instructions for Use and ensure infusion pressure does not exceed the recommendations.
- The catheter should not be advanced into a vessel having a reference vessel diameter smaller than the catheter outer diameter.
- This catheter is designed and intended for one time use only. **Do not re-sterilize and/or reuse.**
- The catheter should not be advanced through an area of resistance unless the source of resistance is identified by fluoroscopy and appropriate steps are taken to reduce or remove the obstruction.
- This catheter should only be used by physicians qualified to perform percutaneous vascular interventions.

## POTENTIAL ADVERSE EVENTS

Vascular catheterization and/or vascular interventions may result in complications including, but not limited to:

- Vessel dissection
- Perforation
- Total occlusion
- Embolism
- Vessel spasm
- Death

## INSTRUCTIONS FOR USE

1. Remove the catheter spiral holder from its sealed packaging.
2. Activate hydrophilic coating by wetting the catheter with heparinized saline solution or sterile water. To ensure hydrophilic activation, wet the entire surface of the catheter.

**NOTE:** The surface of the catheter may become dry if not used immediately after activation. Additional wetting with heparinized saline solution or sterile water will reactivate the hydrophilic coating.

3. Under fluoroscopic guidance, introduce the catheter into the vascular system using standard techniques.

## HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

## REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd