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OSTEO-SITE® VERTEBROPLASTY AND BONE ACCESS NEEDLE SETS

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Osteo-Site Vertebroplasty and Bone Access Needle Sets consist of a stainless steel outer cannula with an accompanying inner stylet and obturator. Sets with the suffix -S in the order number include a stainless steel biopsy needle. A plastic, radiolucent needle holder is also available to aid in needle placement.

INTENDED USE

The Osteo-Site Vertebroplasty and Bone Access Needle Sets are intended for vertebral body access, biopsy and infusion during a vertebroplasty procedure. The products are intended for use by physicians trained and experienced in vertebroplasty and percutaneous bone access techniques. Standard techniques for vertebral body access, biopsy, and infusion should be employed.

CONTRAINDICATIONS

None known

WARNINGS

None known

PRECAUTIONS

Illustrations

- Assure that inner stylet is in proper alignment with outer cannula for appropriate bevel orientation. (**Fig.1**)
- Assure all needles are at the desired location within the vertebral body before cement mixing begins.

INSTRUCTIONS FOR USE

NOTE: The handling and curing characteristics of bone cement vary with temperature, humidity, mixing technique, injection pressure, and procedure time. Physicians should be familiar with the properties, handling characteristics, and application of bone cement before the procedure begins.

1. Under fluoroscopic guidance using both anterior/posterior and lateral views, advance the outer cannula (with inner stylet in place) to the desired location within the vertebral body.
2. Remove the inner stylet from the outer cannula.
3. Perform biopsy or infusion.

NOTE: Upon completion of infusion procedure, the inner stylet may be re-inserted into the outer cannula to remove remaining infusion material from the lumen of the cannula.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.