COOK

MEDICAL

- EN** **Pakter Curved Needle Set**
Instructions For Use
- 2**
- DA** **Pakter buet nålesæt**
Brugsanvisning
- 3**
- DE** **Pakter gebogenes Kanülenset**
Gebrauchsanweisung
- 4**
- EL** **ΣΕΤ κυρτών βελόνων Pakter**
Οδηγίες χρήσης
- 6**
- ES** **Equipo de aguja curvada Pakter**
Instrucciones de uso
- 7**
- FR** **Set d'aiguille courbe Pakter**
Mode d'emploi
- 9**
- IT** **Set con ago curvo Pakter**
Istruzioni per l'uso
- 10**
- NL** **Pakter gebogen naaldset**
Gebruiksaanwijzing
- 12**
- NO** **Pakter buet nålesett**
Bruksanvisning
- 13**
- PT** **Conjunto de agulha curva Pakter**
Instruções de utilização
- 14**
- SV** **Pakter böjt nålset**
Bruksanvisning
- 16**

T _ D C H N S _ R E V 2

PAKTER CURVED NEEDLE SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Pakter Curved Needle Set (prefix DCHNS-) is a two-needle coaxial set consisting of a 21 gage stainless steel outer needle and stylet and a 25 gage curved nitinol needle and stylet.

INTENDED USE

The Pakter Curved Needle Set is intended for discography, diagnostic sampling, aspiration and injection.

CONTRAINDICATIONS

None known

WARNINGS

None known

PRECAUTIONS

- Use of imaging technology (fluoroscopy, CT, ultrasound) is recommended when using this product.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

INSTRUCTIONS FOR USE

1. Use imaging technology or anatomic landmarks to determine the external percutaneous entry site.
2. Use standard facility protocol to prepare the targeted entry site.
3. Remove the Pakter Curved Needle Set assembly from the sterile packaging.
4. If desired, use a scalpel to make a skin nick at the external percutaneous entry site.
5. Insert the assembled 21 gage needle cannula and stylet from the Pakter Curved Needle Set through the skin.
6. Using imaging guidance, advance the cannula/stylet assembly to the desired depth, adjacent to the target anatomic location.
7. Remove the stylet from the 21 gage cannula, leaving the cannula in position.
8. Insert the 25 gage curved needle and stylet into the 21 gage cannula, advancing the needle until the distal curved segment is entirely within the cannula.
9. Orient the notch on the hub of the 25 gage curved needle in the direction in which you want the curved needle to track.
10. Advance the 25 gage curved needle slowly. Adjust the orientation of the curved needle as necessary to enter the target anatomic location.
11. Once the desired anatomic position is acquired, remove the stylet from the 25 gage curved needle.
12. You may inject contrast media or aspirate fluid through the 25 gage curved needle to further confirm desired tip location in the target anatomy.

- Once desired location is confirmed, proceed with aspiration and/or injection of contrast media or other solutions as desired.
- Following completion of the procedure, maintain position of the 21 gauge outer needle and withdraw the 25 gauge curved needle through the 21 gauge cannula.
- With the 25 gauge curved needle removed, remove the 21 gauge cannula and apply a sterile dressing.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2022
Bloomington, IN 47404 U.S.A.