COOK

MEDICAL

CE
0088

EN
3

Deutsch Gastrostomy Catheter

Instructions for Use

DA
4

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Brugsanvisning

DE
5

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EL
6

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T _ D E U T S C H _ R E V 4

DEUTSCH GASTROSTOMY CATHETER

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Deutsch Gastrostomy Catheter is a 16 French Ultrathane® catheter with a Cook-Cope locking loop mechanism, flexible stiffener, enteral feeding adapter and trocar stylet for single-stick placement.

INTENDED USE

The Deutsch Gastrostomy Catheter is intended for gastroenteric feeding.

CONTRAINDICATIONS

The product should not be used in patients with coagulopathy.

WARNINGS

This product contains **NATURAL RUBBER LATEX**, which may cause allergic reactions.

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in gastrostomy techniques. Standard techniques for placement and maintenance of gastrostomy catheters should be employed.
- Manipulation of product requires ultrasound, fluoroscopy or other imaging guidance.
- Ensure that the catheter pigtail is appropriately located within the stomach before securing the catheter to the skin.
- A TFE-coated wire guide must be used with this catheter.
- Hold pigtail retention suture during insertion of stiffening cannula to avoid bunching or tangling of suture.
- Activate hydrophilic coating by wetting surface of device with sterile water or saline. For best results, maintain wetted condition of device during placement.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

INSTRUCTIONS FOR USE

Trocar Placement

1. Insufflate the stomach with air.
2. Following the Instructions for Use provided with the Cope Gastrointestinal Suture Anchor Set (sold separately), anchor the anterior wall of the stomach to the abdominal wall.
3. Using the trocar (pre-loaded with catheter), make a puncture at the appropriate location in the upper third of the abdominal wall.
4. Under fluoroscopic guidance, advance the trocar until its tip is positioned within the fundus of the stomach.
5. Withdraw the trocar.
6. Pull suture end tight to form loop configuration in catheter and inject contrast to verify catheter tip position within stomach.
7. Wrap suture around catheter and secure with a knot, then trim off excess suture and slide sleeve over suture to prevent leakage.
8. Attach catheter to skin surface using physician preferred method of attachment.

Seldinger Placement

1. Insufflate the stomach with air.
2. Following the Instructions for Use provided with the Cope Gastrointestinal Suture Anchor Set (sold separately), anchor the anterior wall of the stomach to the abdominal wall.
3. With an 18 gage introducer needle, puncture the stomach at the appropriate location in the upper third of the abdominal wall.
NOTE: The tract through the abdominal wall may be dilated if necessary.
4. Under fluoroscopic guidance, insert a heavy-duty wire guide into the fundus of the stomach.
5. Slide the catheter over the wire guide until it is within the stomach.
6. Remove the wire guide.
7. Pull suture end tight to form loop configuration in catheter and inject contrast to verify catheter tip position within stomach.
8. Wrap suture around catheter and secure with a knot, then trim off excess suture and slide sleeve over suture to prevent leakage.
9. Attach catheter to skin surface using physician preferred method of attachment.

Unlocking Catheter Loop

1. Advance wire guide into catheter.
2. Uncover suture by folding back the latex sleeve, and release suture.
3. Withdraw catheter over wire guide.

NOTE: For catheter exchange, insert new catheter over wire guide. For removal, withdraw wire guide following removal of catheter.

POST-PLACEMENT CARE

- It is not uncommon for a small amount of intraperitoneal air to be present post-procedure.

- Feedings may be instituted immediately, although it is recommended to begin with dilute feedings and advance gradually. **NOTE:** If signs of peritonitis develop, feedings should be discontinued and appropriate antibiotics instituted. Feedings should not be resumed until all parameters are normal.
- If using a Luer tip syringe for flushing purposes, remove feeding adapter and connect Luer converter to the hub of the catheter. Remove Luer converter and connect the feeding adapter after flushing the catheter.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2011
Bloomington, IN 47404 U.S.A.