- EN**
3 **Double Flexible Tipped Wire Guides**
Instructions for Use
- CS**
4 **Vodicí dráty s ohebnými hroty na obou koncích**
Návod k použití
- DA**
5 **Kateterledere med dobbelt fleksibel spids**
Brugsanvisning
- DE**
6 **Führungsdrähte mit beidseits biegsamer Spitze**
Gebrauchsanweisung
- EL**
7 **Συρμάτινοι οδηγοί διπλού εύκαμπτου άκρου**
Οδηγίες χρήσης
- ES**
8 **Guías con las dos puntas flexibles**
Instrucciones de uso
- FR**
9 **Guides à double extrémité souple**
Mode d'emploi
- HU**
10 **Kettős hajlékony végű vezetődrótok**
Használati utasítás
- IT**
11 **Guide con doppia punta flessibile**
Istruzioni per l'uso
- NL**
12 **Voerdraden met twee flexibele tippen**
Gebruiksaanwijzing
- NO**
14 **Ledevaiere med dobbel fleksibel spiss**
Bruksanvisning
- PL**
15 **Prowadniki z dwoma giętkimi końcówkami**
Instrukcja użycia
- PT**
16 **Fios guia com ponta dupla flexível**
Instruções de utilização
- SV**
17 **Ledare med dubbla, böjliga spetsar**
Bruksanvisning

DOUBLE FLEXIBLE TIPPED WIRE GUIDES

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

Double Flexible Tipped Wire Guides are made with a stainless steel inner mandril and stainless steel outer coil. They come in a variety of diameters, lengths, and tip configurations. These wires are provided with two usable ends. Some devices are provided with a green PTFE coating to ease coaxial use with other devices. Refer to the product label for product specifications.

INTENDED USE

Double Flexible Tipped Wire Guides are used to facilitate the placement of devices during diagnostic and interventional procedures.

CONTRAINDICATIONS

Not intended for use in the coronary or neurovascular arteries.

WARNINGS

- This product is a delicate instrument. Avoid forceful angulation.
- Avoid manipulating or withdrawing the wire guide back through a metal needle or cannula. A sharp edge may scrape or shear material from the wire guide.
- Altering the tip's configuration or curve manually may damage the wire guide.
- Excessive tightening of a torque device may abrade the coating on the wire guide.

PRECAUTIONS

- This product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of devices in angiographic procedures should be employed.
- Do not attempt to torque the wire guide.
- Use medical imaging when you manipulate the wire guide. Do not advance or manipulate the wire guide without visual evidence of the corresponding movement of the distal tip.
- When you use the wire guide with another device, consider the end-hole size and the length of the device in order to ensure a proper fit between the wire guide and the device.
- These wire guides come with two usable ends. Either end may be used, depending on the clinical situation as determined by an experienced operator.
- Do not advance or withdraw a wire guide when resistance is encountered, as a perforation could occur.
- Inspect the wire guide for kinks or damage prior to use.

POTENTIAL ADVERSE EVENTS

Using this product may cause the following adverse events.

- Allergic reactions
- Thrombosis and embolization
- Vessel dissection or perforation

INSTRUCTIONS FOR USE

1. Flush the wire guide holder by attaching a syringe with heparinized saline or sterile water to the fitting of the wire guide holder. Inject enough solution to wet the wire guide surface entirely.

NOTE: If flushing through the wire guide holder is not possible, remove the wire guide from the holder and place it in a bowl of heparinized saline or sterile water, or wet the wire guide surface over the entire length by wiping once or twice with gauze moistened with heparinized saline solution.

2. Carefully remove the wire guide from the holder.
3. If needed, insert a wire guide insertion tool (provided) through the valve assembly or hub of the guiding sheath or other interventional device. Insert the tip of the wire guide through the insertion tool.
4. Standard wire guide techniques may now be employed.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for more information.

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© COOK 2019
Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd