COOK
MEDICAL

CE
0123

EN
3

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T - D L U A C - R E V 1

DUAL LUMEN URETERAL ACCESS CATHETER AND FLEXI-TIP® DUAL LUMEN URETERAL ACCESS CATHETER

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner). Please read all instructions before using this device.

DEVICE DESCRIPTION

The Dual Lumen Ureteral Access Catheter and Flexi-Tip Dual Lumen Ureteral Access Catheter are catheters with an independent .050 inch diameter wire guide lumen and a .050 inch diameter injection lumen.

NOTE: The length of the catheter has ink marks spaced 1 cm apart from each other which are 1 mm in width. In addition to this basic pattern, there are variations at increments of 5 cm. From the distal tip, one wide ink mark 3 mm in width denotes 5 and 30 cm, two ink marks spaced 1 mm apart denotes 10 and 35 cm, three ink marks spaced 1 mm apart denotes 15 and 40 cm, four ink marks spaced 1 mm apart denotes 20 and 45 cm, and five ink marks spaced 1 mm apart denotes 25 cm. The spacing of the ink marks is intended as reference for visual aid only.

INTENDED USE

This device is intended for injection of contrast medium and/or safety wire guide placement during percutaneous and transurethral access. The catheters that are 24 cm in length can be used in pediatric patients 12 years and over.

CONTRAINDICATIONS

This device is contraindicated for use in patients unsuitable for catheterization.

WARNINGS

Sterilized by ethylene oxide gas. Do not use if sterile barrier is damaged. The device is intended for one-time use.

PRECAUTIONS

- These devices are intended for use by physicians trained and experienced in urological techniques. Standard techniques should be employed.
- Avoid bending or kinking the catheter prior to placement. Doing so could damage the integrity of the catheter and result in patient injury.
- If you encounter resistance while advancing or withdrawing the catheter, stop. Determine the cause of the resistance before proceeding.
- Do not withdraw catheter while it is deflected in a cystoscope/endoscope.

POTENTIAL ADVERSE EVENTS

Adverse events that may result from ureteral catheter placement include, but are not limited to:

- Extravasation
- Hemorrhage
- Sepsis
- Edema

- Loss of renal function
- Perforation or laceration of the kidney, renal pelvis, ureter, and/or bladder
- Peritonitis
- Urinary tract infection
- Abrasion

INSTRUCTIONS FOR USE

1. Perform standard techniques for gaining access to the ureter utilizing percutaneous or transurethral access, leaving a wire guide in place.
2. Flush both channels of the dual lumen catheter. **NOTE:** If using an AQ®-coated Dual Lumen Access Catheter, activate the AQ coating by briefly immersing or wiping down the catheter in sterile water or saline.
3. Using direct, fluoroscopic or radiographic visualization advance the catheter over the previously placed wire guide to the desired location within the ureter. **NOTE:** If simultaneous pyelogram is needed, affix the syringe with the contrast solution to the blue lumen port (labeled "D") and inject. If placing a second wire, thread the wire through the blue lumen port (labeled "D") and to the desired location within the ureter.
4. To remove the catheter, maintain wire position with one hand and withdraw the catheter.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

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Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd