COOK

MEDICAL

CE
0088

EN 2	Hinck Spinal Needle with Depth Gauge Instructions For Use
DA 3	Hinck-spinalkanyle med dybdemåler Brugsanvisning
DE 4	Hinck Spinalkanüle mit Tiefenmarkierung Gebrauchsanweisung
EL 6	Ενδοραχιαία βελόνα Hinck με μετρητή βάθους Οδηγίες χρήσης
ES 7	Aguja espinal Hinck con medidor de profundidad Instrucciones de uso
FR 9	Aiguille rachidienne Hinck avec jauge de profondeur Mode d'emploi
IT 10	Ago spinale con indicatore di profondità Hinck Istruzioni per l'uso
NL 11	Hinck spinale naald met dieptemeter Gebruiksaanwijzing
NO 13	Hinck spinalnål med dybdemåler Bruksanvisning
PT 14	Agulha espinal Hinck com medidor de profundidade Instruções de utilização
SV 16	Hinck spinalnål med djupmätare Bruksanvisning

T _ D P N H _ R E V O

DISPOSABLE HINCK DEPTH GUIDE NEEDLE

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Disposable Hinck Depth Guide Needle (prefix DPNH-) consists of a lancet-bevel cannula, a matching beveled stylet and a blunt-tipped stylet. The blunt-tipped stylet is designed to extend 5 millimeters past the proximal tip of the lancet bevel when inserted through the outer cannula lumen. The proximal hub of the cannula is designed with a notch to indicate the directional positioning of the lancet bevel and to secure a fit to the proximal hub of the stylet.

INTENDED USE

The Disposable Hinck Depth Guide Needle is intended for diagnostic sampling of cerebrospinal fluid and for the introduction of contrast medium for myelography.

CONTRAINDICATIONS

None known

WARNINGS

None known

PRECAUTIONS

Use of imaging technology (fluoroscopy, CT, ultrasound) is recommended when using this product.

INSTRUCTIONS FOR USE

1. Use imaging technology or anatomic landmarks and palpation to determine the external percutaneous entry site.
2. Use standard facility protocol to prepare the targeted entry site.
3. Remove the Disposable Hinck Depth Guide Needle components from the sterile packaging.
4. If desired, use a scalpel to make a skin nick at the external percutaneous entry site.
5. Insert the assembled needle cannula and beveled stylet from the Disposable Hinck Depth Guide Needle through the skin and advance to the subarachnoid space.
6. Remove the beveled stylet.
7. Gently insert the blunt-tipped stylet into the needle cannula. The pin on the stylet hub should engage the slot on the cannula hub.
8. If the blunt-tipped stylet does not advance fully into the cannula, retract the cannula while maintaining position with the blunt-tipped stylet.
9. If the blunt-tipped stylet advances fully into the cannula, advance the cannula and stylet until the stylet stops on the floor of the spinal canal.
10. Maintaining position of the needle cannula, remove the blunt-tipped stylet.
11. Proceed with aspiration of spinal fluid and/or injection of contrast media or other solutions as desired.
12. Once the procedures are completed, withdraw the needle cannula and apply a sterile dressing to the puncture site.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

Hinck VC. Spinal Needle with Depth Guide: A New Design for Simplifying Subarachnoid Injections. *Am J Roentgenol* 127:354 1976

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© COOK 2013
Bloomington, IN 47404 U.S.A.