EN 4	Filiform Double Pigtail Stent Set Instructions for Use
CS 6	Souprava filiformního double pigtailového stentu Návod k použití
DA 8	Trådagtigt stentsæt med dobbelt grisehale Brugsanvisning
DE 11	Filiform-Doppelpigtail-Stentset Gebrauchsanweisung
EL 13	Σετ ενδοπρόσθεσης νηματοειδούς διπλού σπειροειδούς άκρου Οδηγίες χρήσης
ES 16	Equipo de stent filiforme de doble pigtail Instrucciones de uso
FR 19	Set d'endoprothèse double pigtail filiforme Mode d'emploi
HU 21	Filiform kettős pigtail sztentkészlet Használati utasítás
IT 24	Set con stent a doppio pigtail e guida filiforme Istruzioni per l'uso
NL 27	Dubbele pigtail stentset met draadvormige voerdraadtip Gebruiksaanwijzing
NO 29	Filiformt stentsett med dobbel pigtail Bruksanvisning
PL 32	Zestaw stentu nitkowatego z dwiema końcówkami pigtail Instrukcja użycia
PT 34	Conjunto de stent de duplo pigtail com fio guia com ponta filiforme Instruções de utilização
SV 37	Filiformstentset med dubbla pigtails Bruksanvisning


Fig. 1


Fig. 2


Fig. 3


Fig. 4


Fig. 5

FILIFORM DOUBLE PIGTAIL STENT SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

- Radiopaque double pigtail stent
- Stent positioner with lockable fitting
- Wire guide (*size and configuration may vary*)

NOTE: Set components may vary.

INTENDED USE

The Filiform Double Pigtail Stent is used for temporary internal drainage from the ureteropelvic junction to the bladder.

The wire guide is used to gain ureteral access, to establish a tract, and to assist in the placement, replacement, and exchange of medical devices during urological and gynecologic procedures.

CONTRAINDICATIONS

None known

WARNINGS

None known

PRECAUTIONS

- Manipulation of the wire guide requires appropriate imaging control. Use caution not to force or over manipulate the wire guide when gaining access.
- When using a wire guide through a metal cannula/needle, use caution as damage may occur to the outer coating.
- When exchanging or withdrawing an instrument over the wire guide, secure and maintain the wire guide in place under fluoroscopy in order to avoid unexpected wire guide displacement.
- The included wire guide is not intended for PTCA use.
- Hydrophilically coated wires are very slippery when wet. Always maintain control of the wire guide when manipulating it through any device.
- Complications of ureteral stent placement are documented. Use of this device should be based upon consideration of risk-benefit factors as they apply to your patient. Informed consent should be obtained to maximize patient compliance with follow-up procedures.
- It is imperative that the filiform flexible end of the wire guide be introduced into the ureter first.
- Do not force set components during placement, replacement, or removal. Carefully remove the set components if any resistance is encountered.
- The stent must not remain indwelling more than twelve months. If the patient's status permits, the stent may be replaced with a new stent.
- The included stent is not intended as a permanent indwelling device.
- A pregnant patient must be more closely monitored for possible stent encrustation due to calcium supplements.
- Improper handling can seriously weaken the stent. Acute bending or overstressing during placement may result in subsequent separation of the stent at the point of stress after a prolonged indwelling period. Angulation of the wire guide or stent should be avoided. Use of a 0-degree scope lens is recommended. Scopes larger than 21.0 French are suggested.

- Individual variations of interaction between stents and the urinary system are unpredictable.
- Periodic evaluation via cystoscopic, radiographic, or ultrasonic means is suggested. The stent must be replaced if encrustation hampers drainage.
- If problems occur using this device, please call your Cook Urological sales representative or contact our Customer Service department at the address/phone listed at www.cookmedical.com.

INSTRUCTIONS FOR USE

Wire Guide Preparation

1. Using aseptic technique, remove the wire from its outer packaging and place in the sterile field.
2. If using a hydrophilically coated wire guide, fill a syringe with sterile water or sterile saline solution and attach it to the flushing port on the wire guide holder. Inject enough solution to wet the wire guide surface entirely. This will activate the hydrophilic coating. **NOTE:** For optimal performance, rehydrate the hydrophilically coated wire guide after exposure to ambient environment or, after extended use, replace it with a new hydrophilically coated wire guide.
3. Remove the wire guide from the holder.

Assembly for Filiform Placement

1. Using a baseline pyelogram, estimate the length of stent needed by taking the ureteral length and adding 1 cm. Accurate measurement enhances drainage efficiency and patient comfort.
2. Introduce the flexible end of the wire guide into the fitting end of the stent positioner and advance it until the flexible end of the wire guide extends half its length beyond the blunt end of the positioner.
NOTE: Wire guides without SlipCoat hydrophilic coating should be well lubricated with a sterile water-soluble lubricant.
NOTE: The lockable fitting on the distal end of the stent positioner must be loosened to allow the wire guide to be introduced.
3. Carefully introduce the flexible end of the wire guide into the blunt end of the stent and advance it until the stent is straightened and the flexible tip of the wire guide extends beyond the tapered end of the stent.
4. While maintaining the position of the stent and the wire guide, advance the stent positioner until its blunt end rests against the blunt end of the stent.
5. Check to ensure the flexible tip of the wire guide extends beyond the tip of the straightened stent and then tighten the lockable fitting to secure the position of the stent positioner on the wire guide.

Endoscopic Placement

1. Advance the flexible end of the stent assembly, under cystoscopic vision, into the ureteral orifice. Advancement of the stent up the ureter is accomplished under fluoroscopic control.
2. When the flexible tip of the wire guide and the proximal end of the stent are fluoroscopically confirmed to be in the renal pelvis, stabilize the positioner, loosen the lockable fitting, and gently withdraw the wire guide to allow the retentive coil to form in the renal pelvis (**Fig. 1**).
3. Maintaining the position of the positioner, continue to gently withdraw the wire guide until the intravesicle coil forms.
4. Completely remove the positioner and the wire guide (**Fig. 2**).
5. Fluoroscopically confirm the position of the stent coils in the bladder and renal pelvis.
6. When necessary, final positioning of the stent can be performed with endoscopic forceps, or a 2-0 loop tether placed through the distal sideport prior to stent placement.

Intraoperative Placement

1. Insert the flexible end of the wire guide into one end of the stent and straighten the entire stent.
2. Grasping the stent and wire guide between the thumb and index finger, advance the stent through the incision in the ureter and into the bladder (**Fig. 3**).

3. Remove the wire guide.
4. Using syringe suction at the sideport, aspirate urine to confirm the stent's proper position in the bladder.
5. Maintaining the position of the stent, reinsert the flexible tip of the wire guide into a sideport near the incision site and advance the wire guide to straighten the second coil.
6. Grasping the stent and wire guide, advance the stent through the same incision and into the renal pelvis (**Fig. 4**).
7. Remove the wire guide, allowing the stent coil to re-form (**Fig. 5**).

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2012
Bloomington, IN 47404 U.S.A.
COOK IRELAND LTD.