COOK

MEDICAL

CE
0088

EN
3

Dotter Intravascular Retriever Set

Instructions for Use

DA
3

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T _ D R S _ R E V 3

DOTTER INTRAVASCULAR RETRIEVER SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Dotter Intravascular Retriever Set consists of a 4-wire helical loop basket (7 cm long, 3 cm wide) with a pin vise handle, a catheter and a Check-Flo® Introducer Set.

INTENDED USE

The Dotter Intravascular Retriever Set is intended to snare a foreign body and withdraw it to a peripheral vascular location.

CONTRAINDICATIONS

None known

WARNINGS

None known

PRECAUTIONS

- The product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of vascular access sheaths, angiographic catheters and wire guides should be employed.
- Manipulation of products requires fluoroscopic control.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

INSTRUCTIONS FOR USE

1. Insert the introducer sheath using standard percutaneous technique, then introduce the catheter into the desired position.
2. Slide the folding sleeve over the helical loop basket and introduce it into the catheter hub.
3. Advance the basket through the catheter until it emerges from the distal tip. **NOTE:** The basket will expand upon emergence from catheter.
4. Position the basket so that the foreign body is caught within it, then maintain basket position and slide the catheter forward to snare the foreign body.
5. Maintaining the snare on the foreign body, loosen the pin vise handle of the basket and advance it to the hub of the catheter.
6. Tighten the pin vise to secure the foreign body.
7. Withdraw the catheter and basket to a peripheral location, and retrieve the foreign body.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2011
Bloomington, IN 47404 U.S.A.