COOK

MEDICAL

CE
0088

EN

Spinal Needle

2

Instructions For Use

DA

Spinalkanyle

3

Brugsanvisning

DE

Spinalkanüle

4

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EL

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T _ D S N _ R E V O

SPINAL NEEDLE

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Disposable Spinal Needle (prefix DSN-) is a two-part needle consisting of a lancet beveled needle cannula and beveled needle stylet.

INTENDED USE

The Disposable Spinal Needle is intended for diagnostic sampling of cerebrospinal fluid and for the introduction of contrast medium for myelography.

CONTRAINDICATIONS

None known

WARNINGS

None known

PRECAUTIONS

Use of imaging technology (fluoroscopy, CT, ultrasound) is recommended when using this product.

INSTRUCTIONS FOR USE

1. Use imaging technology or anatomic landmarks to determine the external percutaneous entry site.
2. Use standard facility protocol to prepare the targeted entry site.
3. Remove the Disposable Spinal Needle assembly from the sterile packaging.
4. If desired, use a scalpel to make a skin nick at the external percutaneous entry site.
5. Insert the assembled Disposable Spinal Needle through the skin.
6. Using imaging guidance, advance the assembled needle to the desired depth or anatomic location.
7. Remove the beveled stylet.
8. You may inject contrast media or aspirate fluid to further confirm desired tip location in target anatomy.
9. Once desired location is confirmed, proceed with aspiration and/or injection of contrast media as desired.
10. Once the procedure is complete, remove the cannula and apply a sterile dressing.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2013
Bloomington, IN 47404 U.S.A.