COOK

MEDICAL

- EN**
3 Disposable Reuter Tip Deflecting Wire Guides
Instructions for Use
- DA**
4 Reuter flekterbare engangskateterledere
Brugsanvisning
- DE**
5 Reuter Einweg-Führungsdrähte mit steuerbarer Spitze
Gebrauchsanweisung
- EL**
6 Αναλώσιμοι συρμάτινοι οδηγοί κάμψης άκρου Reuter
Οδηγίες χρήσης
- ES**
7 Guías con punta deflectora Reuter desechables
Instrucciones de uso
- FR**
9 Guides à extrémité orientable Reuter jetables
Mode d'emploi
- IT**
10 Guide monouso con punta deflettibile Reuter
Istruzioni per l'uso
- NL**
11 Disposable Reuter tipdeflexievoerdraden
Gebruiksaanwijzing
- PT**
12 Fios guia descartáveis com ponta defletora reuter
Instruções de utilização
- SV**
14 Reuter engångsledare med styrbar spets
Bruksanvisning

T _ D T D W _ R E V 6

DISPOSABLE REUTER TIP DEFLECTING WIRE GUIDES

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

DEVICE DESCRIPTION

The Disposable Reuter Tip Deflecting Wire Guide is a handle/wire guide assembly that allows for controlled deflection of the wire guide tip within a catheter lumen. The handle is constructed of plastic, with a plastic side-arm adapter. The wire guides are constructed of stainless steel, and have a tip-deflected curve radius of 5 or 10 mm.

INTENDED USE

Disposable Reuter Tip Deflecting Wire Guides are intended for use in curving or deflecting catheter tips for pulmonary angiography, selective angiography, translumbar aortography, bronchography, repositioning of central venous catheter tips and other vascular and non-vascular applications.

CONTRAINDICATIONS

None known

WARNINGS

- Exercise caution when using a Disposable Reuter Tip Deflecting Wire Guide. This wire guide does not have a soft tip; therefore, to avoid damaging the vessel, **do not advance the distal tip of the wire guide beyond the distal tip of the catheter.**
- If the side-arm fitting is not used, exercise extra caution to ensure that wire guide tip does not extend beyond catheter tip.
- Do not use excessive force when deflecting a Disposable Reuter Tip Deflecting Wire Guide. Keep fingers and thumb in rings while deflecting wire guide. Throw length is designed to help limit excessive use of force in order to prevent damage to handle and wire guide.
- Do not attempt deflection of wire guide until it is inserted into a catheter. Such deflection may result in damage to tip of wire.
- When inserted in catheter, wire guide tip should be near the endhole, but should not be lodged in the catheter tip orifice. If the wire guide tip is lodged in the tip orifice of the catheter, flushing of the catheter tip is ineffective.
- Table-top test the handle, wire guide and catheter assembly prior to use. **Confirm that the wire guide does not protrude through or obstruct the catheter endhole. Ensure that the wire guide, when deflected, cannot exit through a catheter sideport.**
- Always monitor manipulation using fluoroscopic guidance.

PRECAUTIONS

- The product is intended for use by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for placement of vascular access sheaths, angiographic catheters and wire guides should be employed.
- Always use a Disposable Reuter Tip Deflecting Wire Guide with a cross-sectional diameter at least one size larger than the endhole diameter of the catheter with which it is paired (e.g., use a .038 inch wire guide with a .035 inch endhole catheter). This will help prevent wire guide tip from inadvertently exiting catheter tip during manipulation.
- A catheter with sideports may be used with the Disposable Reuter Tip Deflecting Wire Guide. However, care should be taken to ensure that the wire guide cannot accidentally exit through a catheter sideport.
- The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

INSTRUCTIONS FOR USE

NOTE: Throughout duration of procedure, periodically use side-arm fitting to flush system with heparinized saline solution to prevent clot formation on catheter or deflecting wire guide surface.

1. After catheter is properly positioned in vessel, perform a wire guide exchange for the Disposable Reuter Tip Deflecting Wire Guide/Handle assembly. Attach the male Luer lock on the side-arm fitting to the catheter hub.
2. To deflect the wire guide and catheter, insert thumb and fingers in proper rings and slowly but firmly pull fingers back to deflect tip.

WARNING: Exercise caution when using a Reuter Tip Deflecting Wire Guide. The distal tip of the wire guide should not be advanced beyond the distal tip of the catheter. Damage to the vessel wall can occur.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use

the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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Bloomington, IN 47404 U.S.A.