COOK
MEDICAL

CE
0123

EN
3

**Disposable Two-Part Trocar
Needle**

Instructions for Use

CS
4

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T _ D T P T N _ R E V O

DISPOSABLE TWO-PART TROCAR NEEDLE

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner). Please read all instructions before using this device.

DEVICE DESCRIPTION

The Disposable Two-Part Trocar Needle is a rigid needle used for percutaneous nephrostomy access. Gage and length will vary by configuration.

EchoTip® needle design provides enhanced visualization of needle tip when used with ultrasonic imaging equipment.

INTENDED USE

Used under fluoroscopy or ultrasound for initial location and positioning within the renal collecting system during a percutaneous nephrostomy.

CONTRAINDICATIONS

This device is contraindicated in the presence of conditions which create unacceptable risk during percutaneous nephrostomy.

WARNINGS

Do not force the device during the procedure; proceed with caution if excessive resistance is encountered.

PRECAUTIONS

This device is intended for use by physicians trained and experienced in the use of urologic and percutaneous access procedures. Standard techniques should be employed.

POTENTIAL ADVERSE EVENTS

Complications associated with percutaneous nephrostomy access procedures may include, but are not limited to:

- Bleeding
- Renal perforation
- Edema
- Hemorrhage
- Tissue injury

INSTRUCTIONS FOR USE

1. Remove device from its package using sterile technique.
2. Perform fluoroscopy or ultrasound to provide visualization of the target.
3. After prepping the skin in the usual manner, using fluoroscopy and/or ultrasound, make a needle puncture to provide direct access to the intended area of interest.
4. Remove the inner stylet to allow aspiration of urine and ensure correct placement of the needle.
5. Place an appropriately sized wire guide through the cannula of the needle to achieve access into the collecting system and allow placement of subsequent devices.
6. Withdraw the needle by grasping the needle at the skin so the depth of the collecting system can be gauged.

Upon completion of the procedure, dispose of the device per institutional guidelines for biohazardous medical waste.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened and undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

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Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd