EN
3

Dewan Suprapubic Urodynamic Catheter Set

Instructions for Use

CS
4

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T _ D U C S _ R E V 1

DEWAN SUPRAPUBIC URODYNAMIC CATHETER SET

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

DEVICE DESCRIPTION

Set Contains:

- 7.5 French dual lumen polyurethane urodynamic catheter, 25 cm long
- 18 gage stainless steel needle, 7 cm long
- 8.0 French TFE Peel-Away® sheath, 15 cm long
- 8.0 French radiopaque polyethylene dilator, 20 cm long
- .038 inch (0.97 mm) diameter stainless steel wire guide, 60 cm long

INTENDED USE

Used for suprapubic monitoring of bladder pressures. Suprapubic urodynamic studies avoid the effects of urethral catheterization and eliminate the need for a double bladder puncture.

CONTRAINDICATIONS

None known

WARNINGS

None known

PRECAUTIONS

The potential effects of phthalates on pregnant/nursing women or children have not been fully characterized and there may be concern for reproductive and developmental effects.

INSTRUCTIONS FOR USE

1. The bladder is filled transurethraly to a palpable state.
2. Unless scarring exists from previous surgery, the 18 gage needle is passed into the bladder about 2 finger breadths above the symphysis pubis and at a right angle to the abdomen. If scarring is present, then a 30 degree to 45 degree approach cephalad is suggested. **NOTE:** An initial scalpel nick in the skin may be helpful.
3. The floppy end of the wire guide is fed into the needle and advanced into the bladder about half way. The needle is removed.
4. The dilator/Peel-Away® sheath combination is advanced over the wire guide into the bladder. An oscillating motion may be helpful.
5. The dilator is removed, being careful to leave the wire guide in place.
6. The open-end tip of the dual lumen catheter is advanced over the wire guide and the catheter is positioned into the bladder through the Peel-Away® sheath.
7. The sheath is then peeled away and discarded.
8. The integral butterfly fitting is then sutured to the skin to stabilize the catheter externally and the wire guide is removed.
9. A bolus injection of sterile saline is then injected, via the "FILL" lumen, to fill the bladder.
10. Connect the urodynamic measurement equipment to the small lumen extension tube. Proceed with the urodynamic study.
11. To remove, cut stay sutures. Withdraw and dispose of the catheter.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

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© COOK 2019
Bloomington, IN 47404 U.S.A.
Cook Medical Europe Ltd