EN 16	Zenith® Low Profile AAA/Zenith® Alpha® Abdominal Ancillary Components Instructions for Use
CS 22	Přídavné komponenty Zenith® Low Profile AAA/Zenith® Alpha® Abdominal Návod k použití
DA 28	Zenith® AAA reduceret diameter/Zenith® Alpha® Abdominal hjælpekomponenter Brugsanvisning
DE 34	Zenith® Low Profile AAA/Zenith® Alpha® Abdominal Hilfskomponenten Gebrauchsanweisung
EL 41	Βοηθητικά εξαρτήματα χαμηλού προφίλ AAA Zenith®/Zenith® Alpha® Abdominal Οδηγίες χρήσης
ES 48	Componentes auxiliares de la endoprótesis vascular para AAA de perfil bajo Zenith®/endoprótesis vascular Zenith® Alpha® Abdominal Instrucciones de uso
FR 55	Composants auxiliaires Zenith® AAA bas profil/Zenith® Alpha® Abdominal Mode d'emploi
HU 61	Zenith® keskeny AAA/Zenith® Alpha® Abdominal kiegészítő komponensek Használati utasítás
IT 67	Componenti ausiliari per endoprotesi addominale Zenith® Low Profile AAA/Zenith® Alpha® Abdominal Istruzioni per l'uso
NL 74	Hulpcomponenten voor de Zenith® AAA met laag profiel/Zenith® Alpha® Abdominal Gebruiksaanwijzing
NO 81	Hjelpekomponenter for Zenith® lavprofil AAA / Zenith® Alpha® Abdominal Bruksanvisning
PL 87	Elementy pomocnicze do niskoprofilowego stent-graftu wewnątrznaczyniowego Zenith® AAA/stent-graftu wewnątrznaczyniowego Zenith® Alpha® Abdominal Instrukcja użycia
PT 94	Componentes auxiliares AAA Zenith® de baixo perfil e Zenith® Alpha® Abdominal Instruções de utilização
SV 100	Kompletterande komponenter till Zenith® AAA med låg profil/Zenith® Alpha® Abdominal Bruksanvisning

TABLE OF CONTENTS

Illustrations	10-15
1 DEVICE DESCRIPTION	16
1.1 Main Body Extensions	16
1.2 Converters	16
2 INDICATIONS FOR USE	16
3 CONTRAINDICATIONS	16
4 WARNINGS AND PRECAUTIONS	16
4.1 General	16
4.2 Patient Selection, Treatment and Follow-Up	16
4.3 Pre-Procedure Measurement Techniques and Imaging	17
4.4 Device Selection	17
4.5 Implant Procedure	17
4.6 Molding Balloon Use	17
4.7 MRI Information	17
5 POTENTIAL ADVERSE EVENTS	17
6 PATIENT SELECTION AND TREATMENT	18
6.1 Individualization of Treatment	18
7 PATIENT COUNSELING INFORMATION	18
8 HOW SUPPLIED	18
9 CLINICAL USE INFORMATION	18
9.1 Physician Training	18
9.2 Inspection Prior to Use	18
9.3 Materials Required	19
9.4 Materials Recommended	19
9.5 Device Diameter Sizing Guidelines	19
Table 9.5.1 Main Body Extension, ZLBE (Zenith Low Profile/Zenith Alpha Abdominal Body Extension), Graft Diameter Sizing Guide	19
Table 9.5.2 Converter, ZLC (Zenith Low Profile/Zenith Alpha Abdominal Converter), Graft Diameter Sizing Guide	19
10 DIRECTIONS FOR USE	19
Anatomical Requirements	19
General Use Information	19
Pre-Implant Determinants	19
Patient Preparation	19
10.1 Zenith Low Profile AAA/Zenith Alpha Abdominal Ancillary Components	19
General Use Information	19
10.1.1 Converter	19
10.1.2 Main Body Extensions	20
11 IMAGING GUIDELINES AND POSTOPERATIVE FOLLOW-UP	20
11.1 General	20
Table 11.1.1 Imaging Schedule	21
11.2 Contrast and Non-Contrast CT Recommendations	21
Table 11.2.1 Acceptable Imaging Protocols	21
11.3 Abdominal Radiographs	21
11.4 MRI Information	21
11.5 Additional Surveillance and Treatment	21
12 REFERENCES	21

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


ZENITH® LOW PROFILE AAA/ZENITH® ALPHA® ABDOMINAL ANCILLARY COMPONENTS

Read all instructions carefully. Failure to properly follow the instructions, warnings and precautions may lead to serious surgical consequences or injury to the patient.

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).

CAUTION: All contents of the outer pouch (including the introduction system and endovascular grafts) are supplied sterile, for single use only.

For the Zenith product line, there are several applicable Suggested Instructions for Use (IFU). This IFU describes the Suggested Instructions for Use for the Zenith Low Profile AAA/Zenith Alpha Abdominal Ancillary Components. For information on other Zenith components, please refer to the following Suggested Instructions for Use:

- Zenith Alpha® Spiral-Z Endovascular Leg;
- Zenith® Low Profile AAA Endovascular Graft;
- Zenith Alpha® Abdominal Endovascular Graft; and
- Coda® Balloon Catheter.

1 DEVICE DESCRIPTION

Ancillary endovascular components (main body extensions and converters) are available. (Fig. 1) Ancillary components are constructed from woven polyester fabric, self-expanding nitinol stents, and braided polyester and polypropylene suture, providing a conduit that is intended to exclude the aneurysm from blood flow.

The aortic main body extensions can be used to provide additional length to the proximal portion of the endovascular graft. The converters can be used to convert a bifurcated graft into an aortouni-iliac graft, if necessary (e.g., cases of Type III endoleak, limb occlusion or unattainable contralateral limb cannulation).

1.1 Main Body Extensions

Main body extensions utilize 16 French (6.0 mm OD) or 17 French (6.5 mm OD) H&L-B One-Shot Introduction Systems. (Fig. 14) The main body extension introduction system contains a single trigger-wire release mechanism. Deployment of the main body extension is achieved by sheath retraction and removing the distal trigger-wire.

To facilitate fluoroscopic visualization of the stent grafts, the main body extensions have 4 gold radiopaque markers positioned in a circumferential orientation within 2 mm of the most superior aspect of the graft material.

For added hemostasis, the Captor Hemostatic Valve can be loosened or tightened for the introduction and/or removal of ancillary devices into and out of the sheath. The main body extension delivery systems feature a Flexor introducer sheath which resists kinking and is hydrophilically coated. Both features are intended to enhance trackability in the iliac arteries and abdominal aorta.

1.2 Converters

Converters utilize 16 French (6.0 mm OD) or 17 French (6.5 mm OD) H&L-B One-Shot Introduction Systems. (Fig. 6) The converter introduction system contains a single trigger-wire release mechanism. Deployment of the converter is achieved by sheath retraction and removing the distal trigger-wire.

The converter has a single gold radiopaque marker positioned within 2 mm of the most superior aspect of the graft material. For added hemostasis, the Captor Hemostatic Valve can be loosened or tightened for the introduction and/or removal of ancillary devices into and out of the sheath. The converter delivery system features a Flexor introducer sheath, which resists kinking and is hydrophilically coated. Both features are intended to enhance trackability in the iliac arteries and abdominal aorta.

2 INDICATIONS FOR USE

The Zenith Low Profile AAA/Zenith Alpha Abdominal Ancillary Components are indicated for use with the Zenith Low Profile AAA Graft/Zenith Alpha Abdominal Endovascular Graft during either a primary or a secondary procedure in patients who have adequate iliac/femoral access compatible with the required introduction systems.

3 CONTRAINDICATIONS

The Zenith Low Profile AAA/Zenith Alpha Abdominal Ancillary Components are contraindicated in:

- Patients with known sensitivities or allergies to stainless steel, nitinol, polyester, solder (tin, silver), polypropylene, urethane, PTFE, nylon or gold.
- Patients with a systemic infection who may be at increased risk of endovascular graft infection.

4 WARNINGS AND PRECAUTIONS

4.1 General

- Read all instructions carefully. Failure to properly follow the instructions, warnings and precautions may lead to serious consequences or injury to the patient.
- Always have a qualified surgery team available during implantation or reintervention procedures in the event that conversion to open surgical repair is necessary.
- The Zenith Low Profile AAA/Zenith Alpha Abdominal Ancillary Components with the H&L-B One-Shot Introduction System should only be used by physicians and teams trained in vascular interventional techniques (catheter-based and surgical) and in the use of this device. Specific training expectations are described in **Section 9.1, Physician Training**.

- Additional endovascular interventions or conversion to standard open surgical repair following initial endovascular repair should be considered for patients experiencing enlarging aneurysms, unacceptable decrease in fixation length (vessel and component overlap) and/or endoleak. An increase in aneurysm size and/or persistent endoleak or migration may lead to aneurysm rupture.
- Patients experiencing reduced blood flow through the graft limb and/or leaks may be required to undergo secondary interventions or surgical procedures.
- The Zenith Low Profile AAA/Zenith Alpha Abdominal Ancillary Components have not been evaluated for use with grafts other than Zenith Low Profile AAA/Zenith Alpha Abdominal Stent Grafts.

4.2 Patient Selection, Treatment and Follow-Up

- The Zenith Low Profile AAA and Zenith Alpha Abdominal Endovascular main body extensions are designed to treat aortic neck diameters no smaller than 18 mm and no larger than 32 mm. The Zenith Low Profile AAA and Zenith Alpha Abdominal Endovascular main body extensions are designed to treat proximal aortic necks (distal to the lowest renal artery) of at least 15 mm in length.
- Key anatomical elements that may affect successful exclusion of the aneurysm include severe proximal neck angulation (>60 degrees for infrarenal neck to axis of AAA or >45 degrees for suprarenal neck relative to the immediate infrarenal neck); short proximal aortic neck (<15 mm); an inverted funnel shape (greater than 10% increase in diameter over 15 mm of proximal aortic neck length); and circumferential thrombus and/or calcification at the arterial implantation sites, specifically the proximal aortic neck and distal iliac artery interface. In the presence of anatomical limitations, a longer neck may be required to obtain adequate sealing and fixation. Irregular calcification and/or plaque may compromise the attachment and sealing at the fixation sites. Necks exhibiting these key anatomical elements may be more conducive to graft migration or endoleak.
- Adequate iliac or femoral access is required to introduce the device into the vasculature. Access vessel diameter (measured inner wall to inner wall) and morphology (minimal tortuosity, occlusive disease and/or calcification) should be compatible with vascular access techniques and delivery systems of the profile of a 16 French (6.0 mm OD) or 17 French (6.5 mm OD) vascular introducer sheath. Vessels that are significantly calcified, occlusive, tortuous or thrombus-lined may preclude placement of the endovascular graft and/or may increase the risk of embolization or thromboses. A vascular conduit technique may be necessary to achieve success in some patients.
- Pre-existing regions of stenosis/narrowing (less than approximately 20 mm ID in the aorta or 7 to 8 mm ID in the iliacs) have been shown to increase the risk of a thromboembolic event (e.g., graft limb occlusion). The potential for this increased risk in these patients may preclude placement of an endovascular graft. Dilatation of these regions with a noncompliant balloon and/or stent placement may be necessary to help assure maintained graft patency and to reduce the risk for a thromboembolic event. Additionally, the completion angiogram (with stiff wire guides removed) should be reviewed carefully to determine if further treatment in these regions is necessary (e.g., adjunctive ballooning or stenting). Failure to remove the stiff wire guide prior to the angiogram could mask any limb kinking or narrowing that might occur when the wire guide is removed.
- Follow-up imaging should be carefully reviewed for narrowing within the graft leg. Patients with a graft leg lumen less than approximately 5 mm ID may be at increased risk of a thromboembolic event (e.g., graft limb occlusion). Reintervention (e.g., noncompliant ballooning or stenting in these regions) should be considered to help assure maintained graft patency and to reduce the risk of a thromboembolic event.
- Patients with poor outflow or a hypercoagulable state (e.g., cancer) may be at an increased risk of a thromboembolic event.
- The Zenith Low Profile AAA/Zenith Alpha Abdominal Ancillary Components are not recommended in patients who cannot tolerate contrast agents necessary for intraoperative and postoperative follow-up imaging. All patients should be monitored closely and checked periodically for a change in the condition of their disease and the integrity of the endoprosthesis.
- The Zenith Low Profile AAA/Zenith Alpha Abdominal Ancillary Components are not recommended in patients exceeding weight and/or size limits which compromise or prevent the necessary imaging requirements.
- Inability to maintain patency of at least one internal iliac artery or occlusion of an indispensable inferior mesenteric artery may increase the risk of pelvic/bowel ischemia.
- Multiple large, patent lumbar arteries, mural thrombus and a patent inferior mesenteric artery may all predispose a patient to Type II endoleaks. Patients with uncorrectable coagulopathy may also have an increased risk of Type II endoleak or bleeding complications.
- The Zenith Low Profile AAA/Zenith Alpha Abdominal Ancillary Components have not been evaluated in the following patient populations:
 - traumatic aortic injury
 - leaking, pending rupture or ruptured aneurysms
 - mycotic aneurysms
 - pseudoaneurysms resulting from previous graft placement
 - revision of previously placed endovascular grafts
 - uncorrectable coagulopathy
 - indispensable mesenteric artery
 - genetic connective tissue disease (e.g., Marfans or Ehlers-Danlos Syndromes)
 - concomitant thoracic aortic or thoracoabdominal aneurysms
 - patients with active systemic infections
 - pregnant or nursing females
 - less than 18 years of age

- Successful patient selection requires specific imaging and accurate measurements; please see **Section 4.3 Pre-Procedure Measurement Techniques and Imaging**.
- All lengths and diameters of the devices necessary to complete the procedure should be available to the physician, especially when preoperative case planning measurements (treatment diameters/lengths) are not certain. This approach allows for greater intraoperative flexibility to achieve optimal procedural outcomes.

4.3 Pre-Procedure Measurement Techniques and Imaging

- Lack of non-contrast CT imaging may result in failure to appreciate iliac or aortic calcification, which may preclude access or reliable device fixation and seal.
- Pre-procedure imaging reconstruction thicknesses >3 mm may result in sub-optimal device sizing, or in failure to appreciate focal stenoses from CT.
- Clinical experience indicates that contrast-enhanced spiral computed tomographic angiography (CTA) with 3-D reconstruction is the strongly recommended imaging modality to accurately assess patient anatomy prior to treatment with the Zenith Low Profile AAA/Zenith Alpha Abdominal Ancillary Components. If contrast-enhanced spiral CTA with 3-D reconstruction is not available, the patient should be referred to a facility with these capabilities.
- Clinicians recommend positioning the x-ray C-arm during procedural angiography such that the origins of the renal arteries, and particularly the lowest patent renal artery, are well demonstrated prior to deployment of the proximal edge of the graft material (sealing stent) of the main body extension.

Diameters

Utilizing CT, diameter measurements should be determined from the outer wall to outer wall vessel diameter (not lumen measurement) to help with proper device sizing and device selection. The contrast-enhanced spiral CT scan must start 1 cm superior to the celiac axis and continue through the femoral heads at an axial thickness slice of 3 mm or less.

Lengths

Utilizing CT, length measurements should be determined to accurately assess infrarenal proximal neck length as well as planning main body sizes, iliac legs and ancillary components for the Zenith Low Profile AAA and Zenith Alpha Abdominal Endovascular Graft. These reconstructions should be performed in sagittal, coronal, and 3-D.

- **The long-term performance of endovascular grafts has not yet been established. All patients should be advised that endovascular treatment requires lifelong, regular follow-up to assess their health and the performance of their endovascular graft.** Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysm or changes in the structure or position of the endovascular graft) should receive enhanced follow-up. Specific follow-up guidelines are described in **Section 11, IMAGING GUIDELINES AND POSTOPERATIVE FOLLOW-UP**.
- The Zenith Low Profile AAA/Zenith Alpha Abdominal Ancillary Components are not recommended in patients unable to undergo, or who will not be compliant with, the necessary preoperative and postoperative imaging and implantation studies as described in **Section 11, IMAGING GUIDELINES AND POSTOPERATIVE FOLLOW-UP**.
- After endovascular graft placement, patients should be regularly monitored for perigraft flow, aneurysm growth or changes in the structure or position of the endovascular graft. At a minimum, annual imaging is required, including: 1) abdominal radiographs to examine device integrity (separation between components, stent fracture or barb separation) and 2) contrast and non-contrast CT to examine aneurysm changes, perigraft flow, patency, tortuosity and progressive disease. If renal complications or other factors preclude the use of image contrast media, abdominal radiographs and duplex ultrasound may provide similar information.

4.4 Device Selection

Strict adherence to the Zenith Low Profile AAA/Zenith Alpha Abdominal Ancillary Components IFU sizing guide is strongly recommended when selecting the appropriate device size (**Tables 9.5.1** thru **9.5.2**). Appropriate device oversizing has been incorporated into the IFU sizing guide. Sizing outside of this range can result in endoleak, fracture, migration, device infolding or compression.

4.5 Implant Procedure

(Refer to **Section 10, DIRECTIONS FOR USE**)

- Appropriate procedural imaging is required to successfully position the Zenith Low Profile AAA/Zenith Alpha Abdominal Ancillary Components and assure accurate apposition to the aortic wall.
- Do not bend or kink the delivery system. Doing so may cause damage to the delivery system and the Zenith Low Profile AAA/Zenith Alpha Abdominal Ancillary Components.
- To avoid any twist in the endovascular graft, during any rotation of the delivery system, be careful to rotate all of the components of the system together (from outer sheath to inner cannula).
- To avoid damage to the sheath, be careful to advance all components of the system together (from outer sheath to inner cannula).
- Do not continue advancing any portion of the delivery system if resistance is felt during advancement of the wire guide or delivery system. Stop and assess the cause of resistance; vessel, catheter or graft damage may occur. Exercise particular care in areas of stenosis, intravascular thrombosis or in calcified or tortuous vessels.
- Inadvertent partial deployment or migration of the endoprosthesis may require surgical removal.
- Unless medically indicated, do not deploy the Zenith Low Profile AAA/Zenith Alpha Abdominal Ancillary Components in a location that will occlude arteries necessary to supply blood flow to organs or extremities. Do not cover significant renal or mesenteric arteries (exception is the inferior mesenteric artery) with the endoprosthesis. Vessel occlusion may occur. During the clinical study, this device was not studied in patients with two occluded internal iliac arteries.
- Do not attempt to re-sheath the graft after partial or complete deployment.
- Repositioning the stent graft distally after partial deployment of the covered proximal stent may result in damage to the stent graft and/or vessel injury.

- Inaccurate placement and/or incomplete sealing of the Zenith Low Profile AAA/Zenith Alpha Abdominal Ancillary Components within the vessel may result in increased risk of endoleak, migration or inadvertent occlusion of the renal or internal iliac arteries. Renal artery patency must be maintained to prevent/reduce the risk of renal failure and subsequent complications.
- Inadequate fixation of the Zenith Low Profile AAA/Zenith Alpha Abdominal Ancillary Components may result in increased risk of migration of the stent graft. Incorrect deployment or migration of the endoprosthesis may require surgical intervention.
- Systemic anticoagulation should be used during the implantation procedure based on hospital and physician preferred protocol. If heparin is contraindicated, an alternative anticoagulant should be considered.
- To activate the hydrophilic coating on the outside of the Flexor introducer sheath, the surface must be wiped with sterile gauze pads soaked in saline solution. Always keep the sheath hydrated for optimal performance.
- Minimize handling of the constrained endoprosthesis during preparation and insertion to decrease the risk of endoprosthesis contamination and infection.
- Maintain wire guide position during delivery system insertion.
- Fluoroscopy should be used during introduction and deployment to confirm proper operation of the delivery system components, proper placement of the graft, and desired procedural outcome.
- The use of the Zenith Low Profile AAA/Zenith Alpha Abdominal Ancillary Components require administration of intravascular contrast. Patients with pre-existing renal insufficiency may have an increased risk of renal failure postoperatively. Care should be taken to limit the amount of contrast media used during the procedure and to observe preventative methods of treatment to decrease renal compromise.
- As the sheath and/or wire guide is withdrawn, anatomy and graft position may change. Constantly monitor graft position and perform angiography to check position as necessary.
- The converter and main body extension delivery systems cannot be introduced through a main body or iliac leg introducer sheath.
- Use caution during manipulation of catheters, wire and sheaths within an aneurysm. Significant disturbances may dislodge fragments of thrombus, which can cause distal embolization or rupture of the aneurysm.
- Take care not to displace or damage the main body graft during placement or deployment.

4.6 Molding Balloon Use

- Do not inflate the balloon in vessel outside of the graft, as doing so may cause damage to the vessel. Use the balloon in accordance with its labeling.
- Use care when inflating the balloon within the graft in the presence of calcification, as excessive inflation may cause damage to the vessel.
- Confirm complete deflation of the balloon prior to repositioning.
- For added hemostasis, the Captor Hemostatic Valve can be loosened or tightened to accommodate the insertion and subsequent withdrawal of a molding balloon.

Main Body Extension and Converter

- Care should be taken not to displace the main body graft during the placement and deployment of the main body extension or converter.
- To activate the hydrophilic coating on the outside of the Flexor introducer sheath, the surface must be wiped with sterile gauze pads soaked in saline solution. Always keep the sheath hydrated for optimal performance.

4.7 MRI Information

NOTE: If using this device in conjunction with another endovascular graft from the Zenith family, refer to the appropriate device's IFU for additional MRI information.

Nonclinical testing has demonstrated that the Zenith Low Profile AAA/Zenith Alpha Abdominal Ancillary Components are **MR Conditional** according to ASTM F2503. A patient with this device may be safely scanned under the following conditions:

Static Magnetic Field

- Static magnetic field of 3.0 tesla or 1.5 tesla only
- Maximum spatial magnetic gradient of 1600 gauss/cm (16.0 T/m) or less
- Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of ≤ 2.0 W/kg (Normal Operating Mode) for 15 minutes of continuous scanning.

Under the scan conditions provided above, the Zenith Low Profile AAA/Zenith Alpha Abdominal Ancillary Components are not expected to result in a temperature rise of more than 1.3°C after 15 minutes of continuous scanning.

Image Artifact

The image artifact extends approximately 9.5 mm from Zenith Low Profile AAA/Zenith Alpha Abdominal Ancillary Components as found during nonclinical testing when imaged with a gradient echo pulse sequence and 3.0 Tesla MRI system. The image artifact obscures the device lumen.

For US Patients Only

Cook recommends that the patient register the MR conditions disclosed in this IFU with the MedAlert Foundation. The MedAlert Foundation can be contacted in the following manners:

Mail: MedAlert Foundation International
2323 Colorado Avenue
Turlock, CA 95382

Phone: 888-633-4298 (toll free)
209-668-3333 from outside the US
Fax: 209-669-2450

Web: www.medicalert.org

5 POTENTIAL ADVERSE EVENTS

Adverse events that may occur and/or require intervention include:

- Amputation
- Anesthetic complications and subsequent attendant problems (e.g., aspiration)

- Aneurysm enlargement
- Aneurysm rupture and death
- Aortic damage, including perforation, dissection, bleeding, rupture and death
- Aortic fistulation to adjacent organs/anatomic structures
- Arterial or venous thrombosis and/or pseudoaneurysm
- Arteriovenous fistula
- Bleeding, hematoma or coagulopathy
- Bowel complications (e.g., ileus, transient ischemia, infarction, necrosis)
- Cardiac complications and subsequent attendant problems (e.g., arrhythmia, myocardial infarction, congestive heart failure, hypotension, hypertension)
- Claudication (e.g., buttock, lower limb)
- Death
- Edema
- Embolization (micro and macro) with transient or permanent ischemia or infarction
- Endoleak
- Endoprosthesis: improper component placement; incomplete component deployment; component migration; suture break; occlusion; infection; stent fracture; graft material wear; dilatation; erosion; puncture; perigraft flow; barb separation and corrosion
- Fever and localized inflammation
- Genitourinary complications and subsequent attendant problems (e.g., ischemia, erosion, fistula, incontinence, hematuria, infection)
- Graft or native vessel occlusion
- Hepatic failure
- Impotence
- Infection of the aneurysm, device or access site, including abscess formation, transient fever and pain
- Limb stenosis
- Lymphatic complications and subsequent attendant problems (e.g., lymph fistula)
- Neurologic local or systemic complications and subsequent attendant problems (e.g., stroke, transient ischemic attack, paraplegia, paraparesis, paralysis)
- Pulmonary/respiratory complications and subsequent attendant problems (e.g., pneumonia, respiratory failure, prolonged intubation)
- Radiation injury and/or late malignancy
- Renal complications and subsequent attendant problems (e.g., artery occlusion, contrast toxicity, insufficiency, failure)
- Surgical conversion to open repair
- Vascular access site complications, including infection, pain, hematoma, pseudoaneurysm, arteriovenous fistula
- Vascular spasm or vascular trauma (e.g., iliofemoral vessel dissection, bleeding, rupture, death)
- Vessel damage
- Wound complications and subsequent attendant problems (e.g., dehiscence, infection)

6 PATIENT SELECTION AND TREATMENT

(Refer to **Section 4, WARNINGS AND PRECAUTIONS**)

6.1 Individualization of Treatment

Cook recommends that the Zenith Low Profile AAA/Zenith Alpha Abdominal Ancillary Component component diameters be selected as described in **Tables 9.5.1 and 9.5.2** within **Section 9, CLINICAL USE INFORMATION**. All lengths and diameters of the devices necessary to complete the procedure should be available to the physician, especially when preoperative case planning measurements (treatment diameters/lengths) are not certain. This approach allows for greater intraoperative flexibility to achieve optimal procedural outcomes. Additional considerations for patient selection include, but are not limited to:

- Patient's age and life expectancy
- Co-morbidities (e.g., cardiac, pulmonary or renal insufficiency prior to surgery, morbid obesity)
- Patient's suitability for open surgical repair
- Patient's anatomical suitability for endovascular repair
- The risk of aneurysm rupture compared to the risk of treatment with the Zenith Low Profile AAA /Zenith Alpha Abdominal Ancillary Component
- Ability to tolerate general, regional or local anesthesia
- Ilio-femoral access vessel size and morphology (minimal thrombus, calcium and/or tortuosity) should be compatible with vascular access techniques and accessories of the delivery profile of a 16 French (6.0 mm OD) or 17 French (6.5 mm OD) vascular introducer sheath.
- For main body extensions, nonaneurysmal infrarenal aortic segment (neck) proximal to the aneurysm:
 - with a diameter measured outer wall to outer wall of no greater than 32 mm and no less than 18 mm,
 - with an angle less than 60 degrees relative to the long axis of the aneurysm, and
 - with an angle less than 45 degrees relative to the axis of the suprarenal aorta.
- Freedom from significant femoral/iliac artery occlusive disease that would impede flow through the endovascular graft.

The final treatment decision is at the discretion of the physician and patient.

7 PATIENT COUNSELING INFORMATION

The physician and patient (and/or family members) should review the risks and benefits when discussing this endovascular device and procedure, including:

- Risks and differences between endovascular repair and surgical repair
- Potential advantages of traditional open surgical repair
- Potential advantages of endovascular repair
- The possibility that subsequent interventional or open surgical repair of the aneurysm may be required after initial endovascular repair

In addition to the risks and benefits of an endovascular repair, the physician should assess the patient's commitment and compliance to postoperative follow-up as necessary to ensure continuing safe and effective results. Listed below are additional topics to discuss with the patient as to expectations after an endovascular repair:

- **The long-term performance of endovascular grafts has not yet been established. All patients should be advised that endovascular treatment requires lifelong, regular follow-up to assess their health and the performance of their endovascular graft.** Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms or changes in the structure or position of the endovascular graft) should receive enhanced follow-up. Specific follow-up guidelines are described in **Section 11, IMAGING GUIDELINES AND POSTOPERATIVE FOLLOW-UP**.
- Patients should be counseled on the importance of adhering to the follow-up schedule, both during the first year and at yearly intervals thereafter. Patients should be told that regular and consistent follow-up is a critical part of ensuring the ongoing safety and effectiveness of endovascular treatment of AAAs. At a minimum, annual imaging and adherence to routine postoperative follow-up requirements is required and should be considered a lifelong commitment to the patient's health and well-being.
- The patient should be told that successful aneurysm repair does not arrest the disease process. It is still possible to have associated degeneration of vessels.
- Physicians must advise all patients that it is important to seek prompt medical attention if he/she experiences signs of limb occlusion, aneurysm enlargement or rupture. Signs of graft limb occlusion include pain in the hip(s) or leg(s) during walking or at rest or discoloration or coolness of the leg. Aneurysm rupture may be asymptomatic, but usually presents as: pain; numbness; weakness in the legs; any back, chest, abdominal or groin pain; dizziness; fainting; rapid heartbeat or sudden weakness.
- Due to the imaging required for successful placement and follow-up of endovascular devices, the risks of radiation exposure to developing tissue should be discussed with women who are or suspect they are pregnant. Men who undergo endovascular or open surgical repair may experience impotence.

Physicians should refer patients to the *Patient Guide* regarding risks occurring during or after implantation of the device. Procedure-related risks include occlusion, endoleak, aneurysm enlargement, fracture, potential for reintervention and open surgical conversion, rupture and death (See **Section 5, POTENTIAL ADVERSE EVENTS**). The physician should complete the *Patient ID Card* and give it to the patient so that he/she can carry it with him/her at all times. The patient should refer to the card anytime he/she visits additional health practitioners, particularly for any additional diagnostic procedures (e.g., MRI).

8 HOW SUPPLIED

- **STERILE - DO NOT RESTERILIZE - SINGLE USE ONLY.**
- The Zenith Low Profile AAA/Zenith Alpha Abdominal Ancillary Components are supplied sterile (100% Ethylene Oxide) and pre-loaded in peel-open packages.
- The device is intended for single use only. Do not re-sterilize the device.
- Inspect the device and packaging to verify that no damage has occurred as a result of shipping. Do not use this device if damage has occurred or if the sterilization barrier has been damaged or broken. If damage has occurred, do not use the product and return it to Cook.
- Prior to use, verify correct devices (quantity and size) have been supplied for the patient by matching the device to the order prescribed by the physician for that particular patient.
- The main body extensions, and converter devices are loaded into a 16 French (6.0 mm OD) or 17 French (6.5 mm OD) Flexor introducer sheath. The sheath's surface is treated with a hydrophilic coating that, when hydrated, enhances trackability. To activate the hydrophilic coating, the surface must be wiped with a sterile gauze pad soaked in saline solution.
- Do not use after the expiration date printed on the label.
- Store in a cool, dry place.

9 CLINICAL USE INFORMATION

9.1 Physician Training

CAUTION: Always have a qualified surgery team available during implantation or reintervention procedures in the event that conversion to open surgical repair is necessary.

CAUTION: The Zenith Low Profile AAA/Zenith Alpha Abdominal Ancillary Components should only be used by physicians and teams trained in vascular interventional techniques and in the use of this device. The recommended skill/knowledge requirements for physicians using the Zenith Low Profile AAA/Zenith Alpha Abdominal Ancillary Components are outlined below.

Patient selection:

- Knowledge of the natural history of abdominal aortic aneurysms (AAA) and comorbidities associated with AAA repair
- Knowledge of radiographic image interpretation, device selection, planning and sizing

A multi-disciplinary team that has combined procedural experience with:

- Femoral cutdown, arteriotomy and repair
- Percutaneous access and closure techniques
- Non-selective and selective wire guide and catheter techniques
- Fluoroscopic and angiographic image interpretation
- Embolization
- Angioplasty
- Endovascular stent placement
- Snare techniques
- Appropriate use of radiographic contrast material
- Techniques to minimize radiation exposure
- Expertise in necessary patient follow-up modalities

9.2 Inspection Prior to Use

Inspect the device and packaging to verify that no damage has occurred as a result of shipping. Do not use this device if damage has occurred or if the sterilization barrier has been damaged or broken. If damage has occurred, do not use the product and return it to Cook. Prior to use, verify correct devices (quantity and size) have been supplied for the patient by matching the device to the order prescribed by the physician for that particular patient.

9.3 Materials Required

- Fluoroscope with digital angiography capabilities (C-arm or fixed unit)
- Contrast media
- Syringe
- Heparinized saline solution
- Sterile gauze pads

9.4 Materials Recommended

- .035 inch (0.89 mm) extra stiff wire guide, 260 cm; for example:
 - Cook Lunderquist Extra Stiff Wire Guides (LES)
- .035 inch (0.89 mm) standard wire guide; for example:
 - Cook .035 inch wire guides
 - Cook Nimble™ Wire Guides
- Molding Balloons; for example:
 - Cook Coda® Balloon Catheter

- Introducer sets; for example:
 - Cook Check-Flo® Introducer Sets
 - Cook Extra Large Check-Flo® Introducer Sets
 - Cook Flexor® Balkin Up & Over® Contralateral Introducers
- Sizing catheter; for example:
 - Cook Auros® Centimeter Sizing Catheters
- Angiographic radiopaque tip catheters; for example:
 - Cook Beacon® Tip Angiographic Catheters
 - Cook Beacon® Tip Royal Flush Catheters
- Entry needles; for example:
 - Cook single wall entry needles

9.5 Device Diameter Sizing Guidelines

The choice of diameter should be determined from the outer wall to outer wall vessel diameter and not the lumen diameter. Undersizing or oversizing may result in incomplete sealing or compromised flow.

Table 9.5.1 Main Body Extension, ZLBE (Zenith Low Profile/Zenith Alpha Abdominal Body Extension), Graft Diameter Sizing Guide*

Intended Aortic Vessel Diameter ^{1,2} (mm)	Main Body Extension Diameter ³ (mm)	Main Body Extension Length (mm)	Introducer Sheath	
			French Size	ID, OD (mm)
18-19	22	45, 58	16	5.3, 6.0
20-21	24	45, 58	16	5.3, 6.0
22	26	45, 58	16	5.3, 6.0
23-24	28	45, 58	16	5.3, 6.0
25-26	30	45, 58	16	5.3, 6.0
27-28	32	45, 58	16	5.3, 6.0
29-32	36	45, 58	17	5.6, 6.5

¹ Maximum diameter along the proximal fixation site.

² Round measured aortic diameter to nearest mm.

³ Additional considerations may affect choice of diameter.

*All dimensions are nominal.

Table 9.5.2 Converter, ZLC (Zenith Low Profile/Zenith Alpha Abdominal Converter), Graft Diameter Sizing Guide*

Main Body Diameter (mm)	Converter Diameter ¹ (mm)	Converter Length (mm)	Introducer Sheath	
			French Size	ID, OD (mm)
22 & 24	24	66	16	5.3, 6.0
26 & 28	28	66	16	5.3, 6.0
30 & 32	32	66	16	5.3, 6.0
36	36	66	17	5.6, 6.5

¹ Additional considerations may affect choice of diameter.

*All dimensions are nominal.

10 DIRECTIONS FOR USE

Anatomical Requirements

- Iliofemoral access vessel size and morphology (minimal thrombus, calcium and/or tortuosity) should be compatible with vascular access techniques and accessories. Arterial conduit techniques may be required.
- For use of the main body extension, proximal aortic neck lengths should be a minimum of 15 mm with a diameter measured outer wall to outer wall of 18 to 32 mm.

Prior to use of the Zenith Low Profile AAA/Zenith Alpha Abdominal Ancillary Component, review this Suggested Instructions for Use booklet. The following instructions embody a basic guideline for device placement. Variations in the following procedures may be necessary. These instructions are intended to help guide the physician and do not take the place of physician judgment.

General Use Information

- Standard techniques for placement of arterial access sheaths, guiding catheters, angiographic catheters and wire guides should be employed during use of the Zenith Low Profile AAA/Zenith Alpha Abdominal Ancillary Components. The Zenith Low Profile AAA/Zenith Alpha Abdominal Ancillary Components are compatible with .035 inch diameter wire guides.
- Endovascular stent grafting is a surgical procedure, and blood loss from various causes may occur, infrequently requiring intervention (including transfusion) to prevent adverse outcomes. It is important to monitor blood loss from the hemostatic valve throughout the procedure, but is specifically relevant during and after manipulation of the gray positioner. After the gray positioner has been removed, if blood loss is excessive, consider placing an uninflated molding balloon or an introduction system dilator within the valve, restricting flow.

Pre-Implant Determinants

Verify from pre-implant planning that the correct device has been selected. Determinants include:

1. Femoral artery selection for introduction of the main body system (i.e., define respective contralateral and ipsilateral iliac arteries).
2. Angulation of aortic neck, aneurysm and iliac arteries.
3. Quality of the aortic neck.
4. Diameters of infrarenal aortic neck and distal iliac arteries.
5. Distance from renal arteries to the aortic bifurcation.
6. Length from the aortic bifurcation to the internal iliac arteries/attachment site(s).
7. Aneurysm(s) extending into the iliac arteries may require special consideration in selecting a suitable graft/artery interface site.
8. Consider the degree of vascular calcification.

Patient Preparation

1. Refer to institutional protocols relating to anesthesia, anticoagulation and monitoring of vital signs.
2. Position patient on imaging table allowing fluoroscopic visualization from the aortic arch to the femoral bifurcations.
3. Both common femoral arteries should be prepared using standard techniques for either surgical or percutaneous access.

10.1 Zenith Low Profile AAA/Zenith Alpha Abdominal Ancillary Components

General Use Information

Inaccuracies in device size selection or placement, changes or anomalies in patient anatomy or procedural complications can require placement of additional endovascular grafts, extensions and converters. It is vital to maintain wire guide access.

Standard techniques for placement of arterial access sheaths, guiding catheters, angiographic catheters and wire guides should be employed during use of the Zenith Low Profile AAA/Zenith Alpha Abdominal Ancillary Components.

The Zenith Low Profile AAA/Zenith Alpha Abdominal Ancillary Components are compatible with .035 inch wire guides.

10.1.1 Converter

Converters can be used to convert a bifurcated graft into an aortouni-iliac graft, if necessary (e.g., cases of Type III endoleak, limb occlusion or unattainable contralateral limb cannulation). (Fig. 6)

Converter Preparation/Flush

1. Remove inner stylet (from the inner cannula), cannula protector (from the inner cannula) and dilator tip protector (from the dilator tip). Remove Peel-Away sheath from back of the hemostatic valve. (Fig. 7) Elevate distal tip of the system and flush through the stopcock on the hemostatic valve until fluid emerges from the flushing groove near the tip of the introducer sheath. (Fig. 3) Continue to inject a full 20 cc of flushing solution through the device. Discontinue injection and close stopcock.

NOTE: Graft flushing solution of heparinized saline is often used.

2. Attach syringe with heparinized saline to the hub of the inner cannula. Flush until fluid exits the dilator tip. (Fig. 4)

NOTE: When flushing system, elevate distal end of the system to facilitate removal of air.

3. Soak sterile gauze pads in saline solution and use to wipe Flexor introducer sheath to activate the hydrophilic coating. Hydrate both sheath and dilator liberally.

Converter Placement and Deployment

1. Remove main body delivery sheath. Use the main body graft wire guide to introduce the converter into the main body.

NOTE: The converter delivery system cannot be introduced through the main body or iliac leg introducer sheath.

2. Advance slowly until the converter is at the site of the required intervention. (Fig. 8) Verify appropriate stent graft overlap to ensure proper sealing and resistance to migration. The proximal two stents should be positioned in the main body graft and the distal two stents should be positioned in the ipsilateral leg.
3. Ensure the Captor Hemostatic Valve is turned to the open position.
4. Deploy the device by withdrawing the sheath while stabilizing the gray positioner of the delivery system. (Figs. 9 and 10)
5. Remove the safety lock from the black trigger-wire release mechanism. Withdraw and remove the trigger-wire by sliding the black trigger-wire release mechanism off the handle, and then remove via the slot over the inner cannula. (Fig. 11)
6. Continue to deploy the device until the distal stent is uncovered.
7. Withdraw the tapered tip of the introducer back through the converter graft and delivery system while maintaining wire guide position. Ensure the endovascular graft is not displaced during the withdrawal of delivery system.
8. Close the Captor Hemostatic Valve by turning it in a clockwise direction until it stops. (Fig. 12)

Converter Molding Balloon Insertion

NOTE: For information on the use of recommended products, refer to the individual product's Instructions for Use.

1. Prepare molding balloon as follows:
 - Flush wire lumen with heparinized saline.
 - Remove all air from balloon.

CAUTION: Captor Hemostatic Valve must be open prior to repositioning of molding balloon.

2. Advance the molding balloon over the wire guide and through the hemostatic valve to the proximal segment of the converter.
3. Tighten the Captor Hemostatic Valve around the molding balloon with gentle pressure by turning it clockwise.
4. Expand the molding balloon within the proximal segment and the distal segment of the converter using diluted contrast media (as recommended by the manufacturer). (Fig. 13)

CAUTION: Confirm complete deflation of the balloon prior to repositioning.

5. Completely deflate and remove the molding balloon, replace it with an angiographic catheter and perform completion angiograms.
6. If no other endovascular maneuvers are necessary, remove any sheaths, wires and catheters. Repair vessels and close in standard surgical fashion.

10.1.2 Main Body Extensions

Main body extensions are used for extending the proximal body of an *in situ* endovascular graft. (Fig. 14)

Main Body Extension Preparation/Flush

1. Remove inner stylet (from the inner cannula), cannula protector (from the inner cannula) and dilator tip protector (from the dilator tip). Remove Peel-Away sheath from back of the hemostatic valve. (Fig. 7) Elevate distal tip of the system and flush through the stopcock on the hemostatic valve until fluid emerges from the flushing groove near the tip of the introducer sheath. (Fig. 3) Continue to inject a full 20 cc of flushing solution through the device. Discontinue injection and close stopcock.

NOTE: Graft flushing solution of heparinized saline is often used.

2. Attach syringe with heparinized saline to the hub of the inner cannula. Flush until fluid exits the dilator tip. (Fig. 4)

NOTE: When flushing system, elevate distal tip of the system to facilitate removal of air.

3. Soak sterile gauze pads in saline solution and use to wipe Flexor introducer sheath to activate the hydrophilic coating. Hydrate both sheath and dilator liberally.

Main Body Extension Placement and Deployment

1. Remove main body delivery sheath. Use the main body graft wire guide to introduce the main body extension into the main body.

NOTE: The main body extension delivery system cannot be introduced through the main body or iliac leg introducer sheath.

2. Advance slowly until the main body extension is at the site of the required intervention. (Fig. 15)
3. Verify the main body extension position to ensure proper sealing and resistance to migration.
4. Verify the placement with angiography to ensure that the renal arteries remain patent and that proper placement is achieved.

CAUTION: Care should be taken not to displace the main body graft during the placement and deployment of the main body extension.

5. Ensure the Captor Hemostatic Valve is turned to the open position.
6. Deploy the device by withdrawing the sheath while stabilizing the gray positioner of the delivery system. (Figs. 9 and 16) Continue to deploy the device until the most distal stent is uncovered. Stop withdrawing the sheath.
7. Remove the safety lock from the black trigger-wire release mechanism. Withdraw and remove the trigger-wire by sliding the black trigger-wire release mechanism off the handle, and then remove via the slot over the inner cannula. (Fig. 11)
8. Withdraw the tapered tip of the introducer back through the main body extension graft and delivery system while maintaining wire guide position. Ensure the main body extension and endovascular graft are not displaced during withdrawal of the delivery system.
9. Close the Captor Hemostatic Valve by turning it in a clockwise direction until it stops. (Fig. 12)

Main Body Extension Molding Balloon Insertion

NOTE: For information on the use of recommended products, refer to the individual product's Instructions for Use.

1. Prepare molding balloon as follows:
 - Flush wire lumen with heparinized saline.
 - Remove all air from balloon.

CAUTION: Captor Hemostatic Valve must be open prior to repositioning of molding balloon.

2. Advance the molding balloon over the wire guide and through the hemostatic valve of the main body introduction system to the level of the main body extension.
3. Tighten the Captor Hemostatic Valve around the molding balloon with gentle pressure by turning it clockwise.

CAUTION: Do not inflate the balloon in the vessel outside of the graft.

4. Expand the molding balloon within the proximal segment of the main body extension and then the most distal segment of the main body extension using diluted contrast media (as recommended by the manufacturer). (Fig. 17)

CAUTION: Confirm complete deflation of the balloon prior to repositioning.

5. Completely deflate and remove the molding balloon, replace it with an angiographic catheter and perform completion angiograms.
6. If no other endovascular maneuvers are necessary, remove any sheaths, wires and catheters. Repair vessels and close in standard surgical fashion.

11 IMAGING GUIDELINES AND POSTOPERATIVE FOLLOW-UP

11.1 General

- **The long-term performance of endovascular grafts has not yet been established. All patients should be advised that endovascular treatment requires lifelong, regular follow-up to assess their health and the performance of their endovascular graft.** Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms or changes in the structure or position of the endovascular graft) should receive additional follow-up.
- Patients should be counseled on the importance of adhering to the follow-up schedule, both during the first year and at yearly intervals thereafter. Patients should be told that regular and consistent follow-up is a critical part of ensuring the ongoing safety and effectiveness of endovascular treatment of AAAs.
- Physicians should evaluate patients on an individual basis and prescribe their follow-up relative to the needs and circumstances of each individual patient. The recommended imaging schedule is presented in **Table 11.1.1**. This schedule continues to be the minimum requirement for patient follow-up and should be maintained even in the absence of clinical symptoms (e.g., pain, numbness, weakness). Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms or changes in the structure or position of the stent graft) should receive follow-up at more frequent intervals.
- Annual imaging follow-up should include abdominal radiographs and both contrast and non-contrast CT examinations.
- The combination of contrast and non-contrast CT imaging provides information on aneurysm diameter change, endoleak, patency, tortuosity, progressive disease, fixation length and other morphological changes.
- The abdominal radiographs provide information on device integrity (separation between components, stent fracture and barb separation).

Table 11.1.1 lists the minimum imaging follow-up for patients with a Zenith Low Profile AAA Endovascular Graft or a Zenith Alpha Abdominal Endovascular Graft including Ancillary Components. Patients requiring enhanced follow-up should have interim evaluations.

Table 11.1.1 Imaging Schedule

	Pre-op	Intra-op	30-Day	6-Month	12-Month ⁴
CT Scan	X ¹		X ³	X ³	X ³
Abdominal Device X-ray			X	X	X
Angiography	X ²	X			

¹ Imaging should be performed within 6 months before the procedure.

² Required only to resolve any uncertainties in anatomical measurements necessary for graft sizing.

³ Duplex ultrasound may be used for those patients experiencing renal failure or who are otherwise unable to undergo contrast enhanced CT scan.

⁴ Yearly thereafter

11.2 Contrast and Non-Contrast CT Recommendations

- Film sets should include all sequential images at lowest possible slice thickness (≤3 mm). DO NOT perform large slice thickness (>3 mm) and/or omit consecutive CT images/film sets, as it prevents precise anatomical and device comparisons over time.
- All images should include a scale for each film/image. Images should be arranged no smaller than 20:1 images on 14 inch x 17 inch sheets if film is used. It is important to follow acceptable imaging protocols during the CT exam. **Table 11.2.1** lists examples of acceptable imaging protocols

Table 11.2.1 Acceptable Imaging Protocols

	Non-Contrast	Contrast
IV contrast	No	Yes
Acceptable machines	Spiral CT or high performance MDCT capable of >40 seconds	Spiral CT or high performance MDCT capable of >40 seconds
Injection volume	n/a	Per institutional protocol
Injection rate	n/a	>2.5 cc/sec
Injection mode	n/a	Power
Bolus timing	n/a	Test bolus: SmartPrep, C.A.R.E. or equivalent
Coverage – start	Diaphragm	1 cm superior to celiac axis
Coverage – finish	Proximal femur	Profunda femoris origin
Collimation	<3 mm	<3 mm
Reconstruction	2.5 mm throughout – soft algorithm	2.5 mm throughout – soft algorithm
Axial DFOV	32 cm	32 cm
Post-injection runs	None	None

11.3 Abdominal Radiographs

The following views are required:

- Four films: supine-frontal (AP), lateral, 30 degree LPO and 30 degree RPO views centered on umbilicus.
- Record the table-to-film distance and use the same distance at each subsequent examination.

Ensure entire device is captured on each single image format lengthwise.

If there is any concern about the device integrity (e.g., kinking, stent breaks, barb separation, relative component migration), it is recommended to use magnified views. The attending physician should evaluate films for device integrity (entire device length including components) using 2-4X magnification visual aid.

11.4 MRI Information

NOTE: If using this device in conjunction with another endovascular graft from the Zenith family, refer to the appropriate device's IFU for additional MRI information.

Nonclinical testing has demonstrated that the Zenith Low Profile AAA/ Zenith Alpha Abdominal Ancillary Components are **MR Conditional** according to ASTM F2503. A patient with this device may be safely scanned under the following conditions:

Static Magnetic Field

- Static magnetic field of 3.0 tesla or 1.5 tesla only
- Maximum spatial magnetic gradient of 1600 gauss/cm (16.0 T/m) or less
- Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of ≤ 2.0 W/kg (Normal Operating Mode) for 15 minutes of continuous scanning.

Under the scan conditions provided above, the Zenith Low Profile AAA/ Zenith Alpha Abdominal Ancillary Components are not expected to result in a temperature rise of more than 1.3°C after 15 minutes of continuous scanning.

Image Artifact

The image artifact extends approximately 9.5 mm from Zenith Low Profile AAA/Zenith Alpha Abdominal Ancillary Components as found during nonclinical testing when imaged with a gradient echo pulse sequence and a 3.0 Tesla MRI system. The image artifact obscures the device lumen.

For US Patients Only

Cook recommends that the patient register the MR conditions disclosed in this IFU with the Medialert Foundation. The Medialert Foundation can be contacted in the following manners:

- Mail: Medialert Foundation International
2323 Colorado Avenue
Turlock, CA 95382
- Phone: 888-633-4298 (toll free)
209-668-3333 from outside the US
Fax: 209-669-2450
- Web: www.medialert.org

11.5 Additional Surveillance and Treatment

Additional surveillance and possible treatment is recommended for:

- Aneurysms with Type I endoleak
- Aneurysms with Type III endoleak
- Aneurysm enlargement, ≥5 mm of maximum diameter (regardless of endoleak status)
- Migration
- Inadequate seal length

Consideration for reintervention or conversion to open repair should include the attending physician's assessment of an individual patient's co-morbidities, life expectancy and the patient's personal choices. Patients should be counseled that subsequent reinterventions, including catheter based and open surgical conversion are possible following endograft placement.

12 REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

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© COOK 2020