EN
20

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Table of Contents

Illustrations.....	11-19
1 DEVICE DESCRIPTION.....	20
1.1 Zenith AAA Ancillary Components.....	20
1.2 Main Body Extension Delivery System.....	20
1.3 Iliac Leg Delivery System.....	20
1.4 Iliac Leg Extension Delivery System.....	20
1.5 Converter Delivery System.....	20
1.6 Iliac Plug Delivery System.....	20
2 INDICATIONS FOR USE.....	20
3 CONTRAINDICATIONS.....	20
4 WARNINGS AND PRECAUTIONS.....	20
4.1 General.....	20
4.2 Patient Selection, Treatment and Follow-Up.....	20
4.3 Pre-Procedure Measurement Techniques and Imaging.....	20
Diameters.....	20
Lengths.....	20
4.4 Device Selection.....	21
4.5 Implant Procedure.....	21
4.6 Molding Balloon Use.....	21
4.7 Main Body Extension.....	21
4.8 MRI Information.....	21
Static Magnetic Field.....	21
MRI-Related Heating.....	21
Image Artifact.....	21
5 ADVERSE EVENTS.....	21
5.1 Observed Adverse Events.....	21
5.2 Potential Adverse Events.....	21
Device-Related Adverse Event Reporting.....	21
6 SUMMARY OF CLINICAL STUDIES.....	21
7 PATIENT SELECTION AND TREATMENT.....	21
7.1 Individualization of Treatment.....	21
8 PATIENT COUNSELING INFORMATION.....	22
9 HOW SUPPLIED.....	22
10 CLINICAL USE INFORMATION.....	22
10.1 Physician Training.....	22
10.2 Inspection Prior to Use.....	22
10.3 Materials Required.....	22
10.4 Materials Recommended.....	22
10.5 Device Diameter Sizing Guidelines.....	22
Table 10.5.1 Main Body Extension Graft Diameter Sizing Guide.....	22
Table 10.5.2 Flex Iliac Leg Graft Diameter Sizing Guide.....	22
Table 10.5.3 Spiral-Z AAA Iliac Leg Graft Sizing Guide.....	23
Table 10.5.4 Iliac Leg Extension Graft Diameter Sizing Guide.....	23
Table 10.5.5 Converter Graft Diameter Sizing Guide.....	23
Table 10.5.6 Zenith Iliac Plug Graft Diameter Sizing Guide.....	23
11 DIRECTIONS FOR USE.....	23
Anatomical Requirements.....	23
General Use Information.....	23
Pre-Implant Determinants.....	23
Patient Preparation.....	23
11.1 Ancillary Components General Use Information.....	23
11.2 Main Body Extensions.....	23
Main Body Extension Preparation/Flush.....	23
Main Body Extension Placement and Deployment.....	23
Main Body Extension Molding Balloon Insertion.....	24
11.3 Iliac Legs and Iliac Leg Extensions.....	24
Iliac Leg Extension Preparation/Flush.....	24
Iliac Leg Extension Placement and Deployment.....	24
Iliac Leg Extension Molding Balloon Insertion.....	24
11.4 Converter.....	24
Converter Preparation/Flush.....	24
Converter Placement and Deployment.....	24
Converter Molding Balloon Insertion.....	24
11.5 Iliac Plug.....	24
Iliac Plug Preparation/Flush.....	24
Iliac Plug Placement.....	24
Iliac Plug Deployment.....	24
Iliac Plug Molding Balloon Insertion.....	24
Femoral-to-Femoral Crossover.....	25
12 IMAGING GUIDELINES AND POST-OPERATIVE FOLLOW-UP.....	25
12.1 General.....	25
12.2 MRI Information.....	25
Static Magnetic Field.....	25
MRI-Related Heating.....	25
Image Artifact.....	25
12.3 Additional Surveillance and Treatment.....	25
13 PATIENT TRACKING INFORMATION.....	25

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
1. Hub
2. Pin Vise
3. Peel-Away® Sheath
4. Stopcock
5. Connecting Tube
6. Flexor® Introducer Sheath
7. Captor® Hemostatic Valve
8. Gray Positioner
9. Inner Cannula
10. Gripper
11. Trigger-Wire Release
12. Iliac Plug

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
Fig. 31


Fig. 32

ZENITH® AAA ANCILLARY COMPONENTS WITH Z-TRAK™ INTRODUCTION SYSTEM

Read all instructions carefully. Failure to properly follow the instructions, warnings and precautions may lead to serious surgical consequences or injury to the patient.

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

CAUTION: All contents of the outer pouch (including the introduction system and the endovascular grafts) are supplied sterile, for single use only.

For the Zenith product line there are five applicable Suggested Instructions for Use (IFU). This IFU describes the Suggested Instructions for Use for the Zenith AAA Endovascular Graft Ancillary Components (main body extension, iliac leg, iliac leg extension, converter and iliac plug). For information on other Zenith components, please refer to the following Suggested Instructions for Use:

- Zenith AAA Endovascular Graft (Zenith AAA Endovascular Graft main body, iliac legs, iliac leg extension, main body extension, converter and occluder);
- Zenith Flex® AAA Endovascular Graft (Zenith Flex AAA Endovascular Graft main body and iliac legs);
- Zenith Renu® AAA Ancillary Graft (main body extension and converter configurations);
- Coda® Balloon Catheter; and
- Zenith® Spiral-Z™ AAA Iliac Leg

1 DEVICE DESCRIPTION

1.1 Zenith AAA Ancillary Components

Zenith AAA Ancillary Components include main body extensions, iliac legs, iliac leg extensions, converters and iliac plugs. (**Fig. 1**) These components are constructed from the same materials as the principal graft modules. All modules are constructed of full-thickness woven polyester fabric sewn to self-expanding stainless steel Cook-Z® stents with braided polyester and monofilament polypropylene suture. The modules are fully stented to provide stability and the expansile force necessary to open the lumen of the graft during deployment. Additionally, the Cook-Z stents provide the necessary attachment and seal of the graft to the vessel wall. The aortic main body extensions, iliac legs and iliac leg extensions can be used to provide additional length to their respective portions of the endovascular graft. The converters and iliac plugs can be used to convert a bifurcated graft into an aorto-uniiliac graft, if necessary (e.g., cases of Type III endoleak, limb occlusion or unattainable contralateral limb cannulation).

All Zenith AAA Endovascular Graft Ancillary Components feature a Flexor Introducer Sheath which resists kinking and is hydrophilically coated. For added hemostasis, the Captor Hemostatic Valve can be loosened or tightened for the introduction and/or removal of ancillary devices into and out of the sheath. Both features are intended to enhance trackability in the iliac arteries and abdominal aorta.

1.2 Main Body Extension Delivery System

Main body extensions utilize 18 and 20 French Z-Trak Introduction Systems. (**Fig. 2**) The main body extension introduction system contains a single trigger-wire release mechanism. The trigger wires secure the endovascular graft onto the delivery system until released by the physician. Unlike the Zenith AAA main body delivery system, no top cap is included with the dilator tip because the main body extension component does not contain a barbed, uncovered suprarenal stent. Deployment of the main body extension is achieved by sheath retraction and trigger-wire removal.

1.3 Iliac Leg Delivery System

The Zenith AAA iliac legs are shipped preloaded onto the Z-Trak Introduction system. (**Fig. 3**) The delivery system is designed for ease of use with minimal preparation. The iliac leg delivery system uses a 14 or 16 French Z-Trak Introduction System. All systems are compatible with a .035 inch wire guide. Deployment is achieved by sheath retraction.

1.4 Iliac Leg Extension Delivery System

Iliac leg extensions are deployed using 14 and 16 French Z-Trak Introduction Systems similar to those used to deploy the Zenith AAA iliac legs. (**Fig. 4**) These delivery systems have no trigger-wire release mechanism or top cap assembly. Deployment is achieved by sheath retraction.

1.5 Converter Delivery System

Converters utilize 18 and 20 French Z-Trak Introduction Systems. (**Fig. 5**) The converter introduction system has no trigger-wire release mechanism or top cap assembly. Deployment of the converter is achieved by sheath retraction.

1.6 Iliac Plug Delivery System

Iliac plugs are preloaded into 14 and 16 French Z-Trak Introduction Systems. (**Fig. 6**) The iliac plug delivery system contains a single trigger-wire release mechanism. Deployment of the iliac plug is achieved by sheath retraction and trigger-wire removal.

2 INDICATIONS FOR USE

The Zenith AAA Ancillary Components with the Z-Trak Introduction System are indicated for use with the Zenith AAA graft during either a primary or a secondary procedure in patients who have adequate iliac/femoral access compatible with the required introduction systems.

3 CONTRAINDICATIONS

The Zenith AAA Ancillary Components are contraindicated in:

- Patients with known sensitivities or allergies to stainless steel, polyester, solder (tin, silver), nitinol, polypropylene or gold.
- Patients with a systemic infection who may be at increased risk of endovascular graft infection.

4 WARNINGS AND PRECAUTIONS

4.1 General

- Read all instructions carefully. Failure to properly follow the instructions, warnings and precautions may lead to serious consequences or injury to the patient.
- Always have a qualified surgery team available during implantation or reintervention procedures in the event that conversion to open surgical repair is necessary.
- The Zenith AAA Ancillary Components with the Z-Trak Introduction System should only be used by physicians and teams trained in vascular interventional techniques (catheter-based and surgical) and in the use of these devices. Specific training expectations are described in **Section 10.1, Physician Training**.
- Additional endovascular interventions or conversion to standard open surgical repair following endovascular repair should be considered for patients experiencing an enlarging aneurysm, unacceptable decrease in fixation length (vessel and component overlap) and/or endoleak. An increase in aneurysm size and/or persistent endoleak or migration may lead to aneurysm rupture.
- Patients experiencing reduced blood flow through the graft limb and/or leaks may be required to undergo secondary interventions or surgical procedures.
- The Zenith AAA Ancillary Components have not been evaluated for use with grafts other than Zenith AAA stent grafts.

4.2 Patient Selection, Treatment and Follow-Up

- The Zenith AAA main body extension is designed to treat aortic neck diameters no smaller than 18 mm and no larger than 32 mm. The Zenith AAA main body extension is designed to treat proximal aortic necks (distal to the lowest renal artery) of at least 15 mm in length. Iliac artery distal fixation site greater than 10 mm in length and 7.5-20 mm in diameter (measured outer wall to outer wall) is required. These sizing measurements are critical to the performance of the endovascular repair.
- Adequate iliac or femoral access is required to introduce the device into the vasculature.
- Key anatomical elements that may affect successful exclusion of the aneurysm include severe proximal neck angulation (>60 degrees for infrarenal neck to axis of AAA or >45 degrees for suprarenal neck relative to the immediate infrarenal neck); short proximal aortic neck (<15 mm); an inverted funnel shape (greater than 10% increase in diameter over 15 mm of proximal aortic neck length); short distal fixation site (<10 mm); and circumferential thrombus and/or calcification at the arterial implantation sites, specifically the proximal aortic neck and distal iliac artery interface. In the presence of anatomical limitations, a longer neck may be required to obtain adequate sealing and fixation. Irregular calcification and/or plaque may compromise the fixation and sealing of the implantation sites. Necks exhibiting these key anatomic elements may be more conducive to graft migration.
- Access vessel diameter (measured inner wall to inner wall) and morphology (minimal tortuosity, occlusive disease and/or calcification) should be compatible with vascular access techniques and delivery systems of the profile of a 14 French to 20 French vascular introducer sheath. Vessels that are significantly calcified, occlusive, tortuous or thrombus-lined may preclude placement of the ancillary components and/or may increase the risk of embolization. A vascular conduit technique may be necessary to achieve success in some patients.
- Size of the pre-existing Zenith AAA Endovascular Graft should be assessed before selecting a main body extension.
 - Use of a 58 mm length main body extension should only be considered when the length from the lowest renal artery to the bifurcation of the pre-existing graft is >65 mm for TFB (Zenith AAA Endovascular Graft) main body grafts and >73 mm for TFFB (Zenith Flex AAA Endovascular Graft) main body grafts. Failure to use the device in this manner may result in increased risk of inadvertent occlusion of renal arteries.
 - Use of 73 mm length main body extension should only be considered when the length from the lowest renal artery to the bifurcation of the pre-existing graft is >86 mm for TFB (Zenith Flex AAA Endovascular Graft) main body grafts. Failure to use the device in this manner may result in increased risk of inadvertent occlusion of renal arteries.
- The Zenith AAA Ancillary Components with the Z-Trak Introduction System are not recommended in patients who cannot tolerate contrast agents necessary for intraoperative and postoperative follow-up imaging. All patients should be monitored closely and checked periodically for a change in the condition of their disease and the integrity of the endoprosthesis.
- The Zenith AAA Ancillary Components with the Z-Trak Introduction System are not recommended in patients exceeding weight and/or size limits which compromise or prevent the necessary imaging requirements.
- Inability to maintain patency of at least one internal iliac artery or occlusion of an indispensable inferior mesenteric artery may increase the risk of pelvic/ bowel ischemia.
- Multiple large, patent lumbar arteries, mural thrombus and a patent inferior mesenteric artery may all predispose a patient to Type II endoleaks. Patients with uncorrectable coagulopathy may also have an increased risk of Type II endoleak or bleeding complications.
- The safety and effectiveness of the Zenith AAA Ancillary Components with the Z-Trak Introduction System have not been evaluated in the following patient populations:
 - traumatic aortic injury
 - leaking, pending rupture or ruptured aneurysms
 - mycotic aneurysms
 - pseudoaneurysms resulting from previous graft placement
 - revision of previously placed endovascular grafts
 - uncorrectable coagulopathy
 - indispensable mesenteric artery
 - genetic connective tissue disease (e.g., Marfans or Ehlers-Danlos Syndromes)
 - concomitant thoracic aortic or thoracoabdominal aneurysms
 - active systemic infections
 - pregnant or nursing females
 - morbidly obese patients
 - less than 18 years of age
- Successful patient selection requires specific imaging and accurate measurements; please see **Section 4.3 Pre-Procedure Measurement Techniques and Imaging**.
- All lengths and diameters of the devices necessary to complete the procedure should be available to the physician, especially when pre-operative case planning measurements (treatment diameters/lengths) are not certain. This approach allows for greater intraoperative flexibility to achieve optimal procedural outcomes.

4.3 Pre-Procedure Measurement Techniques and Imaging

- Lack of non-contrast CT imaging may result in failure to appreciate iliac or aortic calcification, which may preclude access or reliable device fixation and seal.
- Pre-procedure imaging reconstruction thicknesses >3 mm may result in sub-optimal device sizing, or in failure to appreciate focal stenoses from CT.
- Clinical experience indicates that contrast-enhanced spiral computed tomographic angiography (CTA) with 3-D reconstruction is the strongly recommended imaging modality to accurately assess patient anatomy prior to treatment with the Zenith AAA Ancillary Components. If contrast-enhanced spiral CTA with 3-D reconstruction is not available, the patient should be referred to a facility with these capabilities.
- Clinicians recommend positioning the x-ray C-arm during procedural angiography such that the origins of the renal arteries, and particularly the lowest patent renal artery, are well demonstrated prior to deployment of the proximal edge of the graft material (sealing stent) of the main body extension. Additionally, angiography should demonstrate the iliac artery bifurcations such that the distal common iliacs are well defined relative to the origin of the internal iliac arteries bilaterally, prior to deployment of the iliac leg components or iliac leg extensions.

Diameters

Utilizing CT, diameter measurements should be determined from the outer wall to outer wall vessel diameter (not lumen measurement) to help with proper device sizing and device selection. The contrast-enhanced spiral CT scan must start 1 cm superior to the celiac axis and continue through the femoral heads at an axial thickness slice of 3 mm or less.

Lengths

Utilizing CT, length measurements should be determined to accurately assess infrarenal proximal neck length as well as planning main body sizes, iliac legs and ancillary components for the Zenith (Flex) AAA Endovascular Graft. These reconstructions should be performed in sagittal, coronal and 3-D.

- **The long-term performance of endovascular grafts has not yet been established. All patients should be advised that endovascular treatment requires life-long, regular follow-up to assess their health and the performance of their endovascular graft.** Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysm or changes in the structure or position of the endovascular graft) should receive enhanced follow-up. Specific follow-up guidelines are described in **Section 12, IMAGING GUIDELINES AND POSTOPERATIVE FOLLOW-UP**.

- The Zenith AAA Ancillary Components with the Z-Trak Introduction System are not recommended in patients unable to undergo, or who will not be compliant with, the necessary preoperative and postoperative imaging and implantation studies as described in the Zenith (Flex) AAA Endovascular Graft Instructions for Use. A copy is available online at www.cookmedical.com.
- After endovascular graft placement, patients should be regularly monitored for perigraft flow, aneurysm growth or changes in the structure or position of the endovascular graft. At a minimum, annual imaging is required, including: 1) abdominal radiographs to examine device integrity (separation between components or stent fracture) and 2) contrast and non-contrast CT to examine aneurysm changes, perigraft flow, patency, tortuosity and progressive disease. If renal complications or other factors preclude the use of image contrast media, abdominal radiographs and duplex ultrasound may provide similar information.

4.4 Device Selection

Strict adherence to the Zenith AAA Endovascular Graft Ancillary Components IFU sizing guide is strongly recommended when selecting the appropriate device size (**Tables 10.5.1 through 10.5.6**). Appropriate device oversizing has been incorporated into the IFU sizing guide. Sizing outside of this range can result in endoleak, fracture, migration, device infolding or compression.

4.5 Implant Procedure

(Refer to **Section 11, DIRECTIONS FOR USE**)

- Appropriate procedural imaging is required to successfully position the Zenith AAA Ancillary Components and assure accurate apposition to the aortic wall.
- Do not bend or kink the delivery system. Doing so may cause damage to the delivery system and the Zenith AAA Ancillary Component.
- To avoid any twist in the endovascular graft, during any rotation of the delivery system, be careful to rotate all components of the system together (from outer sheath to inner cannula).
- Do not continue advancing any portion of the delivery system if resistance is felt during advancement of the wire guide or delivery system. Stop and assess the cause of resistance; vessel, catheter or graft damage may occur. Exercise particular care in areas of stenosis, intravascular thrombosis or in calcified or tortuous vessels.
- Inadvertent partial deployment or migration of the endoprosthesis may require surgical removal.
- Unless medically indicated, do not deploy the Zenith AAA Ancillary Components in a location that will occlude arteries necessary to supply blood flow to organs or extremities. Do not cover significant renal or mesenteric arteries (exception is the inferior mesenteric artery) with the endoprosthesis. Vessel occlusion may occur. During the clinical study this device was not studied in patients with two occluded internal iliac arteries.
- Do not attempt to re-sheath the graft after partial or complete deployment.
- Repositioning the stent graft distally after partial deployment of the covered proximal stent may result in damage to the stent graft and/or vessel injury.
- Inaccurate placement and/or incomplete sealing of the Zenith AAA Ancillary Components within the vessel may result in increased risk of endoleak, migration or inadvertent occlusion of the renal or internal iliac arteries. Renal artery patency must be maintained to prevent/reduce the risk of renal failure and subsequent complications.
- Inadequate fixation of the Zenith AAA Ancillary Components may result in increased risk of migration of the stent graft. Incorrect deployment or migration of the endoprosthesis may require surgical intervention.
- Systemic anticoagulation should be used during the implantation procedure based on hospital and physician preferred protocol. If heparin is contraindicated, an alternative anticoagulant should be considered.
- Minimize handling of the constrained endoprosthesis during preparation and insertion to decrease the risk of endoprosthesis contamination and infection.
- Maintain wire guide position during delivery system insertion.
- Always use fluoroscopy for guidance, delivery and observation of any Zenith AAA Ancillary Component within the vasculature.
- The use of the Zenith AAA Ancillary Components with the Z-Trak Introduction System requires administration of intravascular contrast. Patients with preexisting renal insufficiency may have an increased risk of renal failure postoperatively. Care should be taken to limit the amount of contrast media used during the procedure and to observe preventative methods of treatment to decrease renal compromise (e.g., adequate hydration).
- As the sheath and/or wire guide is withdrawn, anatomy and graft position may change. Constantly monitor graft position and perform angiography to check position as necessary.
- Use caution during manipulation of catheters, wires and sheaths within an aneurysm. Significant disturbances may dislodge fragments of thrombus, which can cause distal embolization, or rupture of the aneurysm.
- Care should be taken not to displace or damage the main body graft during placement and deployment.
- The converter and main body extension delivery systems cannot be introduced through a main body or iliac leg introducer sheath.

4.6 Molding Balloon Use

- Do not inflate balloon in the vessel outside of graft, as doing so may cause damage to the vessel. Use the balloon in accordance with its labeling.
- Use care in inflating the balloon within the graft in the presence of calcification, as excessive inflation may cause damage to the vessel.
- Confirm complete deflation of the balloon prior to repositioning.

4.7 Main Body Extension

Care should be taken not to displace the main body graft during the placement and deployment of the main body extension.

4.8 MRI Information

The MR safety and compatibility of the Zenith AAA Ancillary Components is based upon MR safety and compatibility of the Zenith AAA Endovascular Graft, which is made from stents of the same metal.

Non-clinical testing has demonstrated that the Zenith AAA Endovascular Graft is MR Conditional. A patient with this endovascular graft can be scanned safely immediately after placement under the following conditions:

Static Magnetic Field

- Static magnetic field of 3.0 Tesla or less.
- Highest spatial magnetic gradient field of 720 Gauss/cm

Non-clinical evaluation was conducted in a 3 Tesla MR system (General Electric Excite) with a maximum spatial magnetic gradient field of 720 Gauss/cm as measured with a gaussmeter in the position of the static magnetic field pertinent to the patient (i.e., outside of scanner covering, accessible to a patient or individual).

MRI-Related Heating

1.5 Tesla Systems:

- Static magnetic field of 1.5 Tesla
- Maximum whole-body-averaged specific absorption rate (SAR) of 2.0 W/kg for 15 minutes of scanning (i.e., per scanning sequence)

In non-clinical testing, the Zenith AAA Endovascular Graft produced a temperature rise of less than or equal to 1.9 °C at a MR system reported whole-body-averaged specific absorption rate (SAR) of 2.8 W/kg, for 15 minutes of MR scanning in a 1.5 Tesla Magnetom, Siemens Medical Magnetom, Numaris/4 Software, Version Syngo MR 2002B DHHS MR Scanner. The maximum whole-body-averaged specific absorption rate (SAR) was 2.8 W/kg, which corresponds to a calorimetry measured value of 1.5 W/kg.

3.0 Tesla Systems:

- Static magnetic field of 3.0 Tesla
- Maximum whole-body-averaged specific absorption rate (SAR) of 2.0 W/kg for 15 minutes of scanning (i.e., per scanning sequence)

In non-clinical testing, the Zenith AAA Endovascular Graft produced a temperature rise of less than or equal to 1.9 °C at a MR system reported whole-body-averaged specific absorption rate (SAR) of 3.0 W/kg, for 15 minutes of MR scanning in a 3.0 Tesla Excite, GE Electric Healthcare, G3.0-052B Software, MR Scanner. The maximum whole-body-averaged specific absorption rate (SAR) was 3.0 W/kg, which corresponds to a calorimetry measured value of 2.8 W/kg.

Image Artifact

The image artifact extends throughout the anatomical region containing the device, obscuring the view of immediately adjacent anatomical structures within approximately 20 cm of the device, as well as the entire device and its lumen, when scanned in non-clinical testing using the sequence: Fast spin echo, in a 3.0 Tesla, Excite, GE Electric Healthcare, with G3.0-052B Software, MR system with body radiofrequency coil.

For all scanners, the image artifact dissipates as the distance from the device to the area of interest increases. MR scans of the head and neck and lower extremities may be obtained without image artifact. Image artifact may be present in scans of the abdominal region and upper extremities, depending on distance from the device to the area of interest.

Clinical information is available for seventeen patients who received MRI scans after stent-graft implantation. There have been no reported adverse events or device problems in any of these patients as a result of having received an MRI. Additionally, there have been well over 100,000 Zenith AAA Endovascular Grafts implanted worldwide, in which there have been no reported adverse events or device problems as a result of MRI.

Cook recommends that the patient register the MR conditions disclosed in this IFU with the MediAlert Foundation. The MediAlert Foundation can be contacted in the following manners:

Mail:	MediAlert Foundation International 2323 Colorado Avenue Turlock, CA 95382
Phone:	888-633-4298 (toll free) 209-668-3333 from outside the US
Fax:	209-669-2450
Web:	www.medialert.org

5 ADVERSE EVENTS

5.1 Observed Adverse Events

For information on observed adverse events in patients receiving the Zenith (Flex) AAA Endovascular Graft, some of whom received Zenith AAA Ancillary Components, please refer to the Zenith (Flex) AAA Endovascular Graft Instructions for Use. A copy is available online at www.cookmedical.com.

5.2 Potential Adverse Events

Adverse events that may occur and/or require intervention include, but are not limited to:

- Amputation
- Anesthetic complications and subsequent attendant problems (e.g., aspiration)
- Aneurysm enlargement
- Aneurysm rupture and death
- Aortic damage, including perforation, dissection, bleeding, rupture and death
- Arterial or venous thrombosis and/or pseudoaneurysm
- Arteriovenous fistula
- Bleeding, hematoma or coagulopathy
- Bowel complications (e.g., ileus, transient ischemia, infarction, necrosis)
- Cardiac complications and subsequent attendant problems (e.g., arrhythmia, myocardial infarction, congestive heart failure, hypotension, hypertension)
- Claudication (e.g., buttock, lower limb)
- Death
- Edema
- Embolization (micro and macro) with transient or permanent ischemia or infarction
- Endoleak
- Endoprosthesis: improper component placement; incomplete component deployment; component migration; component separation from another graft component; suture break; occlusion; infection; stent fracture; graft material wear; dilatation; erosion; puncture; perigraft flow; and corrosion
- Fever and localized inflammation
- Genitourinary complications and subsequent attendant problems (e.g., ischemia, erosion, fistula, incontinence, hematuria, infection)
- Graft or native vessel occlusion
- Hepatic failure
- Impotence
- Infection of the aneurysm, device or access site, including abscess formation, transient fever and pain
- Lymphatic complications and subsequent attendant problems (e.g., lymph fistula)
- Neurologic local or systemic complications and subsequent attendant problems (e.g., stroke, transient ischemic attack, paraplegia, paraparesis, paralysis)
- Pulmonary/respiratory complications and subsequent attendant problems (e.g., pneumonia, respiratory failure, prolonged intubation)
- Renal complications and subsequent attendant problems (e.g., artery occlusion, contrast toxicity, insufficiency, failure)
- Surgical conversion to open repair
- Vascular access site complications, including infection, pain, hematoma, pseudoaneurysm, arteriovenous fistula
- Vascular spasm or vascular trauma (e.g., iliofemoral vessel dissection, bleeding, rupture, death)
- Vessel damage
- Wound complications and subsequent attendant problems (e.g., dehiscence, infection)

Device-Related Adverse Event Reporting

Any adverse event (clinical incident) involving the Zenith AAA Ancillary Components should be reported to Cook immediately. To report an incident, call the Customer Relations Department at 1-800-457-4500 (24 hour) or 1-812-339-2235. For customers outside the United States, please call your distributor.

6 SUMMARY OF CLINICAL STUDIES

For clinical study information regarding patients receiving the Zenith (Flex) AAA Endovascular Graft, some of whom received Zenith AAA Ancillary Components, please refer to the Zenith (Flex) AAA Endovascular Graft Instructions for Use. A copy is available online at www.cookmedical.com.

7 PATIENT SELECTION AND TREATMENT

(See **Section 4, WARNINGS AND PRECAUTIONS**)

7.1 Individualization of Treatment

Cook recommends that the Zenith AAA Ancillary Component diameters be selected as described in **Tables 10.5.1 through 10.5.6**. All lengths and diameters of the devices necessary to complete the procedure should be available to the physician, especially when preoperative case planning

measurements (treatment diameters/lengths) are not certain. This approach allows for greater intraoperative flexibility to achieve optimal procedural outcomes. The risks and benefits should be carefully considered for each patient before use of the Zenith AAA Ancillary Components. Additional considerations for patient selection include, but are not limited to:

- Patient's age and life expectancy.
- Co-morbidities (e.g., cardiac, pulmonary or renal insufficiency prior to surgery, morbid obesity).
- Patient's suitability for open surgical repair.
- Patient's anatomical suitability for endovascular repair.
- The risk of aneurysm rupture compared to the risk of treatment with the Zenith AAA Ancillary Components.
- Patient's ability to tolerate general, regional or local anesthesia.
- Iliofemoral access vessel size and morphology (minimal thrombus, calcification and/or tortuosity) should be compatible with vascular access techniques and accessories of the delivery profile of a 14 French to 20 French vascular introducer sheath.
- For main body extension, non-aneurysmal infrarenal aortic segment (neck) proximal to the aneurysm:
 - with a diameter measured outer wall to outer wall of no greater than 32 mm and no less than 18 mm,
 - with an angle less than 60 degrees relative to the long axis of the aneurysm, and
 - with an angle less than 45 degrees relative to the axis of the suprarenal aorta.
- Iliac artery distal fixation site greater than 10 mm in length and 7.5 to 20 mm in diameter (measured outer wall to outer wall).
- Freedom from significant femoral/iliac artery occlusive disease that would impede flow through the endovascular graft.

The final treatment decision is at the discretion of the physician and patient.

8 PATIENT COUNSELING INFORMATION

- The physician and patient (and/or family members) should review the risks and benefits when discussing this endovascular device and procedure, including:
- Risks and differences between endovascular repair and surgical repair.
 - Potential advantages of traditional open surgical repair.
 - Potential advantages of endovascular repair.
 - The possibility that subsequent interventional or open surgical repair of the aneurysm may be required after initial endovascular repair.

In addition to the risks and benefits of an endovascular repair, the physician should assess the patient's commitment and compliance to postoperative follow-up as necessary to ensure continuing safe and effective results. Listed below are additional topics to discuss with the patient as to expectations after an endovascular repair:

- **The long-term performance of endovascular grafts has not yet been established. All patients should be advised that endovascular treatment requires life-long, regular follow-up to assess their health and the performance of their endovascular graft.** Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms or changes in the structure or position of the endovascular graft) should receive enhanced follow-up. Specific follow-up guidelines are described in **Section 12, IMAGING GUIDELINES AND POSTOPERATIVE FOLLOW-UP**.
- Patients should be counseled on the importance of adhering to the follow-up schedule, both during the first year and at yearly intervals thereafter. Patients should be told that regular and consistent follow-up is a critical part of ensuring the ongoing safety and effectiveness of endovascular treatment of AAAs. At a minimum, annual imaging and adherence to routine postoperative follow-up requirements are required and should be considered a life-long commitment to the patient's health and well-being.
- The patient should be told that successful aneurysm repair does not arrest the disease process. It is still possible to have associated degeneration of vessels.
- Physicians must advise each patient that it is important to seek prompt medical attention if he/she experiences signs of limb occlusion, aneurysm enlargement or rupture. Signs of graft limb occlusion include pain in the hip(s) or leg(s) during walking or at rest or discoloration or coolness of the leg. Aneurysm rupture may be asymptomatic, but usually presents as: pain; numbness; weakness in the legs; any back, chest, abdominal or groin pain; dizziness; fainting; rapid heartbeat or sudden weakness.
- Due to the imaging required for successful placement and follow-up of endovascular devices, the risks of radiation exposure to developing tissue should be discussed with women who are or suspect they are pregnant.
- Men who undergo endovascular or open surgical repair may experience impotence.

Physicians should refer the patient to the Patient Guide regarding risks occurring during or after implantation of the device. Procedure-related risks include cardiac, pulmonary, neurologic, bowel and bleeding complications. Device-related risks include occlusion, endoleak, aneurysm enlargement, fracture, potential for reintervention and open surgical conversion, rupture and death (See **Section 5, ADVERSE EVENTS**). The physician should complete the Patient I.D. Card and give it to the patient so that he/she can carry it with him/her at all times. The patient should refer to the card anytime he/she visits additional health practitioners, particularly for any additional diagnostic procedures (e.g., MRI).

9 HOW SUPPLIED

- The Zenith AAA Ancillary Components are sterilized by ethylene oxide gas, pre-loaded into the Z-Trak Introduction System, and supplied in peel-open packages.
- The devices are intended for one-time use only. Do not re-sterilize the devices.
- The product is sterile if the package is unopened and undamaged. Inspect the device and packaging to verify that no damage has occurred as a result of shipping. Do not use this device if damage has occurred or if the sterilization barrier has been damaged or broken. If damage has occurred, do not use the product and return to Cook.
- Prior to use, verify correct devices (quantity and size) have been supplied for the patient by matching the device to the order prescribed by the physician for that particular patient.
- Do not use after the "USE BY" (expiration) date printed on the label.
- Store in a cool, dry place.

10 CLINICAL USE INFORMATION

10.1 Physician Training

CAUTION: Always have a qualified surgery team available during implantation or re-intervention procedures in the event that conversion to open surgical repair is necessary.

CAUTION: The Zenith AAA Ancillary Components with the Z-Trak Introduction Systems should only be used by physicians and teams trained in vascular interventional techniques and in the use of this device. The recommended skill/knowledge requirements for physicians using the Zenith AAA Ancillary Components with the Z-Trak Introduction System are outlined below:

Patient selection:

- Knowledge of the natural history of abdominal aortic aneurysms (AAA) and co-morbidities associated with AAA repair.
- Knowledge of radiographic image interpretation, device selection, planning and sizing.

A multidisciplinary team that has combined procedural experience with:

- Femoral cutdown, arteriotomy and repair
- Percutaneous access and closure techniques
- Non-selective and selective wire guide and catheter techniques
- Fluoroscopic and angiographic image interpretation
- Embolization
- Angioplasty
- Endovascular stent placement
- Snare techniques
- Appropriate use of radiographic contrast material
- Techniques to minimize radiation exposure
- Expertise in necessary patient follow-up modalities

10.2 Inspection Prior to Use

Inspect the device and packaging to verify that no damage has occurred as a result of shipping. Do not use this device if damage has occurred or if the sterilization barrier has been damaged or broken. If damage has occurred, do not use the product and return to Cook. Prior to use, verify correct devices (quantity and size) have been supplied for the patient by matching the device to the order prescribed by the physician for that particular patient.

10.3 Materials Required

- Fluoroscope with digital angiography capabilities (C-arm or fixed unit)
- Contrast media
- Syringe
- Heparinized saline solution
- Sterile gauze pads

10.4 Materials Recommended

The following products are recommended. For information on the use of these products, refer to the individual product's instructions for use.

- .035 inch (0.89 mm) extra stiff wire guide, 260 cm; for example:
 - Cook Amplatz Ultra Stiff Wire Guides (AUS2)
 - Cook Lunderquist Extra Stiff Wire Guides (LES)
- .035 inch (0.89 mm) standard wire guide; for example:
 - Cook .035 inch wire guides
 - Cook Nimble™ Wire Guides
- Molding Balloons; for example:
 - Cook Coda Balloon Catheter
- Introducer sets; for example:
 - Cook Check-Flo® Introducer Sets
 - Cook Extra Large Check-Flo Introducer Sets
 - Cook Flexor® Balkin Up & Over® Contralateral Introducers
- Sizing catheter; for example:
 - Cook Auros® Centimeter Sizing Catheters
- Angiographic radiopaque tip catheters; for example:
 - Cook Beacon® Tip Angiographic Catheters
 - Cook Beacon Tip Royal Flush Catheters
- Entry needles; for example:
 - Cook single wall entry needles
- Endovascular Dilators; for example:
 - Cook Endovascular Dilator Sets

10.5 Device Diameter Sizing Guidelines

The choice of diameter should be determined from the outer wall to outer wall vessel diameter and **not** the lumen diameter. Undersizing or oversizing may result in incomplete sealing or compromised flow.

Table 10.5.1 Main Body Extension Graft Diameter Sizing Guide*			
Intended Aortic Vessel Diameter ^{1,2} (mm)	Main Body Extension Diameter ³ (mm)	Main Body Extension Length (mm)	Introducer Sheath (Fr)
18-19	22	39, 58	18
20-21	24	39, 58	18
22-23	26	39, 58	18
24-25	28	39, 58	20
26-27	30	39, 58	20
28	32	39, 58	20
29-32	36	50, 73	20

¹Maximum diameter along the proximal fixation site.
²Round measured aortic diameter to nearest mm.
³Additional considerations may affect choice of diameter.
*All dimensions are nominal.

Table 10.5.2 Flex Iliac Leg Graft Diameter Sizing Guide*			
Intended Iliac Vessel Diameter ^{1,2} (mm)	Iliac Leg Diameter ³ (mm)	Iliac Leg Working Length ⁴ (mm)	Introducer Sheath (Fr)
<8	8	37, 54, 71, 88, 105, 122	14
8-9	10	37, 54, 71, 88, 105, 122	14
10-11	12	39, 56, 73, 90, 107, 124	16
12-13	14	39, 56, 73, 90, 107, 124	16
14-15	16	39, 56, 73, 90	16
16-17	18	39, 56, 73, 90	16
18	20	39, 56, 73, 90	16
19	22	39, 56, 73, 90	16
20	24	39, 56, 73, 90	16

¹Maximum diameter along the distal fixation site.
²Round measured iliac diameter to nearest mm.
³Additional considerations may affect choice of diameter.
⁴Overall leg length = working length + 22 mm docking stent.
*All dimensions are nominal.

Table 10.5.3 Spiral-Z AAA Iliac Leg Graft Sizing Guide*			
Intended Iliac Vessel Diameter ^{1,2} (mm)	Intended Iliac Vessel Diameter ³ (mm)	Iliac Leg Working Length ⁴ (mm)	Introducer Sheath (Fr)
≤8	9	39, 56, 74, 90, 107,122	14
9	11	39, 56, 74, 90, 107,122	14
10-12	13	39, 56, 74, 90, 107,122	14
13-15	16	39, 56, 74, 90	14
16-18	20	39, 56, 74, 90	16
19-20	24	39, 56, 74, 90	16

¹Maximum diameter along the proximal fixation site.
²Round measured iliac diameter to nearest mm.
³Additional considerations may affect choice of diameter.
⁴Overall leg length = working length + 22 mm docking stent.
 *All dimensions are nominal.

Table 10.5.4 Iliac Leg Extension Graft Diameter Sizing Guide*			
Intended Iliac Vessel Diameter ^{1,2} (mm)	Iliac Leg Extension Diameter ³ (mm)	Iliac Leg Extension Length (mm)	Introducer Sheath (Fr)
<8	8	55	14
8-9	10	55	14
10-11	12	55	14
12-13	14	55	14
14-15	16	55	14
16-17	18	55	16
18	20	55	16
19	22	55	16
20	24	55	16

¹Maximum diameter along the distal fixation site.
²Round measured iliac diameter to nearest mm.
³Additional considerations may affect choice of diameter.
 *All dimensions are nominal.

Table 10.5.5 Converter Graft Diameter Sizing Guide*			
Main Body Diameter (mm)	Converter Diameter ¹ (mm)	Converter Length (mm)	Introducer Sheath (Fr)
22	24	80	18
24	24	80	18
26	28	80	20
28	28	80	20
30	32	80	20
32	32	80	20
36	36	82	20

¹Additional considerations may affect choice of diameter.
 *All dimensions are nominal.

Table 10.5.6 Zenith Iliac Plug Graft Diameter Sizing Guide*			
Intended Iliac Vessel Diameter ^{1,2} (mm)	Iliac Plug Diameter ³ (mm)	Iliac Plug Length (mm)	Introducer Sheath (Fr)
8-10	14	30	14
11-12	16	30	14
13-16	20	30	16
17-20	24	30	16

¹Maximum diameter along the distal fixation site.
²Round iliac diameter to nearest mm.
³Additional considerations may affect choice of diameter.
 *All dimensions are nominal.

11 DIRECTIONS FOR USE

Anatomical Requirements

- Iliofemoral access vessel size and morphology (minimal thrombus, calcium and/or tortuosity) should be compatible with vascular access techniques and accessories. Arterial conduit techniques may be required.
- For use of the main body extension, proximal aortic necks should be a minimum of 15 mm with a diameter measured outer wall to outer wall of 18-32 mm. Iliac artery distal fixation site should be greater than 10 mm in length and 7.5-20 mm in diameter (measured outer wall to outer wall).

Prior to use of the Zenith AAA Ancillary Components with the Z-Trak Introduction System, review this *Suggested Instructions for Use* booklet. The following instructions embody a basic guideline for device placement. Variations in the following procedures may be necessary. These instructions are intended to help guide the physician and do not take the place of physician judgement.

General Use Information

- Standard techniques for placement of arterial access sheaths, guiding catheters, angiographic catheters and wire guides should be employed during use of the Zenith AAA Ancillary Components with the Z-Trak Introduction System.
- Endovascular stent grafting is a surgical procedure, and blood loss from various causes may occur, infrequently requiring intervention (including

transfusion) to prevent adverse outcomes. It is important to monitor blood loss from the hemostatic valve throughout the procedure, but is specifically relevant during and after manipulation of the gray positioner. After the gray positioner has been removed, if blood loss is excessive, consider placing an uninflated molding balloon or an introduction system dilator within the valve, restricting flow.

Pre-Implant Determinants

Verify from pre-implant planning that the correct device has been selected. Determinants include:

1. Femoral artery selection for introduction of the delivery system.
2. Angulation of aortic neck, aneurysm and iliac arteries.
3. Quality of the aortic neck, if main body extension is used.
4. Diameters of infrarenal aortic neck and distal iliac arteries.
5. Distance from renal arteries to the bifurcation of a previously placed bifurcated graft, if main body extension is to be used.
6. Length from the bifurcation of a previously placed bifurcated graft to the internal iliac arteries/attachment site(s).
7. Aneurysm(s) extending into the iliac arteries may require special consideration in selecting a suitable graft/artery interface site.
8. Degree of vascular calcification.

Patient Preparation

If using Zenith AAA Ancillary Components as part of a secondary procedure:

1. Refer to institutional protocols relating to anesthesia, anticoagulation and monitoring of vital signs.
 2. Position patient on imaging table allowing fluoroscopic visualization from the aortic arch to the femoral bifurcations.
 3. Expose selected common femoral artery using standard surgical technique. When a femoral-to-femoral crossover will be required as part of the procedure, both common femoral arteries should be exposed.
 4. Establish adequate proximal and distal vascular control of selected femoral vessel.
- NOTE:** The Zenith AAA Ancillary Components are designed to be introduced through an exposed common femoral artery on the selected side of introduction. Angiography at the implant site may be achieved by the use of a straight angiographic catheter inserted from the contralateral side, either by surgical exposure or percutaneous approach.
5. Puncture the selected common femoral arteries using standard technique with an 18 or 19 UT gage arterial needle. Upon vessel entry, insert:
 - Wire guides – standard .035 inch diameter, 145 cm long, J tip or Benton Wire Guide
 - Appropriate size sheaths (e.g., 6 or 8 French)
 - Flush catheter (often radiopaque sizing catheters – e.g., Centimeter Sizing Catheter or straight flush catheter)
 6. Perform angiography to identify level(s) of renals, bifurcation of the previously placed bifurcated graft and iliac bifurcations.

NOTE: If fluoroscope angulation is used with an angulated neck, it may be necessary to perform angiograms using various projections.

11.1 Ancillary Components General Use Information

Inaccuracies in device size selection or placement, changes or anomalies in patient anatomy or procedural complications can require placement of additional endovascular grafts, extensions, iliac plugs and converters. Regardless of the device placed, the basic procedure(s) will be similar to the maneuvers required for Zenith AAA stent grafts. It is vital to maintain wire guide access.

Standard techniques for placement of arterial access sheaths, guiding catheters, angiographic catheters and wire guides should be employed during use of the Zenith AAA Ancillary Components. The Zenith AAA Ancillary Components with the Z-Trak Introduction System are compatible with .035 inch diameter wire guides.

If necessary to increase hemostasis, an uninflated molding balloon or an introduction system dilator may be placed within the hemostatic valve.

11.2 Main Body Extensions

Main body extensions (**Fig. 2**) are used for extending the proximal body of an *in situ* endovascular graft.

Main Body Extension Preparation/Flush

1. Remove gray-hubbed shipping stylet (from the inner cannula) and dilator tip protector (from the dilator tip). Remove Peel-Away sheath from back of the hemostatic valve. (**Fig. 7**) Elevate distal tip of system and flush through the stopcock on the hemostatic valve until fluid emerges from the sideport near the tip of the introducer sheath. (**Fig. 8**) Continue to inject a full 20 cc of flushing solution through the device. Discontinue injection and close stopcock.

NOTE: Graft flushing solution of heparinized saline is often used.

2. Attach syringe with heparinized saline to the hub on the inner cannula. Flush until fluid exits the dilator tip. (**Fig. 9**)

NOTE: When flushing system, elevate distal end of system to facilitate removal of air.

3. Soak sterile gauze pads in saline solution and use to wipe Flexor introducer sheath to activate the hydrophilic coating. Hydrate both sheath and dilator liberally.

Main Body Extension Placement and Deployment

1. Replace J tip wire guide with stiff wire guide (LES), .035 inch, 260 cm long, and advance through catheter and up to the thoracic aorta. Remove flush catheter and sheath. Maintain wire guide position.
2. An angiographic catheter should be introduced via the contralateral femoral artery and positioned at the level of the required implant site.
3. Introduce the main body extension delivery system into the ipsilateral femoral artery.

CAUTION: The main body extension delivery system cannot be introduced through a main body or iliac leg introducer sheath.

4. Advance slowly until the main body extension is at the site of the required intervention. (**Fig. 10**)
5. Verify main body extension position to ensure proper sealing and resistance to migration.
6. Verify placement with angiography to ensure that the renal arteries remain patent and that proper placement is achieved.

CAUTION: Care should be taken not to displace the main body graft during the placement and deployment of the main body extension.

NOTE: Ensure the Captor Hemostatic Valve on the Flexor Introducer sheath is turned to the open position. (**Fig.11**)

7. Use the gripper to stabilize the gray positioner while withdrawing the sheath. (**Figs. 12 and 13**) Continue to deploy the device until the most distal stent is uncovered. (**Fig. 14**) Stop withdrawing the sheath.
8. Remove the safety lock from the trigger-wire release mechanism. Withdraw and remove the trigger-wire by sliding the trigger-wire release mechanism off the handle, and then remove via its slot over the inner cannula. (**Fig. 15**)
9. Withdraw the tapered tip of the introducer back through the main body extension graft and sheath while maintaining wire guide position. Ensure the main body extension and endovascular graft are not displaced during withdrawal of delivery system.
10. Close the Captor Hemostatic Valve on the main body extension introducer sheath by turning it in a clockwise direction until it stops. (**Fig. 16**)

11. Re-check the position of the wire guide. Leave sheath and wire guide in place.

Main Body Extension Molding Balloon Insertion

NOTE: For information on the use of recommended products (listed on page 22) refer to the individual product's instructions for use.

1. Prepare molding balloon as follows:
 - Flush wire lumen with heparinized saline.
 - Remove all air from balloon.
2. In preparation for the insertion of the molding balloon, open the Captor Hemostatic Valve by turning counterclockwise.
3. Advance the molding balloon over the wire guide and through the hemostatic valve of the main body extension introduction system to the level of the main body extension.

CAUTION: Do not inflate balloon in vessel outside of graft.

4. Tighten the Captor Hemostatic Valve around the molding balloon with gentle pressure, by turning it clockwise.
5. Expand the molding balloon within the most proximal segment of the main body extension and then the most distal segment of the main body extension using diluted contrast media (as recommended by the manufacturer). (**Fig. 17**)

CAUTION: Confirm complete deflation of balloon prior to repositioning.

CAUTION: Captor Hemostatic Valve must be open prior to repositioning the molding balloon.

6. Completely deflate and remove molding balloon, replace it with an angiographic catheter and perform completion angiograms.
7. If no other endovascular maneuvers are necessary, remove any sheaths, wires and catheters. Repair vessels and close in standard surgical fashion.

11.3 Iliac Legs and Iliac Leg Extensions

The iliac leg extensions (**Fig. 4**) are used for extending the distal iliac legs and/or bridging an *in situ* endovascular graft. If an iliac leg extension is not available, an iliac leg (**Fig. 3**) can be used. If using an iliac leg/iliac leg extension as part of primary repair, the iliac leg/iliac leg extension delivery system may be introduced through an existing *in situ* 18, 20 or 22 French introduction system sheath.

Iliac Leg Extension Preparation/Flush

1. Remove the gray-hubbed shipping stylet (from the inner cannula) and dilator tip protector (from the dilator tip). Remove Peel-Away sheath from the back of the hemostatic valve. (**Fig. 18**) Elevate distal tip of system and flush through the stopcock on the hemostatic valve until fluid emerges from the flushing groove near the tip of the introducer sheath. (**Fig. 8**) Continue to inject a full 20 cc of flushing solution through the device. Discontinue injection and close stopcock.

NOTE: Graft flushing solution of heparinized saline is often used.

2. Attach syringe with heparinized saline to the hub on the inner cannula. Flush until fluid exits the dilator tip. (**Fig. 9**)

NOTE: When flushing system, elevate distal end of system to facilitate removal of air.

3. Soak sterile gauze pads in saline solution and use to wipe Flexor introducer sheath to activate the hydrophilic coating. Hydrate both sheath and dilator liberally.

Iliac Leg Extension Placement and Deployment

1. Replace J tip wire guide with stiff wire guide (LES), .035 inch, 260 cm long, and advance through catheter and up to the thoracic aorta. Remove flush catheter and sheath. Maintain wire guide position.
2. Introduce the iliac leg extension delivery system into the artery.

NOTE: If the Zenith Flex AAA Endovascular Graft main body delivery system is being utilized to introduce the iliac leg extension, ensure the Captor Hemostatic Valve is in the open position prior to introduction and deployment of the iliac leg extension.

3. Advance slowly until the iliac leg extension is at the site of the required intervention. (**Fig. 19**) Ensure that a proper placement is achieved. Verify appropriate stent graft overlap to ensure proper sealing and resistance to migration.
4. Verify placement with angiography to ensure the internal iliacs will remain patent.

NOTE: Ensure the Captor Hemostatic Valve on the Flexor introducer sheath is turned to the open position. (**Fig. 11**)

5. Use the gripper to stabilize the gray positioner while withdrawing the sheath. (**Figs. 12 and 20**)
6. Continue to deploy the device until the distal stent is uncovered. (**Fig. 21**) Stop withdrawing the sheath.
7. Withdraw the tapered tip of the introducer back through the iliac leg extension graft and sheath while maintaining wire guide position. Ensure the endovascular graft is not displaced during withdrawal of delivery system.

Iliac Leg Extension Molding Balloon Insertion

NOTE: For information on the use of recommended products (listed on page 22) refer to the individual product's instructions for use.

1. Prepare molding balloon as follows:
 - Flush wire lumen with heparinized saline.
 - Remove all air from balloon.
2. In preparation for the insertion of the molding balloon, open the Captor Hemostatic Valve by turning counterclockwise.
3. Advance the molding balloon over the wire guide and through the hemostatic valve to the most proximal segment of the iliac leg extension.

CAUTION: Do not inflate balloon in vessel outside of graft.

4. Tighten the Captor Hemostatic Valve around the molding balloon with gentle pressure, by turning it clockwise.
5. Expand the molding balloon within the most proximal segment and then the most distal segment of the iliac leg extension using diluted contrast media (as recommended by the manufacturer). (**Fig. 22**)

CAUTION: Captor Hemostatic Valve must be open prior to repositioning the molding balloon.

CAUTION: Confirm complete deflation of balloon prior to repositioning.

6. Completely deflate and remove molding balloon, replace it with an angiographic catheter and perform completion angiograms.
7. If no other endovascular maneuvers are necessary, remove any sheaths, wires and catheters. Repair vessels and close in standard surgical fashion.

11.4 Converter

Converters can be used to convert a bifurcated graft into an aorto-uniliac graft, if necessary (e.g., cases of Type III endoleak, limb occlusion or unattainable contralateral limb cannulation). (**Fig. 5**)

Converter Preparation/Flush

1. Remove gray-hubbed shipping stylet (from the inner cannula) and dilator tip protector (from the dilator tip). Remove Peel-Away sheath from back of the hemostatic valve. (**Fig. 18**) Elevate distal tip of system and flush through the stopcock on the hemostatic valve until fluid emerges from the sideport near the tip of the introducer sheath. (**Fig. 8**) Continue to

inject a full 20 cc of flushing solution through the device. Discontinue injection and close stopcock.

NOTE: Graft flushing solution of heparinized saline is often used.

2. Attach syringe with heparinized saline to the hub on the inner cannula. Flush until fluid exits the dilator tip. (**Fig. 9**)

NOTE: When flushing system, elevate distal end of system to facilitate removal of air.

3. Soak sterile gauze pads in saline solution and use to wipe Flexor introducer sheath to activate hydrophilic coating. Hydrate both sheath and dilator liberally.

Converter Placement and Deployment

1. Replace J tip wire guide with stiff wire guide (LES), .035 inch, 260 cm long, and advance through catheter and up to the thoracic aorta. Remove flush catheter and sheath. Maintain wire guide position.
2. Introduce the converter delivery system into the artery.

CAUTION: The converter delivery system cannot be introduced through a main body or iliac leg introducer sheath.

3. Advance slowly until the converter is at the site of the required intervention. (**Fig. 23**) Verify appropriate stent graft overlap to ensure proper sealing and resistance to migration. The proximal two stents should be positioned in the main body graft, and the distal two stents should be positioned in the ipsilateral leg.

NOTE: Ensure the Captor Hemostatic Valve on the Flexor introducer sheath is turned to the open position.

4. Use the gripper to stabilize the gray positioner while withdrawing the sheath. (**Figs. 12 and 24**)
5. Continue to deploy the device until the distal stent is uncovered. (**Fig. 25**)
6. Withdraw the tapered tip of the introducer back through the converter graft and sheath while maintaining wire guide position. Ensure the endovascular graft is not displaced during withdrawal of delivery system.

Converter Molding Balloon Insertion

NOTE: For information on the use of recommended products (listed on page 22) refer to the individual product's instructions for use.

1. Prepare molding balloon as follows:
 - Flush wire lumen with heparinized saline.
 - Remove all air from balloon.
2. In preparation for the insertion of the molding balloon, open the Captor Hemostatic Valve by turning counterclockwise.
3. Advance the molding balloon over the wire guide and through the hemostatic valve to the most proximal segment of the converter.
4. Tighten the Captor Hemostatic Valve around the molding balloon with gentle pressure by turning it clockwise.

CAUTION: Do not inflate balloon in vessel outside of graft.

5. Expand the molding balloon within the most proximal segment and then the most distal segment of the converter using diluted contrast media (as recommended by the manufacturer). (**Fig. 26**)

CAUTION: Captor Hemostatic Valve must be open prior to repositioning the molding balloon.

CAUTION: Confirm complete deflation of balloon prior to repositioning.

6. Completely deflate and remove molding balloon, replace it with an angiographic catheter and perform completion angiograms.
7. If no other endovascular maneuvers are necessary, remove any sheaths, wires and catheters. Repair vessels and close in standard surgical fashion.

11.5 Iliac Plug

The Iliac Plug (**Fig. 6**) is used to occlude an iliac artery, typically in conjunction with a femoral-to-femoral crossover procedure.

Iliac Plug Preparation/Flush

1. Remove the gray-hubbed shipping stylet (from the inner cannula). Remove Peel-Away sheath from the back of the hemostatic valve. (**Fig. 27**) Elevate distal tip of system and flush through the stopcock on the hemostatic valve until fluid exits the dilator tip. (**Fig. 28**) Continue to inject a full 20 cc of flushing solution through the device. Discontinue injection and close stopcock.

NOTE: Graft flushing solution of heparinized saline is often used.

2. Attach syringe with heparinized saline to the hub on the inner cannula. Flush until fluid exits the dilator tip. (**Fig. 29**)

NOTE: When flushing system, elevate distal end of system to facilitate removal of air.

3. Soak sterile gauze pads in saline solution and use to wipe Flexor introducer sheath to activate hydrophilic coating. Hydrate both sheath and dilator liberally.

Iliac Plug Placement

1. Perform angiography to determine proper iliac plug placement (below aortic bifurcation and above iliac bifurcation).
2. Replace J wire with stiff wire guide (LES), .035 inch, at least 145 cm long, and advance through catheter and above the aortic bifurcation, being careful not to disturb any previously placed graft.
3. Insert the delivery system over the wire into the femoral artery. (**Fig. 30**)
4. Advance the delivery system to the intended position in the common iliac artery.
5. Verify placement with angiography of the iliac plug relative to the internal iliac artery.

CAUTION: Maintain wire guide position during delivery system insertion.

NOTE: The iliac plug is tethered to the gray positioner, which allows the iliac plug to be repositioned within the vessel.

Iliac Plug Deployment

NOTE: Ensure the Captor Hemostatic Valve on the Flexor introducer sheath is turned to the open position.

1. Verify position of wire guide in the iliac artery. Deploy the iliac plug by withdrawing the sheath while stabilizing the gray positioner of the delivery system with the gripper. Withdraw the sheath until the iliac plug is uncovered. (**Figs. 31 and 32**) Stop withdrawing sheath.
2. Remove the safety lock from the trigger-wire release mechanism. Withdraw and remove the trigger-wire by sliding the trigger-wire release mechanism off the handle, and then remove via its slot over the inner cannula. (**Fig. 33**)
3. Perform angiography to verify position of the iliac plug with respect to the internal iliac artery.
4. Withdraw the wire guide while maintaining position of the gray positioner. Ensure the endovascular graft is not displaced during withdrawal of wire guide.
5. Slowly withdraw gray positioner. Ensure the endovascular graft is not displaced during withdrawal of gray positioner.

Iliac Plug Molding Balloon Insertion

NOTE: The expandable length of the Zenith AAA Iliac Plug is 30 mm. To avoid inflation of molding balloon outside the distal aspect of the plug, the balloon should fit within the expandable length of the plug.

1. Prepare molding balloon as follows:
 - Flush wire lumen with heparinized saline.
 - Remove all air from balloon.
2. In preparation for insertion of the molding balloon, open the Captor Hemostatic Valve by turning counterclockwise.
3. Re-introduce the J tip wire guide. Advance the wire guide until it curves inside the body of the iliac plug. Advance the molding balloon over the wire guide and through the hemostatic valve into the iliac plug.

CAUTION: Do not inflate balloon in vessel outside of graft.

4. Tighten the Captor Hemostatic Valve around the molding balloon with gentle pressure, by turning it clockwise.
5. Expand the molding balloon within the iliac plug using diluted contrast media (as directed by the manufacturer). (**Fig. 34**)

CAUTION: Confirm complete deflation of balloon prior to repositioning.

CAUTION: Captor Hemostatic Valve must be open prior to repositioning the molding balloon.

6. Completely deflate and remove molding balloon, replace it with an angiographic catheter and perform completion angiograms to verify proper occlusion of the appropriate common iliac artery.

Femoral-to-Femoral Crossover

Perform femoral-to-femoral crossover in standard surgical fashion to revascularize the contralateral limb. Repair vessels and close in standard surgical fashion.

12 IMAGING GUIDELINES AND POSTOPERATIVE FOLLOW-UP

For information on imaging guidelines and postoperative follow-up, please refer to the Zenith Flex AAA Endovascular Graft Instructions for Use. A copy is available online at www.cookmedical.com.

12.1 General

- **The long-term performance of endovascular grafts with secondary endovascular intervention using additional components has not yet been established. All patients should be advised that endovascular treatment requires life-long, regular follow-up to assess their health and the performance of their endovascular graft.** Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms or changes in the structure or position of the endovascular graft) should receive additional follow-up. Patients should be counseled on the importance of adhering to the follow-up schedule, both during the first year and at yearly intervals thereafter. Patients should be told that regular and consistent follow-up is a critical part of ensuring the ongoing safety and effectiveness of endovascular treatment of AAAs.
- Physicians should evaluate patients on an individual basis and prescribe follow-up relative to the needs and circumstances of each individual patient. The minimum requirement for patient follow-up (described in the Instructions for Use for the Zenith Flex AAA Endovascular Graft) should be maintained even in the absence of clinical symptoms (e.g., pain, numbness, weakness). Patients with specific clinical findings (e.g., endoleaks, enlarging aneurysms or changes in the structure or position of the stent graft) should receive follow-up at more frequent intervals.
- Annual imaging follow-up should include abdominal radiographs and both contrast and non-contrast CT examinations. If renal complications or other factors preclude the use of image contrast media, abdominal radiographs, non-contrast CT and duplex ultrasound may be used.
- The combination of contrast and non-contrast CT imaging provides information on aneurysm diameter change, endoleak, patency, tortuosity, progressive disease, fixation length and other morphological changes.
- The abdominal radiographs provide information on device integrity (e.g., separation between components, stent fracture).
- Duplex ultrasound imaging may provide information on aneurysm diameter change, endoleak, patency, tortuosity and progressive disease. In this circumstance, a non-contrast CT should be performed to use in conjunction with the ultrasound. Ultrasound may be a less reliable and sensitive diagnostic method compared to CT.
- The minimum imaging follow-up for patients with Zenith AAA stent grafts is described in the Instructions for Use for the Zenith (Flex) AAA Endovascular Graft, which can be found at www.cookmedical.com. Patients requiring enhanced follow-up should have interim evaluations.

12.2 MRI Information

The MR safety and compatibility of the Zenith AAA Ancillary Components is based upon MR safety and compatibility of the Zenith AAA Endovascular Graft, which is made from stents of the same metal.

Non-clinical testing has demonstrated that the Zenith AAA Endovascular Graft is MR Conditional. A patient with this endovascular graft can be scanned safely immediately after placement under the following conditions:

Static Magnetic Field

- Static magnetic field of 3.0 Tesla or less.
- Highest spatial magnetic gradient field of 720 Gauss/cm

Non-clinical evaluation was conducted in a 3 Tesla MR system (General Electric Excite) with a maximum spatial magnetic gradient field of 720 Gauss/cm as measured with a gaussmeter in the position of the static magnetic field pertinent to the patient (i.e., outside of scanner covering, accessible to a patient or individual).

MRI-Related Heating

1.5 Tesla Systems:

- Static magnetic field of 1.5 Tesla
- Maximum whole-body-averaged specific absorption rate (SAR) of 2.0 W/kg for 15 minutes of scanning (i.e., per scanning sequence)

In non-clinical testing, the Zenith AAA Endovascular Graft produced a temperature rise of less than or equal to 1.4 °C at a MRI system reported whole-body-averaged specific absorption rate (SAR) of 2.8 W/kg, for 15 minutes of MR scanning in a 1.5 Tesla Magnetom, Siemens Medical Magnetom, Numaris/4 Software, Version Syngo MR 2002B DHHS MR Scanner. The maximum whole-body-averaged specific absorption rate (SAR) was 2.8 W/kg, which corresponds to a calorimetry measured value of 1.5 W/kg.

3.0 Telsa Systems

- Static magnetic field of 3.0 Tesla
- Maximum whole-body-averaged specific absorption rate (SAR) of 2.0 W/kg for 15 minutes of scanning (i.e., per scanning sequence)

In non-clinical testing, the Zenith AAA Endovascular Graft produced a temperature rise of less than or equal to 1.9 °C at a MR system reported whole-body-averaged specific absorption rate (SAR) of 3.0 W/kg, for 15 minutes of MR scanning in a 3.0 Tesla Excite, GE Electric Healthcare, G3.0-052B Software, MR Scanner. The maximum whole-body-averaged specific absorption rate (SAR) was 3.0 W/kg, which corresponds to a calorimetry measured value of 2.8 W/kg.

Image Artifact

The image artifact extends throughout the anatomical region containing the device, obscuring the view of immediately adjacent anatomical structures within approximately 20 cm of the device, as well as the entire device and its lumen, when scanned in non-clinical testing using the sequence: Fast spin echo, in a 3.0 Tesla, Excite, GE Electric Healthcare, with G3.0-052B Software, MR system with body radiofrequency coil.

For all scanners, the image artifact dissipates as the distance from the device to the area of interest increases. MR scans of the head and neck and lower extremities may be obtained without image artifact. Image artifact may be present in scans of the abdominal region and upper extremities, depending on distance from the device to the area of interest.

Clinical information is available for seventeen patients who received MRI scans after stent-graft implantation. There have been no reported adverse events or device problems in any of these patients as a result of having received an MRI. Additionally, there have been well over 100,000 Zenith AAA Endovascular Grafts implanted worldwide, in which there have been no reported adverse events or device problems as a result of MRI.

Cook recommends that the patient register the MR conditions disclosed in this IFU with the MedAlert Foundation. The MedAlert Foundation can be contacted in the following manners:

Mail: MedAlert Foundation International
2323 Colorado Avenue
Turlock, CA 95382

Phone: 888-633-4298 (toll free)
209-668-3333 from outside the US

Fax: 209-669-2450

Web: www.medicalert.org

12.3 Additional Surveillance and Treatment

Additional surveillance and possible treatment is recommended for:

- Aneurysms with Type I endoleak
- Aneurysms with Type III endoleak
- Aneurysm enlargement, ≥5 mm of maximum diameter (regardless of endoleak status)
- Migration
- Inadequate seal length

Consideration for reintervention or conversion to open repair should include the attending physician's assessment of an individual patient's co-morbidities, life expectancy and the patient's personal choices. Patients should be counseled that the need for subsequent reinterventions, including catheter-based and open surgical conversion, are possible following endograft placement.

13 PATIENT TRACKING INFORMATION

In addition to these Instructions for Use, the Zenith AAA Ancillary Components with their respective introduction systems are packaged with a Device Tracking Form, which the hospital staff is required to complete and forward to Cook for the purposes of tracking all patients who receive the Zenith AAA Ancillary Components (as required by U.S. Federal Regulation).

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© COOK 2013
Bloomington, IN 47404 U.S.A.